



EZPZ-AP/SP ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

EZPZ AP/SP Mobile Shower Commode Chairs with contoured molded seat, height-adjustable powder-coated aluminum frame, flip-up arm supports, foot supports, dual-locking 5" (125mm) casters and commode pan. **User weight capacity 330 lbs.**

1. EZPZ Attendant Propel / Self Propel Non-Tilting Shower Commode Chairs(19"W seat)

EZ105	EZPZ-AP (Attendant Propel) mobile shower commode chair with 4 dual-locking 5" casters. Will fit over most toilets.	\$1,695
EZ205	EZPZ-SP (Self Propel) mobile shower commode chair with dual-locking 5" front casters and 24" rear wheels. Will fit over most toilets.	\$2,195

2. EZPZ Options and Accessories

EZSS	Soft Seat – removable, contoured foam seat pad	\$315
EZALK	Armrest Locks (Pair)	\$52
EZLTS2PR	Lateral Thoracic Support with Large Pad – Pair	\$740
EZCBASL	Chest Belt – Large	\$69
EZCBICASL	Chest Belt – Infection Control – Large	\$115
EZPBL	Pelvic Belt – Large	\$69
EZPBICL	Pelvic Belt – Infection Control – Large	\$115
EZCSIC19	Snap-on Calf Strap – Infection Control	\$110
ZAHL	Adjustable Heel Loops (Pair)	\$70
ZHS1820	H-Strap	\$105
ZNC	Neoprene Footplate Covers (Each)	\$55
ZIHW	Whizard – winged urine deflector	\$75
EZATSPR	Anti-Tipper for SP (Pair)	\$125
EZQRSP	Quick Release Axles (Pair)	\$55

Specifications	EZPZ-AP (Attendant Propel)	EZPZ-SP (Self Propel)
Weight Capacity	330 lbs	330 lbs
Seat Width	19.6"	19.6"
Seat Depth (front of back post to front of seat)	16.9"	16.9"
Back Width	18.5"	18.5"
Back Height (from top of seat to top of upholstery)	16.9"	16.9"
Toilet Clearance Height	15.3" - 19.2"	15.3" - 19.2"
Toilet Width Clearance	16.9"	16.9"
Arm Support Height Above Seat*	9"	9"
Seat-To-Floor Height* (1" increments)	18.8" - 22.8"	18.8" - 22.8"
Overall Width	22.8"	27.9"
Product Weight	35.2 lbs	39.6 lbs
Frame Material	Aluminum	Aluminum

* without soft seat overlay



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com
T (866) 875-2877