



EZPZ-AT ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

EZPZ-AT – Attendant Tilt Mobile Shower Commode Chair with 35° tilt range, height-adjustable aluminum frame, contoured composite seat, head support fixture with small pad, flip-up arm supports, foot supports, four dual-locking 5” casters and commode pan. **User weight capacity 330 lbs.**

1. EZPZ Tilting Shower Commode Chair (19"W seat)

EZ305	EZPZ-AT w/ 19"W mobile shower commode chair with dual-locking 5" casters. Will fit over most toilets.	\$2,695
-------	---	---------

2. EZPZ Options and Accessories

EZSS	EZPZ Soft Seat Overlay Pad	\$315
EZALK	EZPZ Armrest Locks (Pair)	\$50
EZLTS2PR	Lateral Thoracic Support with Large Pad – Pair	\$740
EZCBTL	Chest Belt – Large	\$65
EZCBICTL	Chest Belt – Infection Control – Large	\$110
EZPBL	Pelvic Belt – Large	\$65
EZPBICL	Pelvic Belt – Infection Control – Large	\$110
EZCSIC	Snap-on Calf Strap – Infection Control	\$110
ZHS1820	H-Strap Calf Support	\$105
ZAHL	Adjustable Heel Loops (Pair)	\$70
ZNC	Neoprene Footplate Covers (Each)	\$55
ZIHW	Whizard – winged urine deflector	\$75
EZATT	Anti-tippers	\$250

Specifications	EZPZ-AT (Attendant Tilt)
Weight Capacity	330 lbs
Tilt Range	35° tilt range including optional 5° of anterior tilt
Seat Width	19.5"
Seat Depth (front of back post to front of seat)	17"
Back Width	18.5"
Back Height (from top of seat to top of upholstery)	17"
Toilet Clearance Height	16.5" - 20.5"
Toilet Clearance Height - with anti-tippers	17.7" - 21.6"
Toilet Width Clearance	17"
Arm Support Height Above Seat*	9"
Seat-To-Floor Height* (1" increments)	20" - 24"
Seat-To-Floor Height* – with anti-tippers	21" - 25"
Overall Width	23"
Product Weight	40 lbs
Frame Material	Aluminum

* without soft seat overlay



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com