



RAZ-SP ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-SP: Self Propel Mobile Shower Commode Chair with height-adjustable stainless-steel frame, Fixed Back, flip-up padded arm supports, foot supports, 22" or 24" rear wheels, front dual-locking 125mm casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 5cm of fore/aft seat adjustment. **User weight capacity 185 Kgs.**

For help with justification, please see our [Evaluation and Medical Justification Form](#).

Need help with this order form? Call customer service at (866) 875-2877 or (519) 942 8441

1. Standard Raz-SP Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]

Z20018	Raz-SP w/ 18"W base frame – Standard base frame width, will fit over most toilets	\$2,805
--------	---	---------

2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.7]

Z20016	Raz-SP w/ 16"W frame width ¹ (base chair price – custom upcharge below will be added)	\$2,805
Z2SW16	Upcharge for custom 16"W frame – will not fit over most toilets	\$570
Z20020	Raz-SP w/ 20"W frame width ² (base chair price – custom upcharge below will be added)	\$2,805
Z2SW20	Upcharge for custom 20"W frame	\$570
Z20022	Raz-SP w/ 22"W frame width (base chair price – custom upcharge below will be added)	\$2,805
Z2SW22	Upcharge for custom 22"W W frame	\$690

Custom Frame Configurations

ZLNGF3	Upcharge for custom 3" Longer Frame	\$460
Z2LOWF4	Upcharge for custom Low Frame ³ (16"-18" STF, 4" lower than standard) – 22" rear wheels only	\$570

¹ 18" W Visco Foam Interface seats can be modified to fit.

² 18" W seats can be selected but must also select Seat Support Kit (ZSSKAS) [Sec.3]

³ Forward Offset Foot Support Receiver ZFOR [Sec. 11] is required when selecting seats that are 20" or deeper.

3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16" D and 19" D, are now designated 17" D and 20" D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Molded Seats (thickness ~6.5cm -2.75")			Visco-Foam Interface Seats (thickness ~3.5")		
ZMS17FRO	18" W x 17" D – Front Access	Standard	ZVFI17BRI	18" W x 17" D – Front Bridged	\$495
ZMS17BRI	18" W x 17" D – Front Bridged	N/C	ZVFI17FRO	18" W x 17" D – Front Access	\$495
ZMS20FRO	18" W x 20" D – Front Access	\$155	ZVFI17LTO	18" W x 17" D – Left Access	\$495
ZMS20BRI	18" W x 20" D – Front Bridged	\$155	ZVFI17RTO	18" W x 17" D – Right Access	\$495
Specialty Seats			ZVFI17RRO	18" W x 17" D – Rear Access	\$495
ZSOM	Shower-Only Seat – Molded 18" W x 16" -21" D	\$460	ZVFI20BRI	18" W x 20" D – Front Bridged	\$495
ZEJRM	E&J Replica Seat – Molded – 18" W x 17" D	\$225	ZVFI20FRO	18" W x 20" D – Front Access	\$570
ZEJRV	E&J Replica Seat – Visco Foam – 18" W x 17" D	\$970	ZVFI20LTO	18" W x 20" D – Left Access	\$570
			ZVFI20RTO	18" W x 20" D – Right Access	\$570
			ZVFI20RRO	18" W x 20" D – Rear Access	\$570
Seat Options					
ZSSKAS	Seat Support Kit for Raz-SP - These are used to support the sides of seats that are narrower than seat frames on which they are mounted. A typical example would be when a 20" W chair is needed but an 18" W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.				\$165
ZSS05	Spacers for Front Seat Brackets – 0.5" (pr) – for adding dump angle (based on 18" D seat)				\$120
ZSS10	Spacers for Front Seat Brackets – 0.75" (pr) – for adding dump angle (based on 18" D seat)				\$140

4. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. See [Custom Seat Assessment Form](#) – Resources tab at razdesigninc.com

ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below) Select if both width and depth are ≤ 21"	\$1,190
ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$1,360
ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if width and/or depth is ≥ 22"	\$1,545

Custom Specifications (for above custom seats)

Base Frame Width: __ 16" / 18" / 20" / 22"
 Custom Seat Width: __ 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than frame width)
 Custom Seat Depth: __ 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Long Frame (ZLNGF3)]
 Access Opening: __ 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)
 Access Opening Width: __ 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required
 Aperture Width:

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.

5. Fixed Back Frames (18" W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (90°/97°) with cane-style posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

ZFB1818	Fixed Back – 18" W x 18" H (Armrest Lock Kit (ZALK) [Sec.11] is not included.)	Standard
ZFBU1818	Tension-Adjustable Back Upholstery for Fixed Back	\$180

6. Adjusta-Back Frames – Back frame widths will be the same as the seat frame widths unless otherwise selected in Sec 8.

Includes height-and angle-adjustable arms, tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (90°/97°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 6.5cm higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. Adjusta-Back – 24"H includes straight posts, which do not allow the back upholstery to be raised. Does not include the Armrest Lock Kit (ZALK).



ZAB __ 15*	Adjusta-Back – 15"H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22")	\$285
ZAB __ 17*	Adjusta Back 17" * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22")	\$285
ZAB __ 20*	Adjusta-Back – 20"H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22")	\$370
ZAB __ 24*	Adjusta-Back – 24"H * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$420

Reclining Backs – 30° Reclining back frame width must be the same as the seat frame width. 2-Point Arms (Z2PLAS / Z2PSAS) are included with Reclining Backs.

ZRBS __ 20*	Reclining Adjusta-Back with 20"H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$1090
ZRBS __ 24*	Reclining Adjusta-Back with 24"H H back (anti-tippers included) * insert back frame width to complete the code (16" / 18" / 20" / 22")	\$1145

7. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.7] (configurations below are not available for recliner backs)

NOTE: For back frame width variances on recliner frames, please request custom quotation.

ZNB__ __ *	Custom Back Frame 2" Narrower than base frame	\$395	Base Frame Width	16"	18"	20"	56	
			Back Frame Width	14" ¹	16" ²	18" ³	20" ⁴	
* Enter base and back frame widths to complete code. eg. ZNB2018 specifies a 20"W chair frame with an 18"W back frame								
ZNB__ __ *	Custom Back Frame 4" Narrower than base frame – requires custom seat	\$395	Base Frame Width (cm)	18"	20"	22"	24"	
			Back Frame Width (cm)	14" ¹	16" ²	18" ³	20" ⁴	
* Enter seat and back frame widths to complete code. eg. ZNB2016 specifies a 18"W seat frame with a 20"W back frame								

¹ Requires Custom Seat

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit ZSSKAS [Sec.3].

⁴ Custom Seat required.

8. Back Support Options – These back supports replace the molded vinyl upholstery on a Fixed Back or tension-adjustable upholstery on an Adjusta-Back. Symphony and Harmony Backs are recommended for infection-control environments.

ZSB14	14" W Symphony Back ¹ – recessed solid back support	\$340
ZSB16	16" W Symphony Back ¹ – recessed solid back support	\$340
ZSB18	18" W Symphony Back ¹ – recessed solid back support	\$340
ZSB20	20" W Symphony Back ¹ – recessed solid back support	\$460
ZSB22	22" W Symphony Back ¹ – recessed solid back support	\$525
ZHB	Harmony Back ^{1,2} (contoured solid back – 18" W x 13" H) Fits 18" W back frames only Fits 17 / 20" back heights only	\$340

¹ Does not fit 15" H backs

² Not recommended for reclining backs.

9. Rear Wheels, Anti-Tippers and Caster Options – Rear Anti-Tippers Are Strongly Recommended

ZRW22	22" Diameter Wheels w/ Textured Handrims	Std	ZCSLSP	User-Actuated Caster Swivel Lock w/ extension handles (pr)	\$570
ZRW24	24" Diameter Wheels w/ Textured Handrims	N/C	ZATSPRR	Anti-Tippers – Rear (pr) standard and long frames ²	\$135
ZQRSP	Quick-Release Axles (pr)	\$55	ZATSPLWRR	Anti-Tippers – (pr) Rear for low frame ²	\$135
ZWLE	Wheel Lock Extension Kit ¹ (pr)	\$90	ZATAS	Anti-Tippers – Front (pr) ²	\$395
ZVH22	Vinyl-Coated Handrims for 22" (pr)	\$455	ZCUSP4	4" Caster Upgrade (2x)	\$70
ZVH24	Vinyl-Coated Handrims for 24" (pr)	\$490	ZCUSP6	6" Caster Upgrade (2x)	\$205
ZPH22	Vertical Projection Handrims – 22" (pr)	\$805			
ZPH24	Vertical Projection Handrims – 24" (pr)	\$765			

¹ Compatible with old-style aluminum Wheel Locks only

² Front Anti-Tippers provided according to caster diameter ordered.

10. Arm Supports and Accessories

ZARMTL	Long Arm Support Tubes (pr)	Standard
ZARMTS	Short Arm Support Tubes (pr) ¹	N/C
ZARMLL	Long Arm Support Pad 13" L – Left	Standard
ZARMLR	Long Arm Support Pad 13" L – Right	Standard
ZARMSL	Short Arm Support Pad 10" L – Left	N/C
ZARMSR	Short Arm Support Pad 10" L – Right	N/C
ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	\$90
Z2PLAS	2-Point Arm Supports ² (pr) with 13" L pads (std with reclining backs)	\$510
Z2PSAS	2-Point Arm Supports ² (pr) with 10" L pads (std with reclining backs)	\$525
ZFAPLL	Long Molded Flat Armpad ³ – Left (13.75" L x 3.15" W – fits ZARM13)	\$95
ZFAPLR	Long Molded Flat Armpad ³ – Right (13.75" L x 3.15" W – fits ZARM13)	\$95
ZFAPSL	Short Molded Flat Armpad ³ – Left (10" L x 2.75" W – fits ZARM10)	\$95
ZFAPSR	Short Molded Flat Armpad ³ – Right (10" L x 2.75" W – fits ZARM10)	\$95
ZTRSL	Arm Trough ⁴ - Small Left (16" L x 4.3" W) NEW	\$130
ZTRSR	Arm Trough ⁴ - Small Right (16" L x 4.3" W) NEW	\$130
ZTRML	Arm Trough ⁴ - Medium Left (47cm L x 10.75" W)	\$130
ZTRMR	Arm Trough ⁴ - Medium Right (47cm L x 10.75" W)	\$130
ZTRSL	Arm Trough ⁴ - Large Left (21" L x 5.1" W) NEW	\$130
ZTRSR	Arm Trough ⁴ - Large Right (21" L x 5.1" W) NEW	\$130
ZPAML	Pivoting Arm Mount ⁵ – Left (for ZFAPLL / ZFAPLR / ZTRSL / ZTRML / ZTRLL)	\$190
ZPAMR	Pivoting Arm Mount ⁵ – Right (for ZFAPLR / ZFAPLR / ZTRSR / ZTRMR / ZTRLR)	\$190
ZASK05	Arm Support Spacer Kit ⁶ – 0.5" per side for all Adjusta-Backs (pr)	\$70
ZASK10	Arm Support Spacer Kit ⁶ – 0.75" per side for all Adjusta-Backs (pr)	\$135
ZPHGL	Pivoting Hand Grips (pr) for 13" L arms ⁷	\$460
ZPHGS	Pivoting Hand Grips (pr) for 10" L arms ⁷	\$460
ZAPSB_*	Anterior Posture Support Bar ^{8,9} * Enter frame width to complete code (14" / 16" / 18" / 20")	\$225
ZAPSB_*	Anterior Posture Support Bar ^{8,9} * Enter frame width to complete code (22" / 24")	\$260
Reclining Back Arm Support Accessories		
Z2PLASR	2-Point Flip-up Arm Supports (Long) (pr)	N/C
Z2PSASR	2-Point Flip-up Arm Supports (Short) (pr)	N/C
ZASKR05	Arm Support Spacers 0.5" per side (pr) for Reclining Back Support	\$70
ZASKR10	Arm Support Spacers 0.75" per side (pr) for Reclining Back Support	\$70
¹ Not compatible with ZPAM / ZFAPL / ZFAPS / ZTRM / ZTRL ² For chairs with same back width and frame width only ³ May require ZPAM for Arm Support to flip up with Flat Arm Pad ⁴ ZPAM is required for Arm Support to flip up with Arm Trough ⁵ Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered. ⁶ ZASK05 or ZASK10 is required for Standard Arm Support to flip up with lateral pad in forward position ⁷ ZPHGL ZPHGS combined with 2 Point Arms requires a custom up charge ⁸ Not compatible with ZFAP / ZTR ⁹ Special length required with ZASK05 and ZASK10 (contact customer service)		

11. Foot and Leg Supports

ZFSL	Foot Support Left – 16" - 20"	Standard
ZFSR	Foot Support Right – 16" - 20"	Standard
ZFET2L	Footrest Extension Tube ¹ Left – 5cm longer	\$30
ZFET2R	Footrest Extension Tube ¹ Right – 5cm longer	\$30
ZMFX8L	MFX8 Foot Support Left – 8" - 12"	\$45
ZMFX8R	MFX8 Foot Support Right – 8" - 12"	\$45
ZMFX12L	MFX12 Foot Support Left – 12" - 16"	\$55
ZMFX12R	MFX12 Foot Support Right – 12" - 16"	\$55
ZADABS18L	Angle/Depth Adj ABS Footplate Left For 18" W frame	\$110
ZADABS18R	Angle/Depth Adj ABS Footplate Right For 18" W frame	\$110
ZADML	Adj Angle/Depth Molded Footplate ² - Left	\$110
ZADMR	Adj Angle/Depth Molded Footplate ² - Right	\$110
ZFET2AL	Footrest Extension Tube Left for ZADML – 5cm longer	\$30
ZFET2AR	Footrest Extension Tube Right for ZADMR – 5cm longer	\$30
ZVS_ _*	V-Style Foot Support (1-piece) 16" - 20" * Enter frame width to complete code (16" / 18")	\$340
ZVS_ _*	V-Style Foot Support (1-piece) 16" - 20" * Enter frame width to complete code (20" / 22")	\$375
ZFBAS	Flip-Back Foot Support (1-piece) 16" - 20"	\$340
ZELRL	Elevating Leg Support ³ – Left (includes ZADML)	\$570
ZELRR	Elevating Leg Support ³ – Right (includes ZADMR)	\$570
ZRLSL	Universal Residual Limb Support – Left	\$225
ZRLSR	Universal Residual Limb Support – Right	\$225
Foot Support Options for 16" W Frames		
ZFS16L	Foot Supports for 16" W frame Left – (16" - 20")	N/C
ZFS16R	Foot Supports for 16" W frame Right – (16" - 20")	N/C
ZADABS16L	Angle/Depth Adj ABS Footplate Left for 16" W frame	\$110
ZADABS16R	Angle/Depth Adj ABS Footplate Right For 16" W frame	\$110
ZMFX816L	MFX8 Foot Support for 16" W Frame Left – 8" - 12"	\$45
ZMFX816R	MFX8 Foot Support for 16" W Frame Right – 8" - 12"	\$45
ZMFX1216L	MFX12 Foot Support for 16" W Frame Left – 12" - 16"	\$55
ZMFX1216R	MFX12 Foot Support for 16" W Frame Right – 12" - 16"	\$55
¹ Not compatible with ZMFX12/ZADM/ZVS/ZFBT		
² Not compatible with Z30016		
³ ZELRL/R compatible with frame widths ≥ 18" and adjusted seat depths ≤ 18" or ≤ 21" w/ 2.75" longer seat frames (Z3LNGSF3)		

12. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$110
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (16" / 18" / 20" / 22")	\$155
ZCSIC_*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (16" / 18" / 20")	\$145
ZCSIC22	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code For 22" W frame	\$165
ZCPL	Calf Panel – Bodypoint Aeromesh ¹ Fits frame widths 16", 18" 20"	\$210
ZHS1820	H-Strap – for 18" / 20" W frames – hook and loop attachment	\$105
ZHS2224	H-Strap – for 22" / 24" W frames – hook and loop attachment	\$130
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$70
ZNC	Neoprene Cover for Standard Footplate (Each)	\$55
ZNCA	Neoprene Cover for Angle/Depth Adjustable (Each)	\$55
ZNC16	Neoprene Covers for Narrow Molded Footplates – 16" W frame (Each)	\$55
ZNCABS16	Neoprene Cover for Narrow Molded ABS Angle/Depth Adjustable (Each)	\$55
ZNCVS	Neoprene Cover for V-Style Foot Support	\$55
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$55
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 5cm laterally or medially. Specify 1 or 2 ___	\$135
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 5cm laterally or medially. Specify 1 or 2 ___	\$135
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 ___	\$140

¹ Not compatible with ZMFX12 / ZFBT / ZVS

13. Positioning Pads – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3" .

Head Supports			
ZHSAS1	Head Support with small pad		\$570
ZHSAS2	Head Support with large pad		\$630
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment		\$90
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 2.75" per side. For up to 4.3" medial offset, specify ZLAEP extension plate.			
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left		\$380
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right		\$380
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left		\$380
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right		\$380
Leg Adductors			
ZLGAS1L	Leg Adductor w/ Small Pad – Left		\$380
ZLGAS1R	Leg Adductor w/ Small Pad – Right		\$380
ZLGAS2L	Leg Adductor w/ Large Pad – Left		\$380
ZLGAS2R	Leg Adductor w/ Large Pad – Right		\$380
Hip Adductors			
ZHP1L	Hip Adductor w/ Small Pad – Left		\$380
ZHP1R	Hip Adductor w/ Small Pad – Right		\$380
ZHP2L	Hip Adductor w/ Large Pad – Left		\$380
ZHP2R	Hip Adductor w/ Large Pad – Right		\$380
Hip Adductors For Reclining Backs			
ZHPR1L	Hip Adductor for Reclining Back w/ Small Pad – Left		\$380
ZHPR1R	Hip Adductor for Reclining Back w/ Small Pad – Right		\$380
ZHPR2L	Hip Adductor for Reclining Back w/ Large Pad – Left		\$380
ZHPR2R	Hip Adductor for Reclining Back w/ Large Pad – Right		\$380
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors			
These plates allow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors			
ZLAEPL	Lateral/Adductor Extension Plate – Left		\$30
ZLAEPR	Lateral/Adductor Extension Plate – Right		\$30
Lateral Thoracic Support	Left	Right	
Leg Adductor	Left	Right	
Hip Adductor	Left	Right	
Total Quantity ____			
Medial Thigh Supports (Pommel)			
ZMTSAS1	Medial Thigh Support w/ Small Pads		\$690
ZMTSAS2	Medial Thigh Support w/ Large Pads		\$690

14. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBS	Pelvic Belt – Small (frame widths 13.75" - 16", 2-piece w/ side-release buckle)	\$50	ZCBS	Chest Belt ² – Small (back widths 13.75" - 16", 1-piece w/ hook and loop closure)	\$40
ZPBM	Pelvic Belt – Medium (frame widths 16" - 18", 2-piece w/ side-release buckle)	\$65	ZCBM	Chest Belt ² – Medium (back widths 16" - 18", 1-piece w/ hook and loop closure)	\$45
ZPBL	Pelvic Belt – Large (frame widths 18" - 22", 2-piece w/ side-release buckle)	\$65	ZCBL	Chest Belt ² – Large (back widths 18" - 22", 1-piece w/ hook and loop closure)	\$55
ZPBXL	Pelvic Belt – XLarge (seat frame widths 22" - 24", 2-piece w/ side-release buckle)	\$75	ZCBXL	Chest Belt – XLarge ² (back frame widths 56 / 24", 1-piece w/ hook and loop closure)	\$70
ZPBICS	Pelvic Belt ¹ – Infection Control, Small (frame widths 13.75" - 16", 2-piece w/ side-release buckle)	\$110	ZCBICS	Chest Belt ^{1,2} – Infection Control, Small 1-piece w/ side-release buckle (back widths 14-16")	\$90
ZPBICM	Pelvic Belt ¹ – Infection Control, Medium (frame widths 16" - 18", 2-piece w/ side-release buckle)	\$110	ZCBICM	Chest Belt ^{1,2} – Infection Control, Medium 1-piece w/ side-release buckle (back widths 16" - 18")	\$90
ZPBICL	Pelvic Belt ¹ – Infection Control, Large (frame widths 18" - 22", 2-piece w/ side-release buckle)	\$110	ZCBICL	Chest Belt ^{1,2} – Infection Control, Large 1-piece w/ side-release buckle (back widths 18" - 22")	\$100
ZPBICXL	Pelvic Belt ¹ – Infection Control, XLarge (seat frame widths 22" - 24", 2-piece w/ side-release buckle)	\$110	ZCBICXL	Chest Belt ^{1,2} – Infection Control, XLarge 1-piece w/ side-release buckle (back frame widths 22" / 24")	\$125
ZPBBPM	Pelvic Belt – Bodypoint Aeromesh, Med	\$295	ZCBBPM	Chest Belt – Bodypoint Aeromesh, Med	\$310
ZPBBPL	Pelvic Belt – Bodypoint Aeromesh, Large	\$295	ZCBBPL	Chest Belt – Bodypoint Aeromesh, Large	\$310
ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh, XLarge	\$225	ZCBBPXL	Chest Belt – Bodypoint Aeromesh, XLarge	\$310

Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPS	N/A	ZPBBPM	16" W base frame	ZPBBPL	18" / 20" W base frame	ZPBBPXL	22" / 24" W base frame
	Chest belt	ZCBBPS	13.75" back frame	ZCBBPM	16" W back frame	ZCBBPL	18" / 20" W back frame	ZCBBPXL	22" / 24" W back frame

¹ Non-absorbent belts – recommended for infection-control environments

² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM / ZCBL / ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

15. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$75
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$170
ZSC17	Splash Curtain ¹ for 17" D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 0.5").	\$345
ZSC20	Splash Curtain ¹ for 20" D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 0.5").	\$390
ZSCFR	Splash Curtain for 17" / 20" D (at front of aperture only, not for Front Open Seats)	\$405
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 0.5"	\$420
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$390
ZBUDBH	Bard urine drainage bag holder (available for all frame widths)	\$190

¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.

16. Other

ZPSK1	Commode Pan Spacer Kit ¹ – 0.75" lower	\$225
ZPSK2	Commode Pan Spacer Kit ¹ – 5cm lower	\$245
ZPSK3	Commode Pan Spacer Kit ¹ – 2.75" lower	\$285
ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$285

¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.

17. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907	Omit Foot Supports - if no other Foot / Leg Support Option is selected	-\$45
------	--	-------	------	--	-------

Some assembly of Raz products is required.

18. Notes:

--

Specifications	Raz-SP (Self Propel)	
	22" wheels	24" wheels
Weight Capacity	400 lbs	400 lbs
Overall Width	Seat width plus 7.75"	Seat width plus 7.75"
Overall Length (foot supports in middle position)	41"	42"
Seat Depth	16" or 21" depending on seat/cushion length	
Clearance Height (under seat)		
- with pan brackets	15.25"-19.25" (1" increments)	
- without pan brackets	16"- 20" (1" increments)	
Clearance Width (to fit around toilet) – 18"W frame	16"	16"
Arm Support Height - above molded seat (1/2" increments)		
- 18"H Fixed Back - standard	8.5"	8.5"
- 15"H Adjusta-Back	5.5"-8.5"	5.5"-8.5"
- 17"H Adjusta-Back	6.5"-9.5"	6.5"-9.5"
- 20"H Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 24"H Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 20"H Reclining Adjusta-Back	7.5"-10.5"	7.5"-10.5"
- 24"H Reclining Adjusta-Back	7.5"-10.5"	7.5"-10.5"
Seat-to-Floor Height (with molded seat)	20"- 24"	20"-24"
Product Weight (with 18"W seat)	47 lbs	49 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)	