



## RAZ-AT

Attendant Propel Mobile Shower Commode Chair



Bathing and Toileting Products for Life

# RAZ-AP FEATURES

The Raz-AP Attendant Propel Mobile Shower Commode Chair comes standard with 4" of height adjustment. The easy-to-clean, stainless steel frame fits over most toilets and is designed with unobstructed side and front access for the user or attendant.



### ALSO AVAILABLE

- Long frame – 3" longer for tall clients
- Low frame - 4" lower seat-to-floor height
- 16", 20" and 22" widths.
- Recliner frame

SPECIFICATIONS	RAZ-AP
WEIGHT CAPACITY	400 lbs
FRAME WIDTH (seat and back frame)	18"
OVERALL WIDTH (with Standard Armpad)	22.5"
OVERALL LENGTH (with footrest)	33.75"
SEAT WIDTH (Commode Seat)	18"
SEAT DEPTH (front of back post to front of seat)	16" - 19"
CLEARANCE HEIGHT UNDER SEAT	15.5" – 19.5" with pan brackets 16.25" – 20.25" without pan brackets
ARMREST HEIGHT – above seat (17" Adjusta-back)	6.5" - 9.5" in ½" increments
SEAT TO FLOOR HEIGHT	20" - 24" in ½" increments
PRODUCT WEIGHT (with foot supports)	38 lbs
FRAME AND HARDWARE MATERIAL	Stainless Steel – Life Time Warranty Against Corrosion

All measurements taken with standard molded seat and standard arm pad.

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice..