



RAZ-AP600

Heavy Duty Attendant Propel
Mobile Shower Commode Chairs



Bathing and Toileting Products for Life

RAZ-AP600 FEATURES

The Raz-AP600 [Attendant-Propel, 600-lb capacity] rehab shower commode chair features the Multi-Position Back with fore/aft center of mass adjustment to accommodate a variety of body contours.



SPECIFICATIONS	RAZ-AP600 (Bariatric Attendant Propel)
WEIGHT CAPACITY	600 lbs
SEAT WIDTH (seat and back frame)	22" - 30"
OVERALL WIDTH (with Standard Armpad)	26.5"
OVERALL LENGTH (with Z709 Footrests)	36" - 39"
SEAT WIDTH (Commode Seat)	23"
SEAT DEPTH (front of back post to front of seat)	18"
CLEARANCE HEIGHT UNDER SEAT	5" casters w/ pan brackets: 14.25" - 18.25" 5" casters w/o pan brackets: 15" - 19" 6" casters w/ pan brackets: 15.25" - 19.25" 6" casters without pan brackets: 16" - 20"
ARMREST HEIGHT ABOVE SEAT	9.5"
SEAT TO FLOOR HEIGHT	19" - 23" (in 1" increments)
PRODUCT WEIGHT	5" casters: 54.5 lbs 6" casters: 56.5 lbs
FRAME AND HARDWARE MATERIAL	Stainless Steel – Life Time Warranty Against Corrosion

All measurements taken with standard molded seat and standard arm pad.

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.