



RAZ-ART


Attendant Rotational Tilt
Mobile Shower Commode Chairs



Bathing and Toileting Products for Life

RAZ-ART FEATURES

The Raz-ART Attendant Rotational Tilt Mobile Shower Commode Chair utilizes a patented, dual-arc, orbital, tilt system that allows for the smallest footprint and lowest tilt effort of any chair on the market. A unique frame design provides two ranges of tool-less height adjustment and a 40° tilt range.



SPECIFICATIONS	RAZ-ART
WEIGHT CAPACITY	350 lbs
OVERALL WIDTH	22.5"
OVERALL LENGTH With head support and standard foot support in middle position	Upright: 35" / Tilted to 40°: 50.5"
SEAT WIDTH	18"
SEAT DEPTH (front of back post to front of seat)	16" - 18" or 19" - 21"
SEAT-TO-FLOOR HEIGHT – 1" (2.5 cm) increments – above standard moulded seat	Lower Range: 20" – 23" Higher Range: 24" – 27"
CLEARANCE HEIGHT UNDER SEAT – 1" (2.5 cm) increments – With pan brackets – Without pan brackets	Lower Range: 15" – 18" Higher Range: 19" – 22" Lower Range: 15.5" – 18.5" Higher Range: 19.5" – 22.5"
ARMREST HEIGHT – above standard moulded seat – 18"H Fixed Back – 15"H Adjusta-Back – 17"H Adjusta-Back – 20"H Adjusta-Back – 24"H Adjusta-Back	7" - 10" in ½" increments 9" 5.5" – 8.5" in 0.5" increments 6.5" – 9.5" in 0.5" increments 7.5" – 10.5" in 0.5" increments 8.5" – 11.5" in 0.5" increment
PRODUCT WEIGHT (with foot supports)	69 lbs
FRAME MATERIAL - base frame and seat frame	Powder-Coated Aluminum – Life Time Warranty Against Corrosion
FRAME MATERIAL – back frame and footrests	Stainless Steel – Life Time Warranty Against Corrosion
HARDWARE MATERIAL	Stainless Steel – Life Time Warranty Against Corrosion

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.