



RAZ-AT600


Attendant Tilt Heavy Duty
Mobile Shower Commode Chair



Bathing and Toileting Products for Life

RAZ-AT600 FEATURES

Raz-AT600 Attendant Tilt – 600 lb Heavy Duty Shower Commode Chair with a 25° tilt range. The specially engineered, Weight Assist Spring System (WASS) is available in two weight ranges and is adjustable to accommodate the user's weight and center of mass for tilt operation with minimal effort by the attendant.



SPECIFICATIONS	RAZ-AT600 (Heavy Duty Attendant Tilt)
WEIGHT CAPACITY	600 lbs
FRAME WIDTH (seat and back frame)	22"
OVERALL WIDTH (with Standard Armpad)	Seat frame width + 4.5"
OVERALL LENGTH (with Z709 Footrests)	38" - 39"
SEAT WIDTH (Commode Seat)	22" - 30"
SEAT DEPTH (front of back post to front of seat)	18"
CLEARANCE HEIGHT UNDER SEAT	18" – 19" with pan brackets 18.75" – 19.75" without pan brackets
ARMREST HEIGHT ABOVE SEAT	9.5"
SEAT TO FLOOR HEIGHT	21.75" - 22.75"
PRODUCT WEIGHT	5" casters: 82 lbs 6" casters: 84 lbs
FRAME AND HARDWARE MATERIAL	Stainless Steel – Life Time Warranty Against Corrosion

All measurements taken with standard molded seat and standard arm pad.

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.