



RAZ-AT

Attendant Tilt Mobile Shower Commode Chairs



Bathing and Toileting Products for Life

RAZ-AT FEATURES

The RAZ-AT (Attendant Tilt) tilts up to 40° to provide pressure reduction and maximum comfort. The unique frame design utilizes lockable gas springs that help lift the seat to an upright position, making the RAZ-AT one of the easiest tilting chairs to operate. The easy-to-clean, stainless steel frame fits over most toilets and is designed with unobstructed side and front access for the user or attendant.

EXTENSIVE RANGE OF OPTIONS AVAILABLE



SPECIFICATIONS	RAZ-AT (Attendant Tilt)
WEIGHT CAPACITY	400 lbs
FRAME WIDTH (seat and back frame)	Standard: 18"
OVERALL WIDTH (with Standard Armpad)	22.5"
OVERALL LENGTH (with footrest)	39"
SEAT WIDTH (Commode Seat)	18"
SEAT DEPTH (front of back post to front of seat)	16" or 19"
CLEARANCE HEIGHT UNDER SEAT	17" – 21" with pan brackets (in 1" increments) 17.75" – 21.75" without pan brackets (in 1" increments)
ARMREST HEIGHT ABOVE SEAT (20" Adjusta-Back)	7.5" - 10.5" in ½" increments
PRODUCT WEIGHT (with foot supports)	61 lbs
SEAT TO FLOOR HEIGHT	21.25" - 25.25" in 1" increments
FRAME AND HARDWARE MATERIAL	Stainless Steel – Life Time Warranty Against Corrosion

All measurements taken with standard molded seat and standard arm pad.

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com
T 1 (866) 650-6555 | F 1 (888) 966-6555