



## RAZ-CAT

Compact Attendant Tilt  
Mobile Shower Commode Chair



Bathing and Toileting Products for Life

## RAZ-CAT FEATURES

The Raz-CAT (Compact Attendant Tilt) has all of the features of our Raz-AT but in a more compact size. For pediatrics, adolescents or smaller adults, the Raz-CAT provides a fit that is better suited to physically smaller clients. To accommodate client growth, the Raz-CAT has adjustable seat depths and is available with 3 back-width options.

### EXTENSIVE RANGE OF OPTIONS AVAILABLE



SPECIFICATIONS	RAZ-CAT (Compact Attendant Tilt)
WEIGHT CAPACITY	300 lbs
FRAME WIDTH (seat and back frame)	18"
BACK FRAME WIDTHS	14", 16" or 18"
OVERALL WIDTH (with Standard Armpads)	21.75", 23" with 18" wide Back Frame
OVERALL LENGTH WHEN UPRIGHT	32"
SEAT WIDTH (Commode Seat)	18"
SEAT DEPTH (front of back post to front of seat)	11" – 16"
BACK HEIGHTS (top of seat to top of upholstery)	15" or 18"
RECLINING BACKS (top of seat to top of upholstery)	15" or 18"
CLEARANCE HEIGHT UNDER SEAT	16.75", 17.75", 18.75", 19.75", 20.75" with pan brackets 17.5", 18.5", 19.5", 20.5", 21.5" without pan brackets
ARMREST HEIGHT ABOVE SEAT	15" High Back: 5.5" – 8" in ½" increments 18" High Back: 6.5" – 9" in ½" increments
SEAT TO FLOOR HEIGHT	21.5" - 25.5" in 1" increments
PRODUCT WEIGHT (with foot supports)	55 lbs
FRAME AND HARDWARE MATERIAL	Stainless Steel – Life Time Warranty Against Corrosion

All measurements taken with standard molded seat and standard arm pad.

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com  
T 1 (866) 650-6555 | F 1 (888) 966-6555