



RAZ-SP

Self Propel Mobile Shower Commode Chair



Bathing and Toileting Products for Life

RAZ-SP FEATURES

The **RAZ-SP** Self Propel Mobile Shower Commode Chair comes standard with 4" of height adjustment, self-propelling rear wheels, 5" dual-locking casters, flip-up padded armrests, and commode pan. The easy-to-clean, stainless steel frame is designed with unobstructed side and front access for the user or attendant.

EXTENSIVE RANGE OF OPTIONS AVAILABLE



ALSO AVAILABLE

- Long frame — 3" longer for tall clients
- Low frame - 4" seat-to-floor height
- Reclining Back
- 16", 20" and 22" widths

| SPECIFICATIONS | RAZ-SP 22" WHEELS | RAZ-SP 24" WHEELS |
|--|---|---|
| WEIGHT CAPACITY | 400 lbs | 400 lbs |
| FRAME WIDTH (seat and back frame) | Standard: 18" | Standard: 18" |
| OVERALL WIDTH (with standard armpad) | 26.5" | 26.5" |
| OVERALL LENGTH (with footrest) | 39" | 39" |
| SEAT WIDTH (commode seat) | 18" | 18" |
| SEAT DEPTH (front of back post to front of seat) | 16" or 19" | 16" or 19" |
| CLEARANCE HEIGHT UNDER SEAT | 16", 17", 18", 20" without pan brackets 15.5", 16.5", 17.5", 19.5" with pan brackets | 15.75", 16.75", 17.75", 19.75" without pan brackets 15.25", 16.25", 17.25", 19.25" with pan brackets |
| ARMREST HEIGHT ABOVE SEAT (17" Adjusta-Back) | 7" - 10" in ½" increments | 7" - 10" in ½" increments |
| SEAT-TO-FLOOR HEIGHT | 20" - 24" in 1" increments | 20" - 24" in 1" increments |
| PRODUCT WEIGHT (with footrests) | 47 lbs | 49 lbs |
| FRAME AND HARDWARE MATERIAL | Stainless Steel (Lifetime warranty against corrosion) | Stainless Steel (Lifetime warranty against corrosion) |

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications, and options are subject to change without prior notice.