



RAZ-SP600

Heavy Duty Self Propel Mobile Shower Commode Chair



Bathing and Toileting Products for Life

RAZ-SP600 FEATURES

The **Raz-SP600** [Self-Propel, 600-lb capacity] rehab shower commode chair features the Multi-Position Back with fore/aft center-of-mass adjustment to accommodate a variety of body contours. Backrest and rear axle plate can be independently adjusted fore/aft to fine-tune center-of-mass position for optimized maneuverability

EXTENSIVE RANGE OF OPTIONS AVAILABLE



- ALSO AVAILABLE**
- Frame widths up to 36"
 - Custom frame depths
 - Custom Seats

SPECIFICATIONS	RAZ-SP600 (Bariatric Self-Propel)
WEIGHT CAPACITY	600 lbs
FRAME WIDTH (seat and back frame)	22" [Custom widths available]
OVERALL WIDTH (with standard arm pads)	30" with Standard Arm Pad
OVERALL LENGTH WHEN UPRIGHT (with Z709 footrests)	40" – 42"
SEAT WIDTH (commode seat)	23"
SEAT DEPTH (front of back post to front of seat)	18"
CLEARANCE HEIGHT UNDER SEAT	5" casters with pan brackets: 15.25" or 16.25" 5" casters without pan brackets: 16" or 17" 6" casters with pan brackets: 16.25", 17.25", 20.25" 6" casters without pan brackets: 17", 18", or 21"
ARMREST HEIGHT ABOVE SEAT	9.5"
SEAT-TO-FLOOR HEIGHT (5" casters)	20" and 21"
SEAT-TO-FLOOR HEIGHT (6" casters)	20", 21 and 24"
PRODUCT WEIGHT (with footrests)	5" casters: 68 lbs 6" casters: 69 lbs
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications, and options are subject to change without prior notice.



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com
T 1 (866) 650-6555 | F 1 (888) 966-6555