



RAZ-ZÜM

Attendant Propel Mobile Shower Commode Chair



The Züm-SP Self Propelled Folding Mobile Shower Commode Chair

ZÜM-AP AND ZÜM-SP FEATURES

The Züm-AP/SP are the first folding frame mobile shower commode chairs that offer the rigidity that you would expect from a premium product. The unique cross brace design allows the Züm to fit fully over a toilet for optimal alignment of the aperture with the toilet bowl.



SPECIFICATIONS	ZÜM-AP	ZÜM-SP
WEIGHT CAPACITY	350 lbs	350 lbs
OVERALL WIDTH – Frame Open	23.25"	25.25"
– Frame Folded	11"	13"
– Frame Folded	7.5" – (back assembly removed)	7.5" (back assembly and wheels removed)
OVERALL LENGTH	34"	40" with 22" wheels
– With standard footrests		41" with 24" wheels
– For SP, rear wheels in centre position		
SEAT WIDTH	18"	18"
CLEARANCE HEIGHT UNDER SEAT		
– With pan brackets	17.25"	17.25"
– Without pan brackets	17.75"	17.75"
ARMREST HEIGHT ABOVE SEAT	8" ¹	8" ¹
SEAT-TO-FLOOR HEIGHT – 1" increments	21.5" ¹	21.5" ¹
– Above standard molded seat		
PRODUCT WEIGHT	44 lbs	57 lbs
FRAME MATERIAL		
– Base frame and seat frame	Aluminum ²	Aluminum ²
– Back frame and footrests	Stainless Steel ²	Stainless Steel ²
– Hardware	Stainless Steel ¹	Stainless Steel ¹

¹ Measurements taken with standard molded seat and standard arm pads

² Lifetime warranty against corrosion

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.