



# TRAKZ TUB SLIDERS

Trakz-AT (Attendant Tilt) and Trakz-AP (Attendant Propel)

ACCESSIBLE BATHING MADE SIMPLE  
– NO RENOVATIONS NECESSARY



RAZ TRAKZ-AT  
Attendant Tilt

**TRAKZ is a tub slider system that enables individuals to shower over an existing bathtub without the need for expensive bathroom renovations.**

- Ideal for seniors and individuals with limited mobility.
- Promotes independence and accessibility in their own homes.
- Versatile and long-term solution for accessible bathing.
- Use the Trakz commode to reduce transfers and to allow for safe toileting and effortless integration onto the bridge for tub transfers.
- Trakz-AT uses a dual-arc tilt system, allowing for the lowest tilt effort of any chair on the market for positioning, ease of transfers and navigating over tub walls.

Bathing and Toileting Products for Life

# TRAKZ-AT FEATURES

**Patented Zbrilok™ Safety System** – Trakz is the only tub slider system with a patented safety feature that prevents the seat and user from moving onto the bridge unless the bridge is securely latched to both the tub base and the mobility base. This patented design, available only from Raz, can prevent serious injuries.



## Rotational Tilt – Dual-Arc Tilt Geometry (Trakz-AT only)

- 35-degree tilt range – for pressure redistribution, trunk positioning and raising client’s feet high enough to clear the tub wall
- Rotational tilt geometry reduces tilt effort and allows the Trakz to be mounted much closer to the rear wall

SPECIFICATIONS	TRAKZ-AP	TRAKZ-AT
WEIGHT CAPACITY		352 lbs
OVERALL WIDTH – with arm supports – with arm supports removed		23" 20.75"
OVERALL LENGTH with standard footrests	35.5"	36"
TILT RANGE	N/A	35°
BACK SUPPORT HEIGHT		20"
BRIDGE LENGTH	7", 8.5", 9.5", 10.5", 12.5", 14.5", 16.5" + custom lengths	
UNDERSIDE OF BRIDGE TO TOP OF SEAT	6.25"	8.5"
UNDERSIDE OF BRIDGE TO FLOOR – base configuration – high range – base configuration – low range		15.75" - 19.75" in 1" increments 11.75" - 15.75" in 1" increments
SEAT WIDTH		18"
SEAT DEPTH		17"
CLEARANCE HEIGHT UNDER SEAT – base configuration – high range – base configuration – low range * without brackets +0.75"		16.25" – 20.25" in 1" increments* 11.25" – 15.25" in 1" increments*
ARMREST HEIGHT ABOVE SEAT		10"
SEAT-TO-FLOOR HEIGHT – base configuration – high range – base configuration – low range	22" – 25" in 1" increments 17" – 21" in 1" increments	24.5" – 27.5" in 1" increments 19.5" – 23.5" in 1" increments
PRODUCT WEIGHT	77.5 lbs	95 lbs
FRAME MATERIALS – base frame – back frame and footrests – hardware		Aluminum <sup>1</sup> Stainless Steel <sup>1</sup> Stainless Steel <sup>1</sup>

<sup>1</sup> Lifetime warranty against corrosion

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.