



160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T.+1-866-650-6555 F. +1-888-966-6555
orders@motioncomposites.com
www.motioncomposites.com



AWARD 2010
HARDING 2013
INNOVATIVE 2014
PRODUCT 2016



5.1

VELOCE
ORDER FORM - CANADA

P.O. number: *	Date:
Company Name: *	Phone: *
Ordered by: *	
Clinician/Sales rep:	
Email:	
Marked for:	

Ship to:*	
Name:	
Address:	
City:	Prov:
Postal Code:	
Phone:	

* *Required*

MA080
English Manual (STD)

MA081
Manuel Francais

Taxes not included

MODEL											
1	WC20	Ultralight VELOCE Wheelchair (265 lb weight limit)						WM0001598 / WAMZ	\$4 495	☐	
2	WC09	HD Heavy Duty package (350 lb weight capacity)						WAMF	\$170	☐	
3	WC10	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber							\$30	☐	
SEAT WIDTH ³											
4	SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ² ☐			NC	☐
5	SW	Seat width:	19" ² ☐	20" ² ☐					WAMA	\$338	☐
SEAT DEPTH											
6	SD	Seat depth:	14" ¹ ☐	15" ¹ ☐	16" ☐	17" ☐	18" ☐			NC	☐
7	SD	Seat depth:	19" ☐	20" ☐					WAMB	\$255	☐


1- Depth achieved with shorter seat sling. 2- Available with HD option 3- Width calculator

FRAME LENGTH						
8	SP01	Standard frame length configuration ¹			STD	☐
9	SP02	Frame length 1" longer than standard configuration			NC	☐
10	SP03	Frame length 2" longer than standard configuration			NC	☐
11	SP04	Frame length 1" shorter than standard configuration			NC	☐
12	SP05	Frame length 2" shorter than standard configuration			NC	☐

1- Distance of 4.37" from seat end to bend

FRONT Seat-to-floor height

Caster sizes	Height available +/- 1/4"	Specify desired height: (inches)
3" Caster	14" to 21"	
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	



CASTERS										
13	CA01	Composite wheel w/ PU tire (1" width)	3"□	4"□	5"□	6"□			NC	□
14	CA10	Pneumatic casters 6" x 1 ¼"							\$160	□
15	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "						\$98	□
16	CA15		5" x 1 ½ "						\$98	□
17	CA16		6" x 1 ½ "						\$98	□
18	CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "						\$175	□
19	CA18		5" x 1 ½ "						\$185	□
20	CA19		6" x 1 ½ "						\$195	□
CASTER OPTIONS										
21	CA13	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)							\$489	□


Veloce - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height: _____
22" wheel ☐	13 1/2 to 18 1/4	
24" wheel ☐	14 1/2 to 19 1/4	
25" wheel ☐	15 to 19 3/4	
26" wheel ☐	15 1/2 to 20 1/4	

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

FRONT FRAME ANGLE


22	FA70	70° Front frame angle			NC	☐
23	FA75	75° Front frame angle			NC	☐
24	FA80	80° Front frame angle			NC	☐
25	FA85	85° Front frame angle			NC	☐
26	FA90	90° Front frame angle			NC	☐

SEAT TO FOOTREST HEIGHT

Set Length to: ¹


Length must be at least 1 3/4" shorter than front seat-to-floor height

1 - + or - 1/2 inches



Footrest width / Taper selection

Seat Width	Straight ¹	1" Taper
Outside Measurement	2 1/2" Narrower than Seat Width	4 1/2" Narrower than Seat Width
14"	11.5"	9.5"
15"	12.5"	10.5"
16"	13.5"	11.5"
17"	14.5"	12.5"
18"	15.5"	13.5"
19"	16.5"	14.5"
20"	17.5"	15.5"



1- May have interference with auto-folding and flip-up footplates

FOOTPLATES					
27	FP06	One piece adjustable angle flip-up (<i>Aluminum</i>)		WAMJ / SEFND0015	STD □
28	FP07	One piece adjustable angle flip-up (<i>Carbon fiber with rubber inserts</i>)		WAMJ / SEFND0015	\$450 □
29	FP28	One piece auto-folding adjustable angle (<i>Aluminum</i>)			\$360 □
30	FP29	One piece auto-folding adjustable angle (<i>Carbon fiber with rubber inserts</i>)			\$440 □
31	FP02	Newton adjustable angle footplate (<i>Composites - Two piece</i>)			\$145 □
32	FP03	Newton adjustable angle (<i>Aluminum - Two piece</i>)			\$165 □
33	FP04	Newton adjustable angle footplate (<i>Carbon fiber with rubber inserts - Two piece</i>)			\$380 □
FOOTPLATE ACCESSORIES					
34	FP08	Length adjustable heel loop (<i>not available with one piece footplate</i>)			NC □
35	FP09	Calf strap		SEFND0010	\$52 □
36	FP30	Bodypoint Aeromesh Padded calf strap			\$129 □

Back canes fixed height without push handle					
STRAIGHT BACK					
37	BC03	Ultralight fixed height carbon fiber back cane ¹			NC □
		9" □ 12" □ 15" □ 18" STD □ 21" □			
38	BC04	8° bend fixed height carbon fiber back cane ¹			\$155 □
		13 ½" □ 15" □ 18" STD □			

1 - Not compatible with push handles

Back canes, height adjustable, with option					
FIXED ANGLE BACK					
39	BC06	Straight carbon fiber back cane			NC □
		9 - 12" □ 12 - 15" □ 15 - 18" □ 18 - 21" □ STD			
40	BC08	8° bend Carbon fiber back cane			\$155 □
		13 ½ - 15" □ 15 - 18" □ 18 - 21" □ STD			
ADJUSTABLE ANGLE BACK - STRAIGHT CANES					
41	BC09	Angle adjustable Straight carbon fiber back cane			\$309 □
		12 - 15" □ 15 - 18" □ 18 - 21" □ STD			
42	BC10	Angle adjustable 8° bend carbon fiber back cane			\$360 □
		13 ½ - 15" □ 15 - 18" □ 18 - 21" □ STD			
43		Select back cane angle: 85° ¹ □ 90° □ 95° STD □ 100° □ 105° □ 110° □			
Set height to: ²					
		(measured from seat sling to top of back upholstery)			
PUSH HANDLE OPTION					
44	PH01	Push handles			STD □
45	PH05	Fold down push handles ^{3,5}			\$150 □
46	PH06	Clamp on height adjustable push handle ^{3, 4, 5} 9 - 20" STD □ 9-23" □			\$325 □
47	PH02	Without push handles			NC □

1 - Not available with a 8 degree bend cane 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit

4 - Not compatible with swing away armrests and 8° bend back canes

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY					
48	BH07	Slip-on Active-Mesh back upholstery			NC □
49	BH08	Adjustable tension Active-Mesh back upholstery		WAMD/WAME	\$245 □
50	BH99	Omit Back upholstery			STD □
RIGID BACK					
51		NXT back series	See pricelist p.7		□

SEAT CUSHION					
52	SC01	2" Cushion			\$175 □
53	SC02	3" Cushion			\$195 □
54	SC99	Omit seat cushion			STD □
55		NXT FIT cushion series	See pricelist p.7		□

Rear wheels												
ONE ARM DRIVE ¹		Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims										
56	HR22	One arm drive ^{2,3}	24"	left ☐	right ☐			WAML	\$1 395	☐		
57	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)							NC	☐		
REAR WHEELS												
58	RW01	Mag ⁴	22" ☐	24" ☐					\$129	☐		
59	RW02	Newton One - spoke wheel ⁵	22" ☐	24" ☐					NC	☐		
60	RW03	Newton Gravity - ultralight wheel	22" ☐	24" ☐	25" ☐	26" ☐			\$360	☐		
61	RW04	Spinerger Spox	22" ☐	24" ☐	25" ☐	26" ☐			\$800	☐		
62		Spinerger Spox spokes colors ⁶ Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐			NC	☐		
63	RW05	Spinerger LX	22" ☐	24" ☐	25" ☐	26" ☐			\$995	☐		
64		Spinerger LX spokes colors ⁶ Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐	NC	☐
65	RW07	Spinerger CLX Carbon Blade		24" ☐	25" ☐				\$1 540	☐		
66		Spinerger CLX spokes colors ⁶ Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐	NC	☐
TIRES						Available						
67	TI01	Soft Urethane 1 1/4 "				22" 24"				NC	☐	
68	TI02	Soft Urethane <i>knobby</i> 1 3/8 "				22", 24", 25", 26"				\$110	☐	
69	TI04	Soft Urethane - <i>Shox</i> 1 "				22", 24", 25", 26"				\$130	☐	
70	TI05	Pneumatic 1 3/8 "				22", 24", 25", 26"				NC	☐	
71	TI08	Pneumatic with Airless insert 1 3/8"				22", 24"				\$215	☐	
72	TI06	Pneumatic HP - 1"				22", 24", 25", 26"				\$110	☐	
73	TI07	Pneumatic HP - <i>Schwalbe Marathon Plus</i> 1"				24", 25", 26"				\$190	☐	

1- Offered with handrims on each side of the wheelchair

2- Only standard handrim position available

3- Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5 - Not compatible with Shox tires

6 - Colour spokes require an additional 10 days lead-time

Handrims										
HANDRIM										
74	HR01	Aluminum anodized							STD	☐
75	HR02	Plastic coated ³						WAMO	\$160	☐
76	HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}						WAMO	\$430	☐
77	HR07	Newton Air Grip ¹ <i>Black</i> ☐						WAMO	\$150	☐
78	HR08	Newton Air Grip ¹ <i>(08) Sunkissed Orange</i> ☐ <i>(09) Burgundy</i> ☐ <i>(10) Sapphire Blue</i> ☐ <i>(11) Monster Green</i> ☐ <i>(12) Charcoal</i> ☐ <i>(13) Fuchsia</i> ☐ <i>(14) Steel Blue</i> ☐ <i>(15) Ferrari Red</i> ☐ <i>(16) Acid Green</i> ☐ <i>(17) White</i> ☐ <i>(18) Sky Blue</i> ☐						WAMO	\$330	☐
79	HR28	Surge ^{1, 2, 3} <i>(1-1/2" Oval handrim w/ Gripton high friction strip)</i>							\$470	☐
80	HR05	Surge LT ^{1, 2, 3} <i>(1-1/8" Oval handrim w/ Gripton high friction strip)</i>							\$500	☐
81	HR04	Natural Fit ^{1, 2, 3}							\$310	☐
82	HR99	Omit handrim							NC	☐

1 - Not available with Mag wheels or OAD

2 - Not available in 20"

3- Only standard handrim position available


Handrims option									
Assembly position for aluminum and Newton AirGrip handrims									
83	HR29	Standard Handrim position <i>Full tab</i>						STD	☐
84	HR30	Narrow Handrim position <i>Short tab 2 positions</i>						\$25	☐
85	HR31	Super Narrow Handrim <i>Short tab 2 positions</i> STD ☐	<i>Short tab 1 position</i> ☐					\$25	☐
Thumb piece ⁴ selection Natural Fit handrims									
86	HR23	Standard grip						\$190	☐
87	HR27	Super grip						\$220	☐
88	HR099	Omit thumb piece						NC	☐

4- Compatible with Surge, Surge LT handrims

AXLE							
89	AX01	Quick release axle			WAMU	STD	☐
90	AX03	Quad release axle			WAMU	\$170	☐

REAR WHEEL CAMBER									
91	CM00	0° camber						STD	☐
92	CM03	3° camber						NC	☐
93	CM06	6° camber						\$80	☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable												Most responsive			
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"			
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK

94	WL08	Newton wheel lock - push to lock				STD	☐
95	WL02	Pull to lock				NC	☐
96	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			WAMR	\$175	☐
97	WL13 WL14	6" removable extension handle				\$67	☐
98	WL05	Scissor Lock				\$175	☐
99	WL99	Omit wheel locks				NC	☐

Armrests

FLIP BACK ARMREST (Convertible to "T" armrest) ¹			Height Adjustment		Removable function		
100	AR01	10" Desk Length	8 - 12" ☐	STD 10 - 14" ☐	Locked ☐	Unlocked ☐	STD \$380 ☐
101	AR02	14" Full Length	8 - 12" ☐	STD 10 - 14" ☐	Locked ☐	Unlocked ☐	STD \$380 ☐
"T" ARMREST ¹			Height Adjustment				
102	AR03	10" Desk Length	8 - 12" ☐	STD 10 - 14" ☐			\$360 ☐
103	AR04	14" Full Length	8 - 12" ☐	STD 10 - 14" ☐			\$360 ☐
SWING AWAY ARMREST ²			Height Adjustment				
104	AR05	Tubular arm pad	8 - 11"				\$295 ☐
105	AR06	10" Desk length arm pad	8 - 11"				\$295 ☐
GEL ARM PAD ³							
106	AR12	14" Full Length Arm Pad - Gel Ovation				SEAND0030	\$150 ☐
107	AR99	Omit armrest					NC ☐

1 - May have a 1/2" adjustment restriction on both limits of the range

2 - Could interfere with some rigid back rests

3 - Not available with swing away armrest

SIDE GUARDS (not available with Flip Back and "T" Armrest)						
108	SG01	Plastic side guards (removable)			WAMV	\$225 ☐
109	SG02	Carbon fiber side guards (fixed) ¹			WAMV	\$255 ☐
110	SG03	Carbon fiber side guards (removable)			WAMV	\$295 ☐
111	SG06	Plastic fender side guards (removable) ²			WAMV	\$320 ☐
112	SG04	Carbon fiber fender side guards (removable)			WAMV	\$495 ☐
113	SG99	Omit side guard				NC ☐

SIDE GUARDS OPTION

114	SG07	Side guards clip (pair) ³			New	NC ☐
-----	------	--------------------------------------	--	--	-----	-----------------------------

1- Not available with back upholstery

2- Not compatible with 20" and 22" rear wheels

3- Not available with fixed carbon fiber side guard

POSITIONING BELT						
115	SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐		\$45 ☐
116	SB02	Auto buckle	Standard 35" ☐	Long 60" ☐		\$55 ☐
117	SB07	Hardware attachment for 2 point seatbelt (1 pair)				NC ☐
118	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)				\$13 ☐
119	SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐	Long 61" ☐		\$82 ☐
120	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0030 \$170 ☐
121	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0020 \$195 ☐
122	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐	Long 24" ☐		\$324 ☐
123	A6SB99	Omit positioning belt				STD ☐

HEADREST ¹

124	AC01	Profiled Headrest	Small (3.5"x6") ☐	Medium (4.5"x7") ☐	Large (5.5"x8") ☐	\$205 ☐
125	AC28	Headrest mounting hardware with stabiliser bar ¹				\$338 ☐

1 - Only available with push handles

ACCESSORIES										
126	AC03	Anti tippers						\$115	☐	
127	AC27	Tip assist ¹						\$60	☐	
128	AC05	WC-19 Transit Tie Down w/Q'straint belt ²		Short belt 20-41" ☐	Long belt 20-46" ☐			\$345	☐	
129	AC06	Cane carrier		Left ☐	Right ☐			\$105	☐	
130	AC08	Spoke guards ³		Available: 22", 24", 26"			WAMP	\$145	☐	
131	AC10	Seat pouch						\$95	☐	
132	AC26	BodyPoint Mobility Bag						\$135	☐	
133	AC23	Newton Backpack						\$120	☐	
134	AC12	Velcro Closing Strap						\$50	☐	
135	AC13	Tool kit						\$50	☐	
136	AC16	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐			\$32	☐	
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	M/XXXXL ☐	

1 - Cannot be combined with Anti Tippers 2 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option
3 - Not available with Mag wheels and Spinergy

CARBON FIBER FINISH WITH DECALS									
137	FC	Highlight color:	(10) Ferrari Red	(11) White ☐	(12) Black ☐	(13) Sunkissed Orange ☐			NC ☐
			(14) Sky Blue ☐	(15) Fuchsia ☐	(16) Acid Green ☐				
TOTAL									

Notes :



BioFit™ CUSHION

	BioFit™	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	595	
	NXCABF1618	15"-16"	18"	595	
	NXCABF1620	15"-16"	20"	595	
	NXCABF1816	17"-18"	16"	595	
	NXCABF1818	17"-18"	18"	595	
	NXCABF1820	17"-18"	20"	595	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	595	
NXCABF2018	19" - 20"	18"	595	
NXCABF2020	19" - 20"	20"	625	
NXCABF2216	21" - 22"	16"	625	
NXCABF2218	21" - 22"	18"	625	
NXCABF2220	21" - 22"	20"	625	

Optima™ THORACIC BACK SUPPORT

	Optima™	Fits Width	Height	\$	✓
	NOACB1613	16"-18"	13"	690	
	NOACB1616	16"-18"	16"	690	
	NOACB1619	16"-18"	19"	690	
	NOACB1813	18"-20"	13"	690	
	NOACB1816	18"-20"	16"	690	
	NOACB1819	18"-20"	19"	690	
	NOACB2016	20"-22"	16"	690	
	NOACB2019	20"-22"	19"	690	
	NOACB2213	22"-24"	13"	720	
	NOACB2216	22"-24"	16"	720	
	NOACB2219	22"-24"	19"	720	

Optima™ DEEP THORACIC BACK SUPPORT

	Optima™ DEEP	Fits Width	Height	\$	✓
	NOPCB1616	16"-18"	16"	790	
	NOPCB1619	16"-18"	19"	790	
	NOPCB1816	18"-20"	16"	790	
	NOPCB1819	18"-20"	19"	790	
	NOPCB2016	20"-22"	16"	790	
	NOPCB2019	20"-22"	19"	790	
	NOPCB2213	22"-24"	13"	820	
	NOPCB2216	22"-24"	16"	820	
	NOPCB2219	22"-24"	19"	820	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAACB2216/20	22"-24"	16"-20"	720	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™ DEEP	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAPCB2216/20	22"-24"	16"-20"	720	

Classic™ PELVIC BACK SUPPORT

	Classic™	Fits Width	Height	\$	✓
	NACB1616	16"-18"	16"	690	
	NACB1618	16"-18"	18"	690	
	NACB1816	18"-20"	16"	690	
	NACB1818	18"-20"	18"	690	
	NACB1820	18"-20"	20"	690	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	720	
	NACB2218	22"-24"	18"	720	
	NACB2220	22"-24"	20"	720	

Classic™ DEEP PELVIC BACK SUPPORT

	Classic™ DEEP	Fits Width	Height	\$	✓
	NPCB1616	16"-18"	16"	690	
	NPCB1618	16"-18"	18"	690	
	NPCB1818	18"-20"	18"	690	
	NPCB1820	18"-20"	20"	690	
	NPCB2020	20"-22"	20"	690	
	NPCB2216	22"-24"	16"	720	
	NPCB2218	22"-24"	18"	720	
	NPCB2220	22"-24"	20"	720	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

		
Model #NIB-0784	\$	✓
10" x 5" Large Rectangular	195	

		
Model # NIB-0818	\$	✓
14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

		
Model #NIB-0747	\$	✓
Multi-Axis Hardware	145	

		
Model # NIB-0823	\$	✓
Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com