



VELOCE
 ORDER FORM - CANADA

1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

2- Dealer/Representative information

Dealer/Representative:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

* Required

MA080 ☐
 English Manual (STD)

MA081 ☐
 Manuel Français

MODEL			Taxes not included		
WC20	☐	Ultralight VELOCE Wheelchair (265 lb weight limit)		WM0001029 / WAMZ	\$4,945
WC09	☐	HD Heavy Duty package (350 lb weight capacity)		WAMF	\$180
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber			\$35
Frame					
SEAT WIDTH ³					
SW	☐	14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ² ☐			NC
	☐	19" ² ☐ 20" ² ☐		WAMA	\$360
SEAT DEPTH					
SD	☐	14" ¹ ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐			NC
	☐	19" ☐ 20" ☐		WAMB	\$275
FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame) ⁴					
SP01	☐	Standard frame length configuration ⁵			STD
SP02	☐	Frame length 1" longer than standard configuration			NC
SP03	☐	Frame length 2" longer than standard configuration			NC
SP04	☐	Frame length 1" shorter than standard configuration			NC
SP05	☐	Frame length 2" shorter than standard configuration			NC
TRANSIT TIE DOWN					
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Short belt 20"-41"			\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Long belt 20"-46"			\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁶			\$335

1- Only available with frame length 1" and 2" longer; 2- Available with HD option


3- Width calculator

4- See chart at page 7; 5- 4.375" from seat sling to front frame

6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

FRONT Seat-to-floor height

Caster sizes	Height available +/- 1/4"	Specify desired height:
3" Caster	14" to 21"	_____ (inches)
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	





CASTERS						
CA01-03	☐	Composite wheel w/ PU tire	3" x 1"			NC
CA01-04	☐		4" x 1"			NC
CA01-05	☐		5" x 1"			NC
CA01-06	☐		6" x 1"			NC
CA10	☐	Pneumatic casters	6" x 1 ¼"			\$170
CA14	☐	Newton UltraCaster Soft-Roll Composite	4" x 1 ½ "			\$100
CA15	☐		5" x 1 ½ "			\$100
CA16	☐		6" x 1 ½ "			\$100
CA17	☐	Newton UltraCaster Soft-Roll Aluminum	4" x 1 ½ "			\$185
CA18	☐		5" x 1 ½ "			\$190
CA19	☐		6" x 1 ½ "			\$205
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				\$510

Veloce - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height: (inches)
22" (501) wheel ☐	13 1/2" to 18 1/4"	
24" (540) wheel ☐	14 1/2" to 18 1/4"	
25" (559) wheel ☐	15" to 19 3/4"	
26" (590) wheel ☐	15 1/2" to 20 1/4"	


Configuration used: 1" seat slope with a center of gravity at 2"
 Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

FRONT FRAME ANGLE						
FA70	☐	70° Front frame angle				NC
FA75	☐	75° Front frame angle				NC
FA80	☐	80° Front frame angle				NC
FA85	☐	85° Front frame angle				NC
FA90	☐	90° Front frame angle				NC

SEAT TO FOOTREST LENGTH	
Set Length to: ¹	Length must be at least 1 ¼" shorter than front seat-to-floor height
1- +/- 1/2" 	

Footrest width / Taper selection		
Seat Width	Straight ¹	1" Taper
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width
14"	11.5"	9.5"
15"	12.5"	10.5"
16"	13.5"	11.5"
17"	14.5"	12.5"
18"	15.5"	13.5"
19"	16.5"	14.5"
20"	17.5"	15.5"

1- May have interference with auto-folding and flip-up footplates



FOOTPLATES								
FP02	☐	Newton adjustable angle footplate - Composites (<i>Two piece</i>)				\$155		
FP03	☐	Newton adjustable angle - Aluminum (<i>Two piece</i>)				\$175		
FP04	☐	Newton adjustable angle footplate - Carbon fiber with rubber inserts (<i>Two piece</i>)				\$395		
FP06	☐	One piece adjustable angle flip-up - Aluminum			WAMJ / SEFND0015	STD		
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts			WAMJ / SEFND0015	\$480		
FP28	☐	One piece auto-folding adjustable angle - Aluminum				\$380		
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				\$470		
FOOTPLATE ACCESSORIES								
FP08	☐	Length adjustable heel loop (<i>not available with one piece footplate</i>)				NC		
FP09	☐	Calf strap			SEFND0010	\$57		
FP30	☐	Bodypoint Aeromesh Padded calf strap				\$140		
Backrest								
FIXED HEIGHT WITHOUT PUSH HANDLE OPTIONS								
BC03-9		Ultralight carbon fiber back cane	9"	*Not available until Fall 2023		NC		
BC03-12		Ultralight carbon fiber back cane	12"	*Not available until Fall 2023		NC		
BC03-15		Ultralight carbon fiber back cane	15"	*Not available until Fall 2023		NC		
BC03-18		Ultralight carbon fiber back cane	18"	*Not available until Fall 2023		NC		
BC03-21		Ultralight carbon fiber back cane	21"	*Not available until Fall 2023		NC		
BC04-13		8° bend carbon fiber back cane	13 ½"	*Not available until Fall 2023		\$165		
BC04-15		8° bend carbon fiber back cane	15"	*Not available until Fall 2023		\$165		
BC04-18		8° bend carbon fiber back cane	18"	*Not available until Fall 2023		\$165		
HEIGHT ADJUSTABLE WITH PUSH HANDLES OPTIONS								
BC06-912		Straight Carbon fiber back cane	9" - 12"	*Not available until Fall 2023		NC		
BC06-1215		Straight Carbon fiber back cane	12" - 15"	*Not available until Fall 2023		NC		
BC06-1518		Straight Carbon fiber back cane	15" - 18"	*Not available until Fall 2023		NC		
BC06-1821		Straight Carbon fiber back cane	18" - 21"	*Not available until Fall 2023		NC		
BC08-1315		8° Bend Carbon fiber back cane	13.5" - 15"	*Not available until Fall 2023		\$165		
BC08-1518		8° Bend Carbon fiber back cane	15" - 18"	*Not available until Fall 2023		\$165		
BC08-1821		8° Bend Carbon fiber back cane	18" - 21"	*Not available until Fall 2023		\$165		
* Replacement options below: 1-Aluminum fixed height and angle, 2- Aluminum fixed height and angle adjustable or 3- Straight carbon fiber angle and height adjustable								
FIXED HEIGHT WITHOUT PUSH HANDLE OPTIONS								
BC30	☐	Straight aluminum back cane	9"			NC		
BC36	☐	Straight aluminum back cane	12"			NC		
BC42	☐	Straight aluminum back cane	15"			NC		
BC63	☐	8° bend aluminum back cane	13 ½"			\$165		
BC64	☐	8° bend aluminum back cane	15"			\$165		
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (<i>MUST select ONE push handle options</i>)								
BC11-912	☐	Straight Aluminum back cane	9" - 12"	increment of 1/2"	Set at ²	NC		
BC11-1215	☐	Straight Aluminum back cane	12" - 15"	increment of 1/2"	Set at ²	NC		
BC11-1518	☐	Straight Aluminum back cane	15" - 18"	increment of 1/2"	Set at ²	NC		
BC11-1821	☐	Straight Aluminum back cane	18" - 21"	increment of 1/2"	Set at ²	NC		
BC12-1315	☐	8° Bend Aluminum back cane	13.5" - 15"	increment of 1/2"	Set at ²	\$80		
BC12-1518	☐	8° Bend Aluminum back cane	15" - 18"	increment of 1/2"	Set at ²	\$80		
BC12-1821	☐	8° Bend Aluminum back cane	18" - 21"	increment of 1/2"	Set at ²	\$80		
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (<i>MUST select ONE push handle options</i>)								
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ²	\$330		
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²	\$330		
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²	\$330		
BC10-1315		8° bend Carbon fiber back cane	13.5" - 15"	*Not available until Fall 2023		\$380		
BC10-1518		8° bend Carbon fiber back cane	15" - 18"	*Not available until Fall 2023		\$380		
BC10-1821		8° bend Carbon fiber back cane	18" - 21"	*Not available until Fall 2023		\$380		
		Select back cane angle:	85° ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐
PUSH HANDLES OPTIONS (<i>ONE option at a time</i>)								
PH01	☐	Push handles				STD		
PH05	☐	Fold down push handles ^{3,5}				\$160		
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" ^{3, 4, 5}				\$345		
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" ^{3, 4, 5}				\$345		
PH02	☐	Omit push handles				NC		

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit
4- Not compatible with swing away armrests and 8° bend back canes; 5- Not compatible with headrest and attendant lock

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery			WAMD/WAME	\$260
BH99	☐	Omit Back upholstery				STD
	☐	NXT back Support	See pricelist			
SEAT CUSHION						
SC01	☐	2" Cushion				\$185
SC02	☐	3" Cushion				\$205
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	

Rear wheels						
ONE ARM DRIVE ¹				Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims, permanent axles, 0° camber		
HR22	☐	One arm drive - aluminum inner handrim ^{2,3}	24"	Left ☐	Right ☐	
HR60	☐	One arm drive - plastic coated inner handrim ^{2,3}	24"	Left ☐	Right ☐	
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS				Available		
RW01	☐	Mag ⁴		22" 24"		\$140
RW02	☐	Newton One - spoke wheel ⁵		22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel		22", 24", 25", 26"		\$380
RW04	☐	Spinergy Spox	Blk STD ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐	22", 24", 25", 26"		\$850
RW05	☐	Spinergy LX	Blk STD ☐ Pur ⁶ ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐ Pink ⁶ ☐ Grn ⁶ ☐ Org ⁶ ☐	22", 24", 25", 26"		\$1,055
RW07	☐	Spinergy CLX ⁷	Blk STD ☐ Pur ⁶ ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐ Pink ⁶ ☐ Grn ⁶ ☐ Org ⁶ ☐	24", 25"		\$1,063
TIRES				Available		
TI01	☐	Soft Urethane - 1 3/8"		22" 24"		NC
TI02	☐	Soft Urethane - knobby 1 3/8 "		22", 24", 25", 26"		\$120
TI04	☐	Soft Urethane - Shox 1"		22", 24", 25", 26"		\$140
TI05	☐	Pneumatic - 1 3/8"		22", 24", 25", 26"		NC
TI08	☐	Pneumatic with Airless insert - 1 3/8" ⁸		22", 24", 26"		\$225
TI06	☐	Pneumatic HP - 1"		22", 24", 25", 26"		\$120
TI07	☐	Pneumatic HP - Puncture resistant, black 1"		24", 25", 26"		\$200

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option

4- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 5- Not compatible with Shox tires;

6- Colour spokes require an additionnal 10 days lead-time; 7- Only available with 1" tires; 8- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims						
HANDRIM				Available		
HR01	☐	Aluminum anodized		22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²		22", 24", 25"		\$80
HR02	☐	Plastic coated ¹		22", 24", 25", 26"	WAMO	\$170
HR03	☐	High friction coated ^{1,2}		22", 24", 25", 26"	WAMO	\$460
HR07 to HR18	☐	Newton Air Grip ²	See page 6 for the color selection and price	22", 24", 25", 26"	WAMO	
HR28	☐	Surge ^{1,2} (1 1/2" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"		\$485
HR05	☐	Surge LT ^{1,2} (1 1/8" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"		\$515
HR04	☐	Natural Fit ^{1,2}		22", 24", 25", 26"		\$320
HR63	☐	Nova H ^{1,2} (1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"		\$515
HR64	☐	Optimum H ^{1,2} (1 1/8" heart-shaped handrim)		22", 24", 25", 26"		\$320
HR99	☐	Omit handrim				NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

Handrims option						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	Full tab			STD
HR030	☐	Narrow Handrim position no cut	Full tab			NC
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions			\$27
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐			\$27
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				\$200
HR27	☐	Super grip				\$230
HR099	☐	Omit thumb piece				NC


3- Only available with Aluminum and Air Grip handrims;

AXLE ¹							
AX01	☐	Quick release axle				WAMU	STD
AX03	☐	Quad release axle				WAMU	\$180

1- Not required with one arm drive

REAR WHEEL CAMBER							
CM00	☐	0° camber					STD
CM03	☐	3° camber					NC
CM06	☐	6° camber					\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
X =	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK							
WL08	☐	Newton wheel lock - push to lock					STD
WL16	☐	Aluminum push to lock with extension					\$73
WL02	☐	Pull to lock					NC
WL04	☐	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)				WAMR	\$185
WL13	☐	6" removable extension handle					\$73
WL14	☐	Scissors Lock					\$185
WL15	☐	Unilateral wheel lock - pull to lock ¹		Left ☐	Right ☐	WAMN	\$335
WL99	☐	Omit wheel locks					NC

1- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests							
HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ¹							
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		\$395
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		\$395
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		\$395
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹							
AR13-610	☐	10" Desk length armpad ⁴	6" - 10"				\$380
AR03-812	☐	10" Desk length armpad	8" - 12"				\$380
AR03-1014	☐	10" Desk length armpad	10" - 14"				\$380
AR14-610	☐	14" Full length armpad ⁴	6" - 10"				\$380
AR04-812	☐	14" Full length armpad	8" - 12"				\$380
AR04-1014	☐	14" Full length armpad	10" - 14"				\$380
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²							
AR05	☐	Tubular armrest - Short pad	8" - 11"				\$315
AR23	☐	Tubular armrest - Long pad	8" - 11"				\$315
AR06	☐	10" Desk length flat armpad	8" - 11"				\$315
AR17	☐	14" Full length flat armpad	8" - 11"				\$315
ARMRESTS OPTIONS							
AR12	☐	14" Full Length armpad ³ - Gel Ovation				SEAND0030	\$160
AR99	☐	Omit armrest					NC

1 - May have a 1/2" adjustment restriction on both limits of the range

2 - Could interfere with some rigid backrests

3 - Not available with swing away and flip back armrest

4 - May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)						
SG01	☐	Plastic side guards (removable)			WAMV	\$235
SG02	☐	Carbon fiber side guards (fixed)			WAMV	\$275
SG03	☐	Carbon fiber side guards (removable)			WAMV	\$315
SG06	☐	Plastic fender side guards - straight bracket (removable)			WAMV	\$340
SG08	☐	Plastic fender side guards - offset bracket (removable) ¹			WAMV	\$340
SG04	☐	Carbon fiber fender side guards (removable)			WAMV	\$510
SG99	☐	Omit side guard				NC
SIDE GUARDS OPTION						
SG07	☐	Side guards clip (pair) ²				NC

1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; 2- Not available with fixed carbon fiber side guard

POSITIONING BELT						
SB01-S	☐	Velcro adjustable - Standard 35"				\$50
SB01-L	☐	Velcro adjustable - Long 60"				\$50
SB02-S	☐	Auto buckle - Standard 35"				\$60
SB02-L	☐	Auto buckle - Long 60"				\$60
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"				\$90
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"				\$90
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"		SERND0030		\$180
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"		SERND0030		\$180
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"		SERND0030		\$180
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"		SERND0020		\$205
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"		SERND0020		\$205
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"		SERND0020		\$205
SB03	☐	Hardware attachment for 2 extra points (1 pair)				\$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"				\$345
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"				\$345
SB99	☐	Omit positioning belt				STD

HEADREST ¹						
AC01-S	☐	Profiled Headrest - Small (3.5"X6")				\$215
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")				\$215
AC01-L	☐	Profiled Headrest - Large (5.5"X8")				\$215
AC28	☐	Headrest mounting hardware with stabilizer bar 1				\$360

1 - Only available with push handles

ACCESSORIES						
AC03	☐	Anti-tippers				\$120
AC29	☐	Swing-away anti-tipper ¹ Left ☐ Right ☐				\$180
AC27	☐	Tip assist ²				\$65
AC06	☐	Cane and crutch holder Left ☐ Right ☐				\$115
AC08	☐	Spoke guards ³ Available: 22", 24", 26"		WAMP		\$155
AC10	☐	Seat pouch				\$100
AC26	☐	BodyPoint Mobility Bag				\$145
AC23	☐	Newton Backpack				\$130
AC30	☐	Newton multi-purpose pouch				\$75
AC12	☐	Velcro Closing Strap				\$55
AC13	☐	Tool kit				\$55
AC16	☐	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐				\$36

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinergy

Colors																				
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																				
AirGrip	☐ \$160	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
																				
Highlight ²	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color Specify PPG Paintit color:	additional 10 days lead-time																	\$422

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.;
Finish between Airgrip and frame might be different;

2- To see different highlight, visit our website: https://www.motioncomposites.com/en_ca/products/folding-wheelchairs/veloce

Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD	N/A	15"	16"	17"	18"	19"	20"
1	15"	16"	17"	18"	19"	20"	N/A
2	16"	17"	18"	19"	20"	N/A	N/A
-1	N/A	N/A	15"	16"	17"	18"	19"
-2	N/A	N/A	N/A	15"	16"	17"	18"

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4

Notes :

BioFit™ CUSHION		WIDTH									
DEPTH		12"	13"	14"	15"	16"	17"	18"	19"	20"	
	12"	☐ \$625	☐ \$625	☐ \$625							
	14"	☐ \$625	☐ \$625	☐ \$625							
	15"		☐ \$625	☐ \$625	☐ \$625						
	16"		☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	
	17"				☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	
	18"				☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	
	19"				☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	
	20"				☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	☐ \$655	
	Option: Change cover for incontinent / infection control outer cover										☐ \$N/C

SEATING AND CUSHION ACCESSORIES

	Model #	Description	Fits Back Width	\$	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	

THORACIC SHAPE			Optima™ THORACIC BACK SUPPORT				CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)			
			Fits Width	Height	Model # CARBON FIBER	\$ CARBON	Model # Aluminum	\$ Aluminum		
	13"-15"	7"			NOACB1307	☐ \$725				
	13"-15"	9"			NOACB1309	☐ \$725				
	13"-15"	13"			NOACB1313	☐ \$725				
	13"-15"	16"			NOACB1316	☐ \$725				
	14"-16"	7"			NOACB1407	☐ \$725				
	14"-16"	9"	NOCFACB1409	☐ \$1150	NOACB1409	☐ \$725				
	14"-16"	13"	NOCFACB1413	☐ \$1150	NOACB1413	☐ \$725				
	14"-16"	16"	NOCFACB1416	☐ \$1150	NOACB1416	☐ \$725				
	15"-17"	7"			NOACB1507	☐ \$725				
	15"-17"	9"	NOCFACB1509	☐ \$1150	NOACB1509	☐ \$725				
	15"-17"	13"	NOCFACB1513	☐ \$1150	NOACB1513	☐ \$725				
	15"-17"	16"	NOCFACB1516	☐ \$1150	NOACB1516	☐ \$725				
	16"-18"	7"			NOACB1607	☐ \$725				
	16"-18"	9"	NOCFACB1609	☐ \$1150	NOACB1609	☐ \$725				
	16"-18"	13"	NOCFACB1613	☐ \$1150	NOACB1613	☐ \$725				
	16"-18"	16"	NOCFACB1616	☐ \$1150	NOACB1616	☐ \$725				
	17"-19"	7"			NOACB1707	☐ \$725				
	17"-19"	9"	NOCFACB1709	☐ \$1150	NOACB1709	☐ \$725				
	17"-19"	13"	NOCFACB1713	☐ \$1150	NOACB1713	☐ \$725				
	17"-19"	16"	NOCFACB1716	☐ \$1150	NOACB1716	☐ \$725				
	18"-20"	7"			NOACB1807	☐ \$725				
	18"-20"	9"	NOCFACB1809	☐ \$1150	NOACB1809	☐ \$725				
	18"-20"	13"	NOCFACB1813	☐ \$1150	NOACB1813	☐ \$725				
	18"-20"	16"	NOCFACB1816	☐ \$1150	NOACB1816	☐ \$725				
	20"-22"	7"			NOACB2007	☐ \$725				
	20"-22"	9"	NOCFACB2009	☐ \$1150	NOACB2009	☐ \$725				
	20"-22"	13"	NOCFACB2013	☐ \$1150	NOACB2013	☐ \$725				
	20"-22"	16"	NOCFACB2016	☐ \$1150	NOACB2016	☐ \$725				

		Optima™ DEEP THORACIC BACK SUPPORT				CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)			
		Fits Width	Height	Model # CARBON FIBER	\$ CARBON	Model # Aluminum	\$ Aluminum		
13"-15"	7"			NOPCB1307	☐ \$825				
13"-15"	9"			NOPCB1309	☐ \$825				
13"-15"	13"			NOPCB1313	☐ \$825				
13"-15"	16"			NOPCB1316	☐ \$825				
14"-16"	7"			NOPCB1407	☐ \$825				
14"-16"	9"	NOCFPCB1409	☐ \$1250	NOPCB1409	☐ \$825				
14"-16"	13"	NOCFPCB1413	☐ \$1250	NOPCB1413	☐ \$825				
14"-16"	16"	NOCFPCB1416	☐ \$1250	NOPCB1416	☐ \$825				
15"-17"	7"			NOPCB1507	☐ \$825				
15"-17"	9"	NOCFPCB1509	☐ \$1250	NOPCB1509	☐ \$825				
15"-17"	13"	NOCFPCB1513	☐ \$1250	NOPCB1513	☐ \$825				
15"-17"	16"	NOCFPCB1516	☐ \$1250	NOPCB1516	☐ \$825				
16"-18"	7"			NOPCB1607	☐ \$825				
16"-18"	9"	NOCFPCB1609	☐ \$1250	NOPCB1609	☐ \$825				
16"-18"	13"	NOCFPCB1613	☐ \$1250	NOPCB1613	☐ \$825				
16"-18"	16"	NOCFPCB1616	☐ \$1250	NOPCB1616	☐ \$825				
17"-19"	7"			NOPCB1707	☐ \$825				
17"-19"	9"	NOCFPCB1709	☐ \$1250	NOPCB1709	☐ \$825				
17"-19"	13"	NOCFPCB1713	☐ \$1250	NOPCB1713	☐ \$825				
17"-19"	16"	NOCFPCB1716	☐ \$1250	NOPCB1716	☐ \$825				
18"-20"	7"			NOPCB1807	☐ \$825				
18"-20"	9"	NOCFPCB1809	☐ \$1250	NOPCB1809	☐ \$825				
18"-20"	13"	NOCFPCB1813	☐ \$1250	NOPCB1813	☐ \$825				
18"-20"	16"	NOCFPCB1816	☐ \$1250	NOPCB1816	☐ \$825				
20"-22"	7"			NOPCB2007	☐ \$825				
20"-22"	9"	NOCFPCB2009	☐ \$1250	NOPCB2009	☐ \$825				
20"-22"	13"	NOCFPCB2013	☐ \$1250	NOPCB2013	☐ \$825				
20"-22"	16"	NOCFPCB2016	☐ \$1250	NOPCB2016	☐ \$825				

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT HARDWARE OPTIONS

All backs come STANDARD with QuickRelease Mounting Hardware. The following are fixed hardware options available at no charge when ordered with a back support.

Model #	Description	Weight of Hardware Only	✓
	QuickFit Quick Release Hardware System	STANDARD (2 LBS)	STANDARD
NIB-1642	Micro Mount 2 point fixed hardware. Back width needs to match chair width	Weight savings of 1.4 lb over STANDARD	☐
NIB-1207	2 point fixed hardware. Width expandable to fit 3 sizes of chair	Weight savings of 1.0 lb over STANDARD	☐

BACK SUPPORT ACCESSORIES – PRIVACY FLAP

	Model #	Description	Fits Back Width	\$	QTY
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	100	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	100	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	100	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	100	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com