


VELOCE
ORDER FORM - CANADA

1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

2- Dealer/Representative information

Dealer/Representative:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

* Required

MA080 ☐
English Manual (STD)

MA081 ☐
Manuel Français

MODEL				Taxes not included		
WC20	☐	Ultralight VELOCE Wheelchair (265 lb weight limit)		WAMZ	WM0001598	\$5,295
WC09	☐	HD Heavy Duty package (350 lb weight capacity)			WAMF	\$180
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber				\$35
Frame						
SEAT WIDTH ³						
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ² ☐
	☐	19" ² ☐	20" ² ☐	21" ⁷ ☐	22" ⁷ ☐	
SEAT DEPTH						
SD	☐	14" ¹ ☐	15" ☐	16" ☐	17" ☐	18" ☐
	☐	19" ☐	20" ☐			
FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame) ⁴						
SP01	☐	Standard frame length configuration ⁵				STD
SP02	☐	Frame length 1" longer than standard configuration				NC
SP03	☐	Frame length 2" longer than standard configuration				NC
SP04	☐	Frame length 1" shorter than standard configuration				NC
SP05	☐	Frame length 2" shorter than standard configuration				NC
TRANSIT TIE DOWN						
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Short belt 20"-41"				\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Long belt 20"-46"				\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁶				\$345

1- Only available with frame length 1" and 2" longer; 2- Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;

4- See chart on page 7; 5- 4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option;

7- Only available with HD option and the one piece adjustable angle flip-up


3. Width calculator


FRONT Seat-to-floor height

Caster sizes	Height available +/- 1/4"	Specify desired height: _____ (inches)
☐ 3" Caster	14" to 21"	
☐ 4" Caster	14.5" to 21"	
☐ 5" Caster	15.5" to 21"	
☐ 6" Caster	16.5" to 21" (N/A with 90° front frame)	

Frog Legs FRONT Seat-to-floor height

Caster Sizes	Height available +/- 1/4"	Specify desired height: _____ (inches)
☐ 3" Caster	15.25" to 21"	
☐ 4" Caster	15.75" to 21"	
☐ 5" Caster	16.75" to 21"	
☐ 6" Caster	17.75" to 21" (N/A with 90° front frame)	



CASTERS				Availability		
CA01-03 to CA01-06	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		NC
CA10	☐	Pneumatic casters (1 ¼")		6"		\$170
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$105
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		\$198
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks				\$525

Veloce - REAR Seat-to-floor Height

WHEEL SIZES		Height available	
☐	22" (501) wheel	13 1/2" to 18 1/4"	Specify desired height: (inches)
☐	24" (540) wheel	14 1/2" to 19 1/4"	
☐	25" (559) wheel	15" to 19 3/4"	
☐	26" (590) wheel	15 1/2" to 20 1/4"	

Configuration used: 1" seat slope with a center of gravity at 2"


Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

FRONT FRAME ANGLE					
FA70	☐	70° Front frame angle			NC
FA75	☐	75° Front frame angle			NC
FA80	☐	80° Front frame angle			NC
FA85	☐	85° Front frame angle			NC
FA90	☐	90° Front frame angle			NC


SEAT TO FOOTREST LENGTH

Set Length to: ¹

Length must be at least 1 ¾" shorter than front seat-to-floor height



Footrest width / Taper selection		
Seat Width	Straight ¹	1" Taper
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width
14"	11.5"	9.5"
15"	12.5"	10.5"
16"	13.5"	11.5"
17"	14.5"	12.5"
18"	15.5"	13.5"
19"	16.5"	14.5"
20"	17.5"	15.5"
21"	18.5"	16.5"
22"	19.5"	17.5"



1- Not compatible with auto-folding footplate

FOOTPLATES						
FP06	☐	One piece adjustable angle flip-up - Aluminum	WAMJ	SEFND0015	STD	
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts	WAMJ	SEFND0015	\$480	
FP28	☐	One piece auto-folding adjustable angle - Aluminum			\$380	
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts			\$470	
FOOTPLATE ACCESSORIES						
FP09	☐	Calf strap		SEFND0010	\$57	
FP30	☐	Bodypoint Aeromesh Padded calf strap			\$140	

Backrest									
HEIGHT ADJUSTABLE WITH PUSH HANDLES OPTIONS (MUST select ONE push handle options)									
BC06-912	☐	Straight Carbon fiber back cane	9" - 12"	increment of 1/2"	Set at ² :				NC
BC06-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ² :				NC
BC06-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :				NC
BC06-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :				NC
BC08-1315	☐	8° Bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ² :				\$175
BC08-1518	☐	8° Bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :				\$175
BC08-1821	☐	8° Bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :				\$175
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ² :				\$330
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :				\$330
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :				\$330
BC10-1315	☐	8° bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ² :				\$380
BC10-1518	☐	8° bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :				\$380
BC10-1821	☐	8° bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :				\$380
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
PUSH HANDLES OPTIONS (ONE option at a time)									
PH01	☐	Push handles							STD
PH05	☐	Fold down push handles ^{3,5}							\$165
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" ^{3, 4, 5}							\$345
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" ^{3, 4, 5}							\$345
PH02	☐	Omit push handles							NC

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit 4- Not compatible with swing away armrests and 8° bend back canes; 5- Not compatible with headrest and attendant lock

BACK UPHOLSTERY									
BH07	☐	Slip-on Active-Mesh back upholstery							NC
BH08	☐	Adjustable tension Active-Mesh back upholstery						WAMD/WAME	\$260
BH99	☐	Omit Back upholstery							STD
	☐	NXT back Support						See price list	
SEAT CUSHION									
SC01	☐	2" Cushion							\$185
SC02	☐	3" Cushion							\$205
SC99	☐	Omit seat cushion							STD
	☐	NXT cushion						See price list	


Rear wheels													
ONE ARM DRIVE ¹			Only available with Newton One rear wheels, permanent axles, 0° camber										
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐					WAML	\$1,480			
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐					WAML	\$1,595			
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)								NC			
REAR WHEELS								Available					
RW01	☐	Mag ³				22" 24"			\$140				
RW02	☐	Newton One - spoke wheel ⁴				22" 24"			SF				
RW03	☐	Newton Gravity - ultralight wheel				22", 24", 25", 26"			\$390				
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	22", 24", 25", 26"	\$875				
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	22", 24", 25", 26"	\$1,095
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	24", 25"	\$1,650
TIRES								Available					
AIRLESS													
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹				22" 24"			NC				
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{7,9}				22", 24", 25", 26"			\$120				
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>				22", 24", 25", 26"			\$140				
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸				22", 24", 26"			\$225				
PNEUMATIC													
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹				22", 24", 25", 26"			NC				
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>				22", 24", 25", 26"			\$120				
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>				24", 25", 26"			\$200				
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰				22", 24", 25"			\$130				
<i>1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available; 3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additionnal 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;</i>													
Handrims													
HANDRIM								Available					
HR01	☐	Aluminum anodized				22", 24", 25", 26"			STD				
HR61	☐	Aluminum black hard anodized ²				22", 24", 25"			\$85				
HR02	☐	Plastic coated ¹				22", 24", 25", 26"	WAMO		\$170				
HR03	☐	High friction coated ^{1,2,3}				22", 24", 25", 26"	WAMO		\$460				
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	<i>See page 6 for the color selection and price</i>			22", 24", 25", 26"	WAMO						
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"			\$485				
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"			\$515				
HR04	☐	Natural Fit ^{1,2,3,4}				22", 24", 25", 26"			\$320				
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)				22", 24", 25", 26"			\$485				
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"			\$515				
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)				22", 24", 25", 26"			\$320				
HR99	☐	Omit handrim							NC				
<i>1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive; 4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece</i>													
Handrims option													
ASSEMBLY POSITION													
HR29	☐	Standard Handrim position	Full tab						STD				
HR030	☐	Narrow Handrim position no cut	Full tab						NC				
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions						\$29				
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐	Short tab 1 position ☐					\$29				
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS													
HR23	☐	Standard grip							\$200				
HR27	☐	Super grip							\$230				
HR099	☐	Omit thumb piece							NC				

AXLE ¹						
AX01	☐	Quick release axle			WAMU	STD
AX03	☐	Quad release axle			WAMU	\$180

1- Not required with one arm drive

REAR WHEEL CAMBER						
CM00	☐	0° camber				STD
CM03	☐	3° camber				NC
CM06	☐	6° camber				\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"	
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK						
WL08	☐	Newton Wheel lock - push to lock			PRODUCT IMPROVED	STD
WL16	☐	Aluminum push to lock with extensions				\$73
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock (<i>works best with pneumatic or high pressure tire</i>)			WAMR	\$185
WL13	☐	6" removable extension handle ¹				\$73
WL05	☐	Scissors Lock				\$185
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐	WAMN	\$335
WL99	☐	Omit wheel locks				NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests						
HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ¹						
AR01-812	☐	10" Desk length armpad 8" - 12"	Locked ☐	Unlocked STD ☐		\$395
AR01-1014	☐	10" Desk length armpad 10" - 14"	Locked ☐	Unlocked STD ☐		\$395
AR02-812	☐	14" Full length armpad 8" - 12"	Locked ☐	Unlocked STD ☐		\$395
AR02-1014	☐	14" Full length armpad 10" - 14"	Locked ☐	Unlocked STD ☐		\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹						
AR13-610	☐	10" Desk length armpad ⁴ 6" - 10"				\$380
AR03-812	☐	10" Desk length armpad 8" - 12"				\$380
AR03-1014	☐	10" Desk length armpad 10" - 14"				\$380
AR14-610	☐	14" Full length armpad ⁴ 6" - 10"				\$380
AR04-812	☐	14" Full length armpad 8" - 12"				\$380
AR04-1014	☐	14" Full length armpad 10" - 14"				\$380
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²						
AR05	☐	Tubular armrest - Short pad 8" - 11"				\$315
AR23	☐	Tubular armrest - Long pad 8" - 11"				\$315
AR06	☐	10" Desk length flat armpad 8" - 11"				\$315
AR17	☐	14" Full length flat armpad 8" - 11"				\$315
ARMRESTS OPTIONS						
AR12	☐	14" Full Length armpad ³ - Gel Ovation			SEAND0030	\$160
AR99	☐	Omit armrest				NC

1- May have a 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with fixed carbon fiber side guards

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)						
SG01	☐	Plastic side guards (<i>removable</i>)	Straight bracket - STD	Offset bracket ¹	WAMV	\$235
SG02	☐	Carbon fiber side guards (<i>fixed</i>)			WAMV	\$285
SG03	☐	Carbon fiber side guards (<i>removable</i>)	Straight bracket - STD	Offset bracket ¹	WAMV	\$320
SG06	☐	Plastic fender side guards (<i>removable</i>)	Straight bracket - STD	Offset bracket ¹	WAMV	\$340
SG04	☐	Carbon fiber fender side guards (<i>removable</i>)	Straight bracket - STD	Offset bracket ¹	WAMV	\$510
SG99	☐	Omit side guard				NC

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1" ;

POSITIONING BELT						
SB01-S	☐	Velcro adjustable - Small 35" ¹				\$50
SB01-L	☐	Velcro adjustable - Long 60"				\$50
SB02-S	☐	Auto buckle - Small 35" ¹				\$60
SB02-M	☐	Auto buckle - Medium 48"				\$60
SB02-L	☐	Auto buckle - Long 60"				\$60
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"				\$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"				\$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"			SERND0030	\$185
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"			SERND0030	\$185
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"			SERND0030	\$185
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"			SERND0020	\$210
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"			SERND0020	\$210
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"			SERND0020	\$210
SB03	☐	Hardware attachment for 2 extra points (1 pair)				\$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"				\$355
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"				\$355
SB99	☐	Omit positioning belt				STD

1- not recommended for seat width of 18" and more

HEADREST ¹						
AC01-S	☐	Profiled Headrest - Small (3.5"X6")				\$215
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")				\$215
AC01-L	☐	Profiled Headrest - Large (5.5"X8")				\$215
AC28	☐	Headrest mounting hardware with stabilizer bar ¹				\$360

1- Only available with push handles

ACCESSORIES						
AC03	☐	Anti-tippers				\$120
AC29	☐	Swing-away anti-tipper ¹ Left ☐ Right ☐				\$185
AC27	☐	Tip assist ²				\$65
AC06	☐	Cane and crutch holder Left ☐ Right ☐				\$115
AC08	☐	Spoke guards ³ Available: 22", 24", 26"			WAMP	\$155
AC10	☐	Seat pouch				\$100
AC26	☐	BodyPoint Mobility Bag				\$145
AC23	☐	Newton Backpack				\$130
AC30	☐	Newton multi-purpose pouch				\$75
AC12	☐	Velcro Closing Strap				\$55
AC13	☐	Tool kit				\$55
AC16	☐	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐				\$36

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinergy

Colors																				
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																				
AirGrip	☐ \$160		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
																				
Highlight ²	☐		☐					☐		☐	☐			☐			☐			
FC99	☐	Custom Color	Specify PPG Paintit color: additional 10 days lead-time																	NC

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- To see different highlight, visit our website: https://www.motioncomposites.com/en_co/products/folding-wheelchairs/veloce

Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD	N/A						
1							N/A
2						N/A	N/A
-1	N/A	N/A					
-2	N/A	N/A	N/A				

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4

Notes :

FOLDING ON CHAIR – CDN ORDER FORM

EXTRA COVERS

Extra Outer smartx^{3D} Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	\$	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	__ __ X __ __	130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	__ __ X __ __	150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	__ __ X __ __	130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	__ __ X __ __	150	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com