


VELOCE
ORDER FORM - CANADA

1- Purchasing information

P.O. number:* _____ Date: _____
Company Name:* _____
Purchasing contact:* _____
Phone:* _____
Email:* _____
Contact preference:* Phone ☐ Email ☐
Marked for: _____

3- Shipping information

Ship to:* _____
Name: _____
Address: _____
City: _____ Province: _____
Postal Code: _____
Phone: _____

2- Dealer/Representative information

Dealer/Representative:* _____
Phone:* _____
Email:* _____
Contact preference:* Phone ☐ Email ☐

* Required

MA080 ☐

English Manual (STD)

MA081 ☐

Manuel Français

MODEL						Taxes not included		
WC20	☐	Ultralight VELOCE Wheelchair (265 lb weight limit)						\$5,295
WC09	☐	HD Heavy Duty package (350 lb weight capacity)						\$180
Frame								
SEAT WIDTH ³								
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ² ☐		NC
	☐	19" ² ☐	20" ² ☐	21" ⁷ ☐	22" ⁷ ☐			\$360
SEAT DEPTH								
SD	☐	14" ¹ ☐	15" ☐	16" ☐	17" ☐	18" ☐		NC
	☐	19" ☐	20" ☐					\$275
FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame) ⁴								
SP01	☐	Standard frame length configuration ⁵						STD
SP02	☐	Frame length 1" longer than standard configuration						NC
SP03	☐	Frame length 2" longer than standard configuration						NC
SP04	☐	Frame length 1" shorter than standard configuration						NC
SP05	☐	Frame length 2" shorter than standard configuration						NC
TRANSIT TIE DOWN								
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Short belt 20"-41"						\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Long belt 20"-46"						\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁶						\$345

1- Only available with frame length 1" and 2" longer; 2- Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;

4- See chart on page 7; 5- 4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option;

7- Only available with HD option and the one piece adjustable angle flip-up


3. Width calculator


FRONT Seat-to-floor height

Caster sizes		Height available +/- 1/4"	Specify desired height: _____ (inches)
☐	3" Caster	14" to 21"	
☐	4" Caster	14.5" to 21"	
☐	5" Caster	15.5" to 21"	
☐	6" Caster	16.5" to 21" (N/A with 90° front frame)	

Frog Legs FRONT Seat-to-floor height

Caster Sizes		Height available +/- 1/4"	Specify desired height: _____ (inches)
☐	3" Caster	15.25" to 21"	
☐	4" Caster	15.75" to 21"	
☐	5" Caster	16.75" to 21"	
☐	6" Caster	17.75" to 21" (N/A with 90° front frame)	





FRONT CASTERS						
CASTERS				Availability		
CA01-03 to CA01-06	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		NC
CA10	☐	Pneumatic casters (1 ¼")		6"		\$170
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$105
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		\$198
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks				\$525

Veloce - REAR Seat-to-floor Height

WHEEL SIZES		Height available	Specify desired height: (inches)
☐	22" (501) wheel	13 1/2" to 18 1/4"	
☐	24" (540) wheel	14 1/2" to 19 1/4"	
☐	25" (559) wheel	15" to 19 3/4"	
☐	26" (590) wheel	15 1/2" to 20 1/4"	


Configuration used: 1" seat slope with a center of gravity at 2"
 Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

FRONT FRAME ANGLE						
FA70	☐	70° Front frame angle				NC
FA75	☐	75° Front frame angle				NC
FA80	☐	80° Front frame angle				NC
FA85	☐	85° Front frame angle				NC
FA90	☐	90° Front frame angle				NC

SEAT TO FOOTREST LENGTH		
Set Length to: ¹	Length must be at least 1 ¾" shorter than front seat-to-floor height	
1- +/- 1/2"		
		


Footrest width / Taper selection		
Seat Width	Straight ¹	1" Taper
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width
14"	11.5"	9.5"
15"	12.5"	10.5"
16"	13.5"	11.5"
17"	14.5"	12.5"
18"	15.5"	13.5"
19"	16.5"	14.5"
20"	17.5"	15.5"
21"	18.5"	16.5"
22"	19.5"	17.5"

1- Not compatible with auto-folding footplate



REAR WHEELS						
ONE ARM DRIVE ¹			Only available with Newton One rear wheels, permanent axles, 0° camber			
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐		\$1,480
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐		\$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS					Available	
RW01	☐	Mag ³			22" 24"	\$140
RW02	☐	Newton One - spoke wheel ⁴			22" 24"	SF
RW03	☐	Newton Gravity - ultralight wheel			22", 24", 25", 26"	\$390
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐	22", 24", 25", 26"	\$875
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	22", 24", 25", 26"	\$1,095
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	24", 25"	\$1,650
TIRES					Available	
AIRLESS						
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹			22" 24"	NC
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{7,9}			22", 24", 25", 26"	\$120
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>			22", 24", 25", 26"	\$140
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸			22", 24", 26"	\$225
PNEUMATIC						
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹			22", 24", 25", 26"	NC
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			22", 24", 25", 26"	\$120
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"	\$200
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰			22", 24", 25"	\$130
1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available; 3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additionnal 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 25" and 26" are not available with Newton Gravity wheels 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;						
HANDRIMS						
HANDRIM					Available	
HR01	☐	Aluminum anodized			22", 24", 25", 26"	STD
HR61	☐	Aluminum black hard anodized ²			22", 24", 25"	\$85
HR02	☐	Plastic coated ¹			22", 24", 25", 26"	\$170
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"	\$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See page 6 for the color selection and price		22", 24", 25", 26"	
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	\$485
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	\$515
HR04	☐	Natural Fit ^{1,2,3,4}			22", 24", 25", 26"	\$320
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"	\$485
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"	\$515
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)			22", 24", 25", 26"	\$320
HR99	☐	Omit handrim				NC
1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive; 4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece						
HANDRIMS OPTION						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	Full tab			STD
HR030	☐	Narrow Handrim position no cut	Full tab			NC
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions			\$29
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐	Short tab 1 position ☐		\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				\$200
HR27	☐	Super grip				\$230
HR099	☐	Omit thumb piece				NC
5- Only available with Aluminum and Air Grip handrims;						
AXLE ¹						
AX01	☐	Quick release axle				STD
AX03	☐	Quad release axle				\$180
1- Not required with one arm drive						
REAR WHEEL CAMBER						
CM00	☐	0° camber				STD
CM03	☐	3° camber				NC
CM06	☐	6° camber				\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
X =	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK						
WL08	☐	Newton Wheel lock - push to lock				STD
WL16	☐	Aluminum push to lock with extensions				\$73
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock <i>(works best with pneumatic or high pressure tire)</i>				\$185
WL13 WL14	☐	6" removable extension handle ¹				\$73
WL05	☐	Scissors Lock				\$185
WL15	☐	Unilateral wheel lock - pull to lock ²		Left ☐	Right ☐	\$335
WL99	☐	Omit wheel locks				NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

https://www.motioncomposites.com/en-ca/support-and-education/how-to-documents

FOOTPLATES						
FP06	☐	One piece adjustable angle flip-up - Aluminum				STD
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				\$480
FP28	☐	One piece auto-folding adjustable angle - Aluminum				\$380
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				\$470
FOOTPLATE ACCESSORIES						
FP09	☐	Calf strap				\$57
FP30	☐	Bodypoint Aeromesh Padded calf strap				\$140

SEATING AND BACKREST

BACK UPHOLSTERY								
BH07	☐	Slip-on Active-Mesh back upholstery					NC	
BH08	☐	Adjustable tension Active-Mesh back upholstery					\$260	
BH99	☐	Omit Back upholstery					STD	
	☐	NXT back Support				See price list		
SEAT CUSHION								
SC01	☐	2" Cushion					\$185	
SC02	☐	3" Cushion					\$205	
SC99	☐	Omit seat cushion					STD	
	☐	NXT cushion				See price list		
HEIGHT ADJUSTABLE WITH PUSH HANDLES OPTIONS (MUST select ONE push handle options)								
BC06-912	☐	Straight Carbon fiber back cane	9" - 12"	increment of 1/2"	Set at ² :		NC	
BC06-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ² :		NC	
BC06-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		NC	
BC06-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		NC	
BC08-1315	☐	8° Bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ² :		\$175	
BC08-1518	☐	8° Bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		\$175	
BC08-1821	☐	8° Bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		\$175	
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)								
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ² :		\$330	
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		\$330	
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		\$330	
BC10-1315	☐	8° bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ² :		\$380	
BC10-1518	☐	8° bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ² :		\$380	
BC10-1821	☐	8° bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ² :		\$380	
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐
PUSH HANDLES OPTIONS (ONE option at a time)								
PH01	☐	Push handles					STD	
PH05	☐	Fold down push handles ^{3,5}					\$165	
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" ^{3, 4, 5}					\$345	
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" ^{3, 4, 5}					\$345	
PH02	☐	Omit push handles					NC	

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit

4- Not compatible with swing away armrests and 8° bend back canes and 9-12" back canes; 5- Not compatible with headrest and attendant lock

ARMRESTS								
HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ¹								
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐			\$395
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐			\$395
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐			\$395
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐			\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹								
AR13-610	☐	10" Desk length armpad ⁴	6" - 10"					\$380
AR03-812	☐	10" Desk length armpad	8" - 12"					\$380
AR03-1014	☐	10" Desk length armpad	10" - 14"					\$380
AR14-610	☐	14" Full length armpad ⁴	6" - 10"					\$380
AR04-812	☐	14" Full length armpad	8" - 12"					\$380
AR04-1014	☐	14" Full length armpad	10" - 14"					\$380
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²								
AR05	☐	Tubular armrest - Short pad	8" - 11"					\$315
AR23	☐	Tubular armrest - Long pad	8" - 11"					\$315
AR06	☐	10" Desk length flat armpad	8" - 11"					\$315
AR17	☐	14" Full length flat armpad	8" - 11"					\$315
ARMRESTS OPTIONS								
AR12	☐	14" Full Length armpad ³ - <i>Gel Ovation</i>						\$160
AR99	☐	Omit armrest						NC

1- May have a 1/2 " adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with fixed carbon fiber side guards

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES								
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>								
SG01	☐	Plastic side guards <i>(removable)</i>		Straight bracket - STD	Offset bracket ¹			\$235
SG02	☐	Carbon fiber side guards <i>(fixed)</i>						\$285
SG03	☐	Carbon fiber side guards <i>(removable)</i>		Straight bracket - STD	Offset bracket ¹			\$320
SG06	☐	Plastic fender side guards <i>(removable)</i>		Straight bracket - STD	Offset bracket ¹			\$340
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>		Straight bracket - STD	Offset bracket ¹			\$510
SG99	☐	Omit side guard						NC

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1" ;

POSITIONING BELT								
SB01-S	☐	Velcro adjustable - Small 35" ¹						\$50
SB01-L	☐	Velcro adjustable - Long 60"						\$50
SB02-S	☐	Auto buckle - Small 35" ¹						\$60
SB02-M	☐	Auto buckle - Medium 48"						\$60
SB02-L	☐	Auto buckle - Long 60"						\$60
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"						\$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"						\$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"						\$185
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"						\$185
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"						\$185
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"						\$210
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"						\$210
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"						\$210
SB03	☐	Hardware attachment for 2 extra points (1 pair)						\$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"						\$355
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"						\$355
SB99	☐	Omit positioning belt						STD

1- not recommended for seat width of 18" and more

HEADREST ¹								
AC01-S	☐	Profiled Headrest - Small (3.5"X6")						\$215
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")						\$215
AC01-L	☐	Profiled Headrest - Large (5.5"X8")						\$215
AC28	☐	Headrest mounting hardware with stabilizer bar ¹						\$360

1- Only available with push handles

ACCESSORIES												
WC10	☐	E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>										\$35
AC03	☐	Anti-tippers										\$120
AC29	☐	Swing-away anti-tipper ¹		Left ☐	Right ☐						\$185	
AC27	☐	Tip assist ²										\$65
AC06	☐	Cane and crutch holder		Left ☐	Right ☐						\$115	
AC08	☐	Spoke guards ³		Available: 22", 24", 25", 26"								\$155
AC10	☐	Seat pouch										\$100
AC26	☐	BodyPoint Mobility Bag										\$145
AC23	☐	Newton Backpack										\$130
AC30	☐	Newton multi-purpose pouch										\$75
AC12	☐	Velcro Closing Strap										\$55
AC13	☐	Tool kit										\$55

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate; 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerdy; the 24" spoke guard will be installed on the 25" wheel

Colors																					
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	☐	\$160	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	\$350
Frame	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC
Highlight ²	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC
FC99	☐	Custom Color Specify PPG Paintit color:																		<i>additional 10 days lead-time</i>	\$422

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- To see different highlight, visit our website: https://www.motioncomposites.com/en_ca/products/folding-wheelchairs/veloce

Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD	N/A						
1							N/A
2						N/A	N/A
-1	N/A	N/A					
-2	N/A	N/A	N/A				

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4



Explore the
EVO Custom Catalog.

CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :




FOLDING ON CHAIR – CDN ORDER FORM

BioFit™ CUSHION		WIDTH											
 Use wheelchair configuration Or use this size: NXCBF _____ (W" x D")	DEPTH	12"	\$625	\$625	\$625		N/A	N/A	N/A	N/A	N/A	N/A	N/A
		13"						N/A	N/A	N/A	N/A	N/A	N/A
		14"	\$625	\$625	\$625								
		15"		\$625	\$625	\$625	\$625						
		16"		\$625	\$625	\$625	\$625	\$625	\$625	\$655	\$655	\$655	\$655
		17"	N/A			\$625	\$625	\$625	\$625	\$655	\$655		
		18"	N/A			\$625	\$625	\$625	\$625	\$655	\$655	\$655	\$655
		19"	N/A			\$625	\$625	\$625	\$625	\$655	\$655		
		20"	N/A			\$625	\$625	\$625	\$655	\$655	\$655	\$655	\$655
		Option: Change cover for incontinent/infection control outer cover											

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.

*Customization (red boxes) is available with an upcharge of \$160.00 and has a 4 to 6 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES

	MODEL #	Description	Fits Back Width	Retail	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$80	

THORACIC SHAPE			Optima™				THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)			
	WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM				
	14"	13"	NOCFACB1413	\$1,150	NOACB1413	\$725				
	14"	16"	NOCFACB1416	\$1,150	NOACB1416	\$725				
	15"	13"	NOCFACB1513	\$1,150	NOACB1513	\$725				
	15"	16"	NOCFACB1516	\$1,150	NOACB1516	\$725				
	16"	13"	NOCFACB1613	\$1,150	NOACB1613	\$725				
	16"	16"	NOCFACB1616	\$1,150	NOACB1616	\$725				
	16"	19"			NOACB1619	\$725				
	17"	13"	NOCFACB1713	\$1,150	NOACB1713	\$725				
	17"	16"	NOCFACB1716	\$1,150	NOACB1716	\$725				
	18"	13"	NOCFACB1813	\$1,150	NOACB1813	\$725				
	18"	16"	NOCFACB1816	\$1,150	NOACB1816	\$725				
	18"	19"			NOACB1819	\$725				
	19"	16"			NOACB1916	\$725				
	19"	19"			NOACB1919	\$725				
	20"	13"	NOCFACB2013	\$1,150	NOACB2013	\$725				
	20"	16"	NOCFACB2016	\$1,150	NOACB2016	\$725				
	20"	19"			NOACB2019	\$725				
	21"	16"			NOACB2116	\$755				
	21"	19"			NOACB2119	\$755				
	22"	13"			NOACB2213	\$755				
	22"	16"			NOACB2216	\$755				
	22"	19"			NOACB2219	\$755				
	ITEM		WIDTH X DEPTH		RETAIL	QTY				
	Replacement Cover & Foam		NOACC		\$130					

		Optima™ DEEP				THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)			
WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM				
14"	13"	NOCFPCB1413	\$1,250	NOPCB1413	\$825				
14"	16"	NOCFPCB1416	\$1,250	NOPCB1416	\$825				
15"	13"	NOCFPCB1513	\$1,250	NOPCB1513	\$825				
15"	16"	NOCFPCB1516	\$1,250	NOPCB1516	\$825				
16"	13"	NOCFPCB1613	\$1,250	NOPCB1613	\$825				
16"	16"	NOCFPCB1616	\$1,250	NOPCB1616	\$825				
16"	19"			NOPCB1619	\$825				
17"	13"	NOCFPCB1713	\$1,250	NOPCB1713	\$825				
17"	16"	NOCFPCB1716	\$1,250	NOPCB1716	\$825				
18"	13"	NOCFPCB1813	\$1,250	NOPCB1813	\$825				
18"	16"	NOCFPCB1816	\$1,250	NOPCB1816	\$825				
18"	19"			NOPCB1819	\$825				
19"	16"			NOPCB1916	\$825				
19"	19"			NOPCB1919	\$825				
20"	13"	NOCFPCB2013	\$1,250	NOPCB2013	\$825				
20"	16"	NOCFPCB2016	\$1,250	NOPCB2016	\$825				
20"	19"			NOPCB2019	\$825				
21"	16"			NOPCB2116	\$860				
21"	19"			NOPCB2119	\$860				
22"	13"			NOPCB2213	\$860				
22"	16"			NOPCB2216	\$860				
22"	19"			NOPCB2219	\$860				
ITEM		WIDTH X DEPTH		RETAIL	QTY				
Replacement Cover & Foam		NOPCC		\$130					

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

THORACIC SHAPE	xtend™	WIDTH	HEIGHT	RETAIL	✓
	NHAACB (L)1213/17	12"	13"-17"	\$725	
	NHAACB (L)1413/17	14"	13"-17"	\$725	
	NHAACB (L)1613/17	16"	13"-17"	\$725	
	NHAACB (L)1813/17	18"	13"-17"	\$755	
	NHAACB (L)2013/17	20"	13"-17"	\$755	
	NHAACB (L)2213/17	22"	13"-17"	\$755	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		N3LHAACC		\$130	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓
	NHAPCB1616/20	16"	16"-20"	\$725	
	NHAPCB1816/20	18"	16"-20"	\$725	
	NHAPCB2016/20	20"	16"-20"	\$725	
	NHAPCB2216/20	22"	16"-20"	\$755	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		N3HAPCC		\$130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

Classic™ PELVIC BACK SUPPORT

PELVIC SHAPE	Classic™	WIDTH	HEIGHT	RETAIL	✓
	NACB1616	16"	16"	\$725	
	NACB1618	16"	18"	\$725	
	NACB1816	18"	16"	\$725	
	NACB1818	18"	18"	\$725	
	NACB1820	18"	20"	\$725	
	NACB2020	20"	20"	\$725	
	NACB2216	22"	16"	\$755	
	NACB2218	22"	18"	\$755	
	NACB2220	22"	20"	\$755	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		NACC		\$130	

Classic™ DEEP PELVIC BACK SUPPORT

	Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓
	NPCB1616	16"	16"	\$830	
	NPCB1618	16"	18"	\$830	
	NPCB1818	18"	18"	\$830	
	NPCB1820	18"	20"	\$830	
	NPCB2020	18"	20"	\$830	
	NPCB2216	20"	16"	\$860	
	NPCB2218	22"	18"	\$860	
	NPCB2220	22"	20"	\$860	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		NPCC		\$130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$70	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

MODEL #	Description	RETAIL	✓
NIB-0784	10" x 5" Large Rectangular	\$205	
NIB-0818	14" x 5" Extra-Large Rectangular	\$205	

MultiFit™ HEAD SUPPORT – FIXTURES

MODEL	Description	RETAIL	✓
NIB-0747	Multi-Axis Hardware	\$150	
NIB-0823	Head Support Receiver Bracket	\$60	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	NXCCBF _____ v	\$130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	NXCCBF _____ v	\$150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF _____ v	\$130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF _____ v	\$150	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com