



**VELOCE**

ORDER FORM - INTERNATIONAL

**1- Purchasing information**

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_  
 Company Name:\* \_\_\_\_\_  
 Purchasing contact:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_  
 Contact preference:\* Phone  Email

**3- Shipping information**

Ship to:\* \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_  
 Zip Code: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**2- ATP/Clinician information**

ATP/Clinician:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_  
 Contact preference:\* Phone  Email  \* Required

- MA083  English Manual (STD)    MA084  German Manual    MA085  Italian Manual    MA086  Dutch Manual    MA087  Spanish Manual    MA088  French Manual    MA089  Portuguese Manual

MODEL							USD (Taxes not included)	
WC20	<input type="checkbox"/>	Ultralight VELOCE Wheelchair (120 kg - 265 lb weight limit)					\$5,295	
WC09	<input type="checkbox"/>	HD Heavy Duty package (159 kg - 350 lb weight capacity)					\$540	
<b>Frame</b>								
<b>SEAT WIDTH<sup>3</sup> - cm (inches)</b>								
SW	<input type="checkbox"/>	35.56 (14") <input type="checkbox"/>	38.10 (15") <input type="checkbox"/>	40.64 (16") <input type="checkbox"/>	43.18 (17") <input type="checkbox"/>	45.72 (18") <sup>2</sup> <input type="checkbox"/>	48.26 (19") <sup>2</sup> <input type="checkbox"/>	NC
	<input type="checkbox"/>	50.80 (20") <sup>2</sup> <input type="checkbox"/>					\$535	
<b>SEAT DEPTH - cm (inches)</b>								
SD	<input type="checkbox"/>	35.56 (14") <sup>1</sup> <input type="checkbox"/>	38.10 (15") <input type="checkbox"/>	40.64 (16") <input type="checkbox"/>	43.18 (17") <input type="checkbox"/>	45.72 (18") <input type="checkbox"/>	48.26 (19") <input type="checkbox"/>	NC
	<input type="checkbox"/>	50.80 (20") <input type="checkbox"/>					\$515	
<b>FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame)<sup>4</sup></b>								
SP01	<input type="checkbox"/>	Standard frame length configuration <sup>5</sup>					STD	
SP02	<input type="checkbox"/>	Frame length 2.54 cm (1") longer than standard configuration					NC	
SP03	<input type="checkbox"/>	Frame length 5.08 cm (2") longer than standard configuration					NC	
SP04	<input type="checkbox"/>	Frame length 2.54 cm (1") shorter than standard configuration					NC	
SP05	<input type="checkbox"/>	Frame length 5.08 cm (2") shorter than standard configuration					NC	
<b>TRANSIT TIE DOWN</b>								
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>6</sup>		Short belt 50-104 cm (20" - 41")		\$365		
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>6</sup>		Long belt 50-116 cm (20" - 46")		\$365		
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>6</sup>					\$335	

1- Only available with frame length 1" and 2" longer; 2- Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;  
 4- See chart on page 7; 5- 4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option;  
 7- Only available with HD option and the one piece adjustable angle flip-up

[3. Width calculator](#)




**FRONT Seat-to-floor height**

Caster sizes	Height available +- 1/4"	Specify desired height:     _____ (inches)
<input type="checkbox"/> 3" Caster	35.56cm to 53.34cm (14" to 21")	
<input type="checkbox"/> 4" Caster	36.83cm to 53.34cm (14.5" to 21")	
<input type="checkbox"/> 5" Caster	39.37cm to 53.34cm (15.5" to 21")	
<input type="checkbox"/> 6" Caster	41.91cm to 53.34cm (16.5" to 21") (N/A with 90° front frame)	

**Frog Legs FRONT Seat-to-floor height**

Caster Sizes	Height available +- 1/4"	Specify desired height:     _____ (inches)
<input type="checkbox"/> 3" Caster	38.7cm to 53.3cm (15.25" to 21")	
<input type="checkbox"/> 4" Caster	40.0cm to 53.3cm (15.75" to 21")	
<input type="checkbox"/> 5" Caster	42.54cm to 53.3cm (16.75" to 21")	
<input type="checkbox"/> 6" Caster	45.08cm to 53.3cm (17.75" to 21") (N/A with 90° front frame)	



## FRONT CASTERS


CASTERS			Availability		
CA01-03 to CA01-06	<input type="checkbox"/>	Composite wheel w/ PU tire (1")	3", 4", 5", 6"		NC
CA10	<input type="checkbox"/>	Pneumatic casters (1 ¼")	6"		\$175
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 ½ ")	4", 5", 6"		\$105
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 ½ ")	4", 5", 6"		\$210
CASTER OPTIONS					
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks			\$515


### REAR Seat-to-floor Height

WHEEL SIZES	Height available		Specify desired height:  <hr style="width: 100px; margin: 0 auto;"/> (inches)
<input type="checkbox"/> 22" (501) wheel	34.29cm to 46.35cm (13.5" to 18.25")		
<input type="checkbox"/> 24" (540) wheel	36.83cm to 48.89cm (14.5" to 19.25")		
<input type="checkbox"/> 25" (559) wheel	38.10cm to 50.16cm (15" to 19.75")		
<input type="checkbox"/> 26" (590) wheel	39.37cm to 51.43cm (15.5" to 20.25")		


Configuration used: 2.5 cm (1") seat slope with a center of gravity at 5cm (2")

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1.25 cm (½") of selected Seat to Floor height"

FRONT FRAME ANGLE					
FA70	<input type="checkbox"/>	70° Front frame angle			NC
FA75	<input type="checkbox"/>	75° Front frame angle			NC
FA80	<input type="checkbox"/>	80° Front frame angle			NC
FA85	<input type="checkbox"/>	85° Front frame angle			NC
FA90	<input type="checkbox"/>	90° Front frame angle			NC

SEAT TO FOOTREST LENGTH		
Set Length to <sup>1</sup> :	Length must be at least 4,5 cm (1 ¼") shorter than front seat-to-floor height	
1 - +/- 1.27 cm (1/2")		
		

Footrest width / Taper selection			
Seat Width	Straight <sup>1</sup>		1" Taper
Outside Measurement	6,4 cm Narrower than Seat Width	11,4 cm Narrower than Seat Width	
35.56	14"	29.21cm/11.5"	24.13cm/9.5"
38.10	15"	31.75cm/12.5"	26.67cm/10.5"
40.64	16"	34.29cm/13.5"	29.21cm/11.5"
43.18	17"	36.83cm/14.5"	31.75cm/12.5"
45.72	18"	39.37cm/15.5"	34.29cm/13.5"
48.26	19"	41.91cm/16.5"	36.83cm/14.5"
50.80	20"	44.45cm/17.5"	39.37cm/15.5"
53.34	21"	46.99cm/18.5"	41.91cm/16.5"
55.88	22"	49.53cm/19.5"	44.45cm/17.5"



<sup>1</sup> - Not compatible with auto-folding footplate

## REAR WHEELS

ONE ARM DRIVE <sup>1</sup>		Only available with Newton One rear wheels, permanent axles, 0° camber				
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$1,510
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$1,595
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS						
						Available
RW01	<input type="checkbox"/>	Mag <sup>3</sup>			22" 24"	\$57
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>4</sup>			22" 24"	NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel			22", 24", 25", 26"	\$350
RW04	<input type="checkbox"/>	Spinergy Spox	Blk STD <input type="checkbox"/>	Wht <sup>5</sup> <input type="checkbox"/> Blue <sup>5</sup> <input type="checkbox"/> Red <sup>5</sup> <input type="checkbox"/> Yel <sup>5</sup> <input type="checkbox"/>	22", 24", 25", 26"	\$895
RW05	<input type="checkbox"/>	Spinergy LX	Blk STD <input type="checkbox"/>	Pur <sup>5</sup> <input type="checkbox"/> Wht <sup>5</sup> <input type="checkbox"/> Blue <sup>5</sup> <input type="checkbox"/> Red <sup>5</sup> <input type="checkbox"/> Yel <sup>5</sup> <input type="checkbox"/> Pink <sup>5</sup> <input type="checkbox"/> Grn <sup>5</sup> <input type="checkbox"/> Org <sup>5</sup> <input type="checkbox"/>	22", 24", 25", 26"	\$1,125
RW07	<input type="checkbox"/>	Spinergy CLX <sup>6</sup>	Blk STD <input type="checkbox"/>	Pur <sup>5</sup> <input type="checkbox"/> Wht <sup>5</sup> <input type="checkbox"/> Blue <sup>5</sup> <input type="checkbox"/> Red <sup>5</sup> <input type="checkbox"/> Yel <sup>5</sup> <input type="checkbox"/> Pink <sup>5</sup> <input type="checkbox"/> Grn <sup>5</sup> <input type="checkbox"/> Org <sup>5</sup> <input type="checkbox"/>	24", 25"	\$1,600
TIRES						
AIRLESS						
TI01	<input type="checkbox"/>	Soft Urethane 3.49cm (1 3/8") - Low Tread, black			22" 24"	NC
TI02	<input type="checkbox"/>	Soft Urethane 3.49cm (1 3/8") - Medium Tread, black <sup>7</sup>			22", 24", 25", 26"	\$105
TI04	<input type="checkbox"/>	Soft Urethane 2.54cm (1") - No tread, black - <i>Shox</i>			22", 24", 25", 26"	\$125
TI08	<input type="checkbox"/>	Pneumatic Airless insert 3.49cm (1 3/8") - Medium Tread <sup>8</sup>			22", 24", 26"	\$90
PNEUMATIC						
TI05	<input type="checkbox"/>	Pneumatic 3.49cm (1 3/8") Medium Tread			22", 24", 25", 26"	NC
TI06	<input type="checkbox"/>	Pneumatic HP 2.54cm (1") - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			22", 24", 25", 26"	\$110
TI07	<input type="checkbox"/>	Pneumatic HP 2.54cm (1") - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"	\$145
TI03	<input type="checkbox"/>	Pneumatic 5.08cm (2") - Large Tread, black			22", 24", 25"	\$115

- 1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available;  
 3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time;  
 6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 25" and 26" are not available with Newton Gravity wheels  
 8 - Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels

## HANDRIMS

HANDRIM					Available	
HR01	<input type="checkbox"/>	Aluminum anodized			22", 24", 25", 26"	STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>			22", 24", 25"	\$80
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>			22", 24", 25", 26"	\$145
HR03	<input type="checkbox"/>	High friction coated <sup>1,2,3</sup>			22", 24", 25", 26"	\$380
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip <sup>2</sup>	<i>See page 6 for the color selection and price</i>		22", 24", 25", 26"	
HR28	<input type="checkbox"/>	Surge <sup>1,2,3,4</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	\$425
HR05	<input type="checkbox"/>	Surge LT <sup>1,2,3,4</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	\$465
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2,3,4</sup>			22", 24", 25", 26"	\$395
HR62	<input type="checkbox"/>	Simi H <sup>1,2,3</sup> (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"	\$395
HR63	<input type="checkbox"/>	Nova H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"	\$425
HR64	<input type="checkbox"/>	Optimum H <sup>1,2,3</sup> (1 1/8" heart-shaped handrim)			22", 24", 25", 26"	\$325
HR99	<input type="checkbox"/>	Omit handrim				NC

- 1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;  
 4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

## HANDRIMS OPTION


ASSEMBLY POSITION					
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab		STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab		NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>5</sup>	Short tab 2 positions		\$32
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>5</sup>	Short tab 2 positions STD <input type="checkbox"/>	Short tab 1 position <input type="checkbox"/>	\$32
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	<input type="checkbox"/>	Standard grip			\$95
HR27	<input type="checkbox"/>	Super grip			\$150
HR099	<input type="checkbox"/>	Omit thumb piece			NC

- 5- Only available with Aluminum and Air Grip handrims;

AXLE <sup>1</sup>					
AX01	<input type="checkbox"/>	Quick release axle			STD
AX03	<input type="checkbox"/>	Quad release axle			\$135

- 1- Not required with one arm drive

REAR WHEEL CAMBER					
CM00	<input type="checkbox"/>	0° camber			STD
CM03	<input type="checkbox"/>	3° camber			NC
CM06	<input type="checkbox"/>	6° camber			\$165

CENTER OF GRAVITY														
	<b>Choose position of rear wheel axle plate:</b> Measure distance X from front of back post to center of axle receiver.													
	X =	2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"	10,2 4"
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
If not indicated, default setting is at 3,8 cm/1.5"														
Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6,3 cm (2.5"). Possible conflict between a 26" rear wheel and a center of gravity over 6,3 cm (2.5"). Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG														

WHEEL LOCK						
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock				STD
WL16	<input type="checkbox"/>	Aluminum push to lock <b>with</b> extensions				\$83
WL02	<input type="checkbox"/>	Pull to lock				NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock <i>(works best with pneumatic or high pressure tire)</i>				\$165
WL13 WL14	<input type="checkbox"/>	6" removable extension handle <sup>1</sup>				\$83
WL05	<input type="checkbox"/>	Scissors Lock				\$155
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock <sup>2</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>		\$335
WL99	<input type="checkbox"/>	Omit wheel locks				NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilities and restriction, refer to the table at page 6 and the "How-to documents" page of our website:

[https://www.motioncomposites.com/en\\_ca/support-and-education/how-to-documents](https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents)

FOOTPLATES						
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum				NC
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				\$470
FP28	<input type="checkbox"/>	One piece auto-folding adjustable angle - Aluminum				\$310
FP29	<input type="checkbox"/>	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				\$380
FOOTPLATE ACCESSORIES						
FP09	<input type="checkbox"/>	Calf strap				\$43
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap				\$125

### SEATING AND BACKREST

BACK UPHOLSTERY						
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery				NC
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery				\$190
BH99	<input type="checkbox"/>	Omit Back upholstery				NC

SEAT CUSHION						
SC01	<input type="checkbox"/>	5.08 cm (2") Cushion				\$194
SC02	<input type="checkbox"/>	7.62 cm (3") Cushion				\$200
SC99	<input type="checkbox"/>	Omit seat cushion				STD

HEIGHT ADJUSTABLE WITH PUSH HANDLES <i>(MUST select ONE push handle options)</i>						
BC06-912	<input type="checkbox"/>	Straight Carbon fiber back cane	22.86 - 30.48cm (9-12")	increment of 1/2"	Set at <sup>2</sup> :	NC
BC06-1215	<input type="checkbox"/>	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	increment of 1/2"	Set at <sup>2</sup> :	NC
BC06-1518	<input type="checkbox"/>	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	increment of 1/2"	Set at <sup>2</sup> :	NC
BC06-1821	<input type="checkbox"/>	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	increment of 1/2"	Set at <sup>2</sup> :	NC
BC08-1315	<input type="checkbox"/>	8° Bend Carbon fiber back cane	34.29 - 38.10cm (13 1/2-15")	increment of 1/2"	Set at <sup>2</sup> :	\$60
BC08-1518	<input type="checkbox"/>	8° Bend Carbon fiber back cane	38.10 - 45.72 (15-18")	increment of 1/2"	Set at <sup>2</sup> :	\$60
BC08-1821	<input type="checkbox"/>	8° Bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	increment of 1/2"	Set at <sup>2</sup> :	\$60

ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>								
BC09-1215	<input type="checkbox"/>	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	increment of 1/2"	Set at <sup>2</sup> :	\$220		
BC09-1518	<input type="checkbox"/>	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	increment of 1/2"	Set at <sup>2</sup> :	\$220		
BC09-1821	<input type="checkbox"/>	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	increment of 1/2"	Set at <sup>2</sup> :	\$220		
BC10-1315	<input type="checkbox"/>	8° bend Carbon fiber back cane	34.29 - 38.10cm (13 1/2-15")	increment of 1/2"	Set at <sup>2</sup> :	\$275		
BC10-1518	<input type="checkbox"/>	8° bend Carbon fiber back cane	38.10 - 45.72 (15-18")	increment of 1/2"	Set at <sup>2</sup> :	\$275		
BC10-1821	<input type="checkbox"/>	8° bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	increment of 1/2"	Set at <sup>2</sup> :	\$275		
		Select back cane angle:	85° <sup>1</sup> <input type="checkbox"/>	90° <input type="checkbox"/>	95° <b>STD</b> <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>

PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>						
PH01	<input type="checkbox"/>	Push handles				NC
PH05	<input type="checkbox"/>	Fold down push handles <sup>3,5</sup>				\$150
PH06-920	<input type="checkbox"/>	Clamp on height adjustable push handles 22.86 - 50.80cm (9" - 20") <sup>3,4,5</sup>				\$330
PH06-923	<input type="checkbox"/>	Clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") <sup>3,4,5</sup>				\$330
PH02	<input type="checkbox"/>	Omit push handles				NC

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat slina to top of back upholstery; 3- 120 kg/265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes and 9-12" back canes; 5- Not compatible with headrest and attendant lock

## ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST (convertible to "T" armrest) <sup>1</sup>							
AR01-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		\$390
AR01-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		\$390
AR02-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		\$390
AR02-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		\$390
HEIGHT ADJUSTABLE "T" ARMREST <sup>1</sup>							
AR13-610	<input type="checkbox"/>	25.4cm (10") Desk length armpad <sup>4</sup>	15.25 - 25cm (6" - 10")				\$360
AR03-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")				\$360
AR03-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")				\$360
AR14-610	<input type="checkbox"/>	35.56cm (14") Full length armpad <sup>4</sup>	15.25 - 25cm (6" - 10")				\$360
AR04-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")				\$360
AR04-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")				\$360
HEIGHT ADJUSTABLE SWING AWAY ARMREST <sup>2</sup>							
AR05	<input type="checkbox"/>	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")				\$250
AR23	<input type="checkbox"/>	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")				\$250
AR06	<input type="checkbox"/>	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")				\$250
AR17	<input type="checkbox"/>	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")				\$250
ARMREST OPTIONS							
AR12	<input type="checkbox"/>	35.56cm (14") Full Length Armpad - Gel Ovation (pair) <sup>3</sup>					\$165
AR99	<input type="checkbox"/>	Omit armrest					NC

1- May have a 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid back rests; Not available with fixed carbon fiber side guards

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- May have some restriction with seat-to-floor height and rear wheel

## ACCESSORIES

SIDE GUARDS (not available with Flip Back and "T" Armrest)							
SG01	<input type="checkbox"/>	Plastic side guards (removable)	Straight bracket - STD	Offset bracket <sup>1</sup>			\$210
SG02	<input type="checkbox"/>	Carbon fiber side guards (fixed)					\$275
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)	Straight bracket - STD	Offset bracket <sup>1</sup>			\$315
SG06	<input type="checkbox"/>	Plastic fender side guards (removable)	Straight bracket - STD	Offset bracket <sup>1</sup>			\$295
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)	Straight bracket - STD	Offset bracket <sup>1</sup>			\$550
SG99	<input type="checkbox"/>	Omit side guard					NC

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1"

POSITIONING BELT							
SB01-S	<input type="checkbox"/>	Velcro adjustable - Small 88.9cm (35") <sup>1</sup>					\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 152.40cm (60")					\$50
SB02-S	<input type="checkbox"/>	Auto buckle - Small 88.9cm (35") <sup>1</sup>					\$95
SB02-M	<input type="checkbox"/>	Auto buckle - Medium 121.92cm (48")					\$95
SB02-L	<input type="checkbox"/>	Auto buckle - Long 152.40cm (60")					\$95
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")					\$95
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")					\$95
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")					\$185
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")					\$185
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")					\$185
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")					\$205
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")					\$205
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")					\$205
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)					\$17
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")					\$310
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")					\$310
SB99	<input type="checkbox"/>	Omit positioning belt					STD

1- Not recommended for seat width of 45.72cm (18") and more

HEADREST <sup>1</sup>							
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small 8.8x15.2cm (3.5"x6")					\$195
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium 11.4x17.8cm (4.5"x7")					\$195
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large 14x2.3cm(5.5"x8")					\$195
AC28	<input type="checkbox"/>	Headrest mounting hardware with tension bar					\$320

1- Only available with push handles

ACCESSORIES							
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber					\$35
AC03	<input type="checkbox"/>	Anti-tippers					\$155
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$165
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>					\$67
AC06	<input type="checkbox"/>	Cane and crutch holder	Left <input type="checkbox"/>	Right <input type="checkbox"/>			\$130
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup>	Available: 22", 24", 25", 26"				\$125
AC10	<input type="checkbox"/>	Seat pouch					\$83
AC30	<input type="checkbox"/>	Newton multi-purpose pouch					\$62
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag					\$135
AC23	<input type="checkbox"/>	Newton Backpack					\$130
AC12	<input type="checkbox"/>	Velcro Closing Strap					\$57
AC13	<input type="checkbox"/>	Tool kit					\$50

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); Not available with amputee axle plate 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerigy; the 24" spoke guard will be installed on the 25" wheel

Colors																					
Options <sup>1</sup>	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$265
Frame		STD	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
Highlight <sup>2</sup>	<input type="checkbox"/>		<input type="checkbox"/>					<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>			<input type="checkbox"/>				NC
FC99	<input type="checkbox"/>	Custom Color	<a href="#">Specify PPG Paintit color:</a>															<i>additional 10 days lead-time</i>	\$390		

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;  
 2- To see different highlight, visit our website: [https://www.motioncomposites.com/en\\_intl/products/folding-wheelchairs/veloce](https://www.motioncomposites.com/en_intl/products/folding-wheelchairs/veloce)

### Frame length chart

Frame length	Seat Depth						
	35.56 (14")	38.10 (15")	40.64 (16")	43.18 (17")	45.72 (18")	48.26 (19")	50.80 (20")
STD	N/A						
1							N/A
2						N/A	N/A
-1	N/A	N/A					
-2	N/A	N/A	N/A				

### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4



### CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The EVO Custom Request Form is the best place to start.

Notes :