


P.O. number: * _____ Date: _____
Company Name: * _____ Phone: * _____
Ordered by: * _____
Clinician/ATP/RTS: _____
Email: _____
Marked for: _____

Ship to: * _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

		MA082 English Manual (STD)	MA089 Spanish Manual	Taxes not included	
* Required				USD	
MODEL					
1	WC20	Ultralight VELOCE Wheelchair (265 lb weight limit)		K0005	\$4 795 ☐
2	WC09	HD Heavy Duty package (350 lb weight capacity)		WAMF	\$500 ☐
3	WC10	E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>			\$30 ☐
SEAT WIDTH ³					
4	SW	Seat width:	14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ² ☐ 19" ² ☐		NC ☐
5	SW	Seat width:	20" ² ☐	E2201	\$495 ☐
SEAT DEPTH					
6	SD	Seat depth:	14" ¹ ☐ 15" ¹ ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐		NC ☐
7	SD	Seat depth:	20" ☐	E2203	\$475 ☐


1 - Depth achieved with shorter seat sling. 2 - Available with HD option [3. Width calculator](#)

FRAME LENGTH					
8	SP01	Standard frame length configuration ¹		STD	☐
9	SP02	Frame length 1" longer than standard configuration		NC	☐
10	SP03	Frame length 2" longer than standard configuration		NC	☐
11	SP04	Frame length 1" shorter than standard configuration		NC	☐
12	SP05	Frame length 2" shorter than standard configuration		NC	☐

1- Distance from seat and front frame is 4.5"

FRONT Seat-to-floor height

Caster sizes	Height available +- 1/4"	
3" Caster	14" to 21"	Specify desired height: _____ (inches)
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	




CASTERS					
13	CA01	Composite wheel w/ PU tire (1" width)	3" ☐ 4" ☐ 5" ☐ 6" ☐		NC ☐
14	CA10	Pneumatic casters 6" x 1 1/4"		\$160	☐
15	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 1/2 "	\$80	☐
16	CA15		5" x 1 1/2 "	\$80	☐
17	CA16		6" x 1 1/2 "	\$80	☐
18	CA17		4" x 1 1/2 "	\$175	☐
19	CA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 1/2 "	\$185	☐
20	CA19		6" x 1 1/2 "	\$195	☐
CASTER OPTIONS					
21	CA13	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)		E1015	\$475 ☐


Veloce - REAR Seat-to-floor Height


WHEEL SIZES	Height available	Specify desired height: _____
22" (501) wheel ☐	13 1/2 to 18 1/4	
24" (540) wheel ☐	14 1/2 to 19 1/4	
25" (559) wheel ☐	15 to 19 3/4	
26" (590) wheel ☐	15 1/2 to 20 1/4	

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

FRONT FRAME ANGLE							
22	FA70	70° Front Frame Angle				NC	☐
23	FA75	75° Front frame angle				NC	☐
24	FA80	80° Front Frame Angle				NC	☐
25	FA85	85° Front frame angle				NC	☐
26	FA90	90° Front Frame Angle				NC	☐

SEAT TO FOOTREST HEIGHT	
Set Length to ¹ :	Length must be at least 1 ¾" shorter than front seat-to-floor height
1 - +/- 1/2"	
	

Footrest width / Taper selection			
Seat Width	Straight ¹	1" Taper	
Outside Measurement	2 1/2" Narrower than Seat Width	4 1/2" Narrower than Seat Width	
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	
17"	14.5"	12.5"	
18"	15.5"	13.5"	
19"	16.5"	14.5"	
20"	17.5"	15.5"	

1- May interfere with auto-folding and flip-up footplates

FOOTPLATES						
27	FP06	One piece adjustable angle flip-up <i>(Aluminum)</i>		K0040	NC	☐
28	FP07	One piece adjustable angle flip-up <i>(Carbon fiber with rubber inserts)</i>		K0040	\$435	☐
29	FP28	One piece auto-folding adjustable angle <i>(Aluminum)</i>		E0951	\$288	☐
30	FP29	One piece auto-folding adjustable angle <i>(Carbon fiber with rubber inserts)</i>		E0951	\$352	☐
31	FP02	Newton adjustable angle footplate <i>(Composites - Two piece)</i>		K0040	\$210	☐
32	FP03	Newton adjustable angle <i>(Aluminum - Two piece)</i>		K0040	\$230	☐
33	FP04	Newton adjustable angle footplate <i>(Carbon fiber with rubber inserts - Two piece)</i>		K0040	\$415	☐
FOOTPLATE ACCESSORIES						
34	FP08	Length adjustable heel loop <i>(not available with one piece footplate)</i>		E0951	\$40	☐
35	FP09	Calf strap		K0038	\$35	☐
36	FP30	Bodypoint Aeromesh Padded calf strap			\$115	☐

1 - Not compatible with push handles

FIXED ANGLE BACK									
39	BC06	Straight carbon fiber back cane						NC	☐
		9 - 12" ☐	12 - 15" ☐	15 - 18" ☐	18 - 21" ☐ STD				
40	BC08	8° bend Carbon fiber back cane						\$50	☐
		13 ½ - 15" ☐	15 - 18" ☐	18 - 21" ☐	STD				
ADJUSTABLE ANGLE BACK									
41	BC09	Angle adjustable Straight carbon fiber back cane						\$200	☐
		12 - 15" ☐	15 - 18" ☐	18 - 21" ☐	STD				
42	BC10	Angle adjustable 8° bend carbon fiber back cane						\$250	☐
		13 ½ - 15" ☐	15 - 18" ☐	18 - 21" ☐	STD				
43		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
Set height to: ²						(measured from seat sling to top of back upholstery)			
		PUSH HANDLE OPTION							
44	PH01	Push handles						NC	☐
45	PH05	Fold down push handles ^{3, 5}						\$120	☐
46	PH06	Clamp on height adjustable push handle ^{3, 4, 5}		9 - 20" STD ☐	9-23" ☐			\$300	☐
47	PH02	Without push handles						NC	☐

BACK UPHOLSTERY							
48	BH07	Slip-on Active-Mesh back upholstery				STD	☐
49	BH08	Adjustable tension Active-Mesh back upholstery			E2611	\$300	☐
50	BH99	Omit Back upholstery				NC	☐
RIGID BACK							
51		NXT back series		See pricelist p.7			☐

Rear wheels									
ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims									
56	HR22	One arm drive ^{2,3}	24"	left ☐	right ☐			\$1 395	☐
57	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)						NC	☐
REAR WHEELS									
58	RW01	Mag ⁴	22" ☐	24" ☐				\$50	☐
59	RW02	Newton One - spoke wheel ⁵	22" ☐	24" ☐				NC	☐
60	RW03	Newton Gravity - ultralight wheel	22" ☐	24" ☐	25" ☐	26" ☐		\$300	☐
61	RW04	Spinergy Spox	22" ☐	24" ☐	25" ☐	26" ☐		\$795	☐
62		Spinergy Spox spokes colors ⁶ Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐		NC	☐
63	RW05	Spinergy LX	22" ☐	24" ☐	25" ☐	26" ☐		\$995	☐
64		Spinergy LX spokes colors ⁶ Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
65	RW07	Spinergy CLX Carbon Blade		24" ☐	25" ☐			\$1 480	☐
66		Spinergy CLX spokes colors ⁶ Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
TIRES					Available				
67	TI01	Soft Urethane 1 1/4 "			22" 24"		E2218	NC	☐
68	TI02	Soft Urethane <i>knobby</i> 1 3/8 "			22", 24", 25", 26"		E2218	\$95	☐
69	TI04	Soft Urethane - <i>Shox</i> 1"			22", 24", 25", 26"		E2218	\$115	☐
70	TI05	Pneumatic 1 3/8 "			22", 24", 25", 26"		E2211/E2212	NC	☐
71	TI08	Pneumatic with Airless insert 1 3/8"			22", 24"		E2211/E2213	\$80	☐
72	TI06	Pneumatic HP - 1"			22", 24", 25", 26"		E2211/E2212	\$95	☐
73	TI07	Pneumatic HP - <i>Schwalbe Marathon Plus</i> 1"			24", 25", 26"		E2211/E2212	\$125	☐

1 - Offered with handrims on each side of the wheelchair

2 - Only standard handrim position available

3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5 - Not compatible with Shox tires

6 - Colour spokes require an additionnal 10 days lead-time

Handrims									
HANDRIM									
74	HR01	Aluminum anodized						STD	☐
75	HR02	Plastic coated ³						\$95	☐
76	HR03	Q-Grip High friction Neoprene coated ^{1,2,3}						\$350	☐
77	HR07	Newton Air Grip ¹ Black ☐						\$95	☐
78	HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐	(09) Burgundy ☐	(10) Sapphire Blue ☐	(11) Monster Green ☐	(12) Charcoal ☐	(13) Fuchsia ☐	(14) Steel Blue ☐	\$240 ☐
		(15) Ferrari Red ☐	(16) Acid Green ☐	(17) White ☐	(18) Sky Blue ☐				
79	HR28	Surge ^{1,2,3} (1-1/2" Oval handrim w/ Gripton high friction strip)						\$375	☐
80	HR05	Surge LT ^{1,2,3} (1-1/8" Oval handrim w/ Gripton high friction strip)						\$405	☐
81	HR04	Natural Fit ^{1,2,3}						\$310	☐
82	HR99	Omit handrim						NC	☐

1 - Not available with Mag wheels or OAD

2 - Not available in 20"


3 - Only standard handrim position available

Handrims option									
Assembly position for aluminum and Newton AirGrip handrims									
83	HR29	Standard Handrim position	Full tab					STD	☐
84	HR30	Narrow Handrim position	Short tab 2 positions					\$25	☐
85	HR31	Super Narrow Handrim	Short tab 2 positions STD ☐	Short tab 1 position ☐				\$25	☐
Thumb piece ⁴ selection for Natural Fit handrims									
86	HR23	Standard grip						\$85	☐
87	HR27	Super grip						\$140	☐
88	HR099	Omit thumb piece						NC	☐

4 - Compatible with Surge, Surge LT handrims

AXLE									
89	AX01	Quick release axle						STD	☐
90	AX03	Quad release axle					K0108	\$125	☐
REAR WHEEL CAMBER ANGLE									
91	CM00	0° camber						NC	☐
92	CM03	3° camber						NC	☐
93	CM06	6° camber						\$150	☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable												Most responsive			
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"			
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK

94	WL08	Newton Wheel lock - push to lock				STD	☐
95	WL02	Pull to lock				NC	☐
96	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			E0974	\$150	☐
97	WL13 WL14	6" removable extension handle			E0961	\$75	☐
98	WL05	Scissor Lock				\$140	☐
99	WL99	Omit wheel locks				NC	☐

Armrests

FLIP BACK ARMREST (Convertible to "T" armrest) ¹			Height Adjustment		Removable function		
100	AR01	10" Desk Length	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐	STD \$350 ☐
101	AR02	14" Full Length	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐	STD \$350 ☐
"T" ARMREST ¹			Height Adjustment				
102	AR03	10" Desk Length	8 - 12" STD ☐	10 - 14" ☐		E0973	\$325 ☐
103	AR04	14" Full Length	8 - 12" STD ☐	10 - 14" ☐		E0973	\$325 ☐
SWING AWAY ARMREST ²			Height Adjustment				
104	AR05	Tubular arm pad	8 - 11"			E0973	\$225 ☐
105	AR06	10" Desk length arm pad	8 - 11"			E0973	\$225 ☐
GEL ARM PAD ³							
106	AR12	14" Full Length Arm Pad - Gel Ovation				SEAND0030	\$150 ☐
107	AR99	Omit armrest					NC ☐

1 - May have a 1/2" adjustment restriction on both limits of the range

2 - Could interfere with some rigid back rests

3 - Not available with swing away armrest

SIDE GUARDS (not available with Flip Back and "T" Armrest)

108	SG01	Plastic side guards (removable)				\$195	☐
109	SG02	Carbon fiber side guards (fixed) ¹				\$250	☐
110	SG03	Carbon fiber side guards (removable)				\$288	☐
111	SG06	Plastic fender side guards (removable) ²			WAMV	\$270	☐
112	SG04	Carbon fiber fender side guards (removable)			WAMV	\$495	☐
113	SG99	Omit side guard				NC	☐
SIDE GUARDS OPTION							
114	SG07	Side guards clip (pair) ³			New	NC	☐

1 - Not available with back upholstery

2 - Not compatible with 20" and 22" rear wheels

3 - Not available with fixed carbon fiber side guard

POSITIONING BELT

115	SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐		E0978	\$45	☐
116	SB02	Auto buckle	Standard 35" ☐	Long 60" ☐		E0978	\$85	☐
117	SB07	Hardware attachment for 2 point seatbelt (1 pair)					NC	☐
118	SB03	Additional attachment hardware for 4 point seatbelt					\$15	☐
119	SB04	BodyPoint Hip Belt - 2 point not padded		Medium 56" ☐	Long 61" ☐		\$80	☐
120	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		\$170	☐
121	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		\$190	☐
122	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps		Medium 18" ☐	Long 24" ☐		\$285	☐
123	SB99	Omit positioning belt					STD	☐

HEADREST ¹

124	AC01	Profiled Headrest	Small (3.5"x6") ☐	Medium (4.5"x7") ☐	Large (5.5"x8") ☐		\$180	☐
125	AC28	Headrest mounting hardware with stabiliser bar ¹					\$295	☐

1 - Only available with push handles

ACCESSORIES											
126	AC03	Anti tippers						E0971	\$140	☐	
127	AC27	Tip assist ¹							\$60	☐	
128	AC05	WC-19 Transit Tie Down w/Q'straint belt ²		Short belt 20-41" ☐		Long belt 20-46" ☐			\$320	☐	
129	AC06	Cane carrier		Left ☐		Right ☐		E2207	\$195	☐	
130	AC08	Spoke guards ³		Available: 22", 24", 26"				K0065	\$115	☐	
131	AC10	Seat pouch							\$75	☐	
132	AC26	BodyPoint Mobility Bag							\$125	☐	
133	AC23	Newton Backpack							\$120	☐	
134	AC12	Velcro Closing Strap							\$50	☐	
135	AC13	Tool kit							\$45	☐	
136	AC16	Motion Composites Polo Shirt		W/S ☐	W/M ☐	W/L ☐				\$32	☐
				M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	M/XXXXL ☐	
137	AC15	Military Patches		Army ☐	Navy ☐	Air Force ☐	Marines ☐	Coast Guard ☐		\$90	☐
138	AC15	Combat Wounded Purple Heart Patch							\$90	☐	

1 - Cannot be combined with Anti Tippers 2 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

3 - Not available with Mag wheels and Spinergy

CARBON FIBER FINISH									
139	FC	Highlight color:		(10) Ferrari Red	(11) White ☐	(12) Black ☐	(13) Sunkissed Orange ☐		NC ☐
				(14) Sky Blue ☐	(15) Fuchsia ☐	(16) Acid Green ☐			

TOTAL	
-------	--

Notes :



BioFit™ CUSHION

 HCPCS CODE: E2607-E2608	BioFit™	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	425	
	NXCABF1618	15"-16"	18"	425	
	NXCABF1620	15"-16"	20"	425	
	NXCABF1816	17"-18"	16"	425	
	NXCABF1818	17"-18"	18"	425	
	NXCABF1820	17"-18"	20"	425	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	425	
NXCABF2018	19" - 20"	18"	425	
NXCABF2020	19" - 20"	20"	425	
NXCABF2216	21" - 22"	16"	455	
NXCABF2218	21" - 22"	18"	455	
NXCABF2220	21" - 22"	20"	455	

Optima™ THORACIC BACK SUPPORT

 HCPCS CODE: E2620-E2621	Optima™	Fits Width	Height	\$	✓
	NOACB1613	16"-18"	13"	660	
	NOACB1616	16"-18"	16"	660	
	NOACB1619	16"-18"	19"	660	
	NOACB1813	18"-20"	13"	660	
	NOACB1816	18"-20"	16"	660	
	NOACB1819	18"-20"	19"	660	
	NOACB2016	20"-22"	16"	660	
	NOACB2019	20"-22"	19"	660	
	NOACB2213	22"-24"	13"	710	
	NOACB2216	22"-24"	16"	710	
	NOACB2219	22"-24"	19"	710	

Optima™ DEEP THORACIC BACK SUPPORT

 HCPCS CODE: E2620-E2621	Optima™ DEEP	Fits Width	Height	\$	✓
	NOPCB1616	16"-18"	16"	660	
	NOPCB1619	16"-18"	19"	660	
	NOPCB1816	18"-20"	16"	660	
	NOPCB1819	18"-20"	19"	660	
	NOPCB2016	20"-22"	16"	660	
	NOPCB2019	20"-22"	19"	660	
	NOPCB2213	22"-24"	13"	710	
	NOPCB2216	22"-24"	16"	710	
	NOPCB2219	22"-24"	19"	710	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

 HCPCS CODE: E2620-E2621	Xtend™	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAACB2216/20	22"-24"	16"-20"	710	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

 HCPCS CODE: E2620-E2621	Xtend™ DEEP	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAPCB2216/20	22"-24"	16"-20"	710	

Classic™ PELVIC BACK SUPPORT

 HCPCS CODE: E2620-E2621	Classic™	Fits Width	Height	\$	✓
	NACB1616	16"-18"	16"	660	
	NACB1618	16"-18"	18"	660	
	NACB1816	18"-20"	16"	660	
	NACB1818	18"-20"	18"	660	
	NACB1820	18"-20"	20"	660	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	710	
	NACB2218	22"-24"	18"	710	
	NACB2220	22"-24"	20"	710	

Classic™ DEEP PELVIC BACK SUPPORT

 HCPCS CODE: E2620-E2621	Classic™ DEEP	Fits Width	Height	\$	✓
	NPCB1616	16"-18"	16"	660	
	NPCB1618	16"-18"	18"	660	
	NPCB1818	18"-20"	18"	660	
	NPCB1820	18"-20"	20"	660	
	NPCB2020	20"-22"	20"	660	
	NPCB2216	22"-24"	16"	710	
	NPCB2218	22"-24"	18"	710	
	NPCB2220	22"-24"	20"	710	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

 HCPCS CODE: E0955	Model #NIB-0784			Model # NIB-0818		
	\$	✓		\$	✓	
	10" x 5" Large Rectangular	195		14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

 Model #NIB-0747	Model #NIB-0747			Model # NIB-0823		
	\$	✓		\$	✓	
	Multi-Axis Hardware	145		Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com