



1- Purchasing information

P.O. number:* _____ Date: _____
 Company Name:* _____
 Purchasing contact:* _____
 Phone:* _____
 Email:* _____
 Contact preference:* Phone ☐ Email ☐
 Marked for: _____

3- Shipping information

Ship to:* _____
 Name: _____
 Address: _____
 City: _____ State: _____
 Zip Code: _____
 Phone: _____

2- ATP/Clinician information

ATP/Clinician:* _____
 Phone:* _____
 Email:* _____
 Contact preference:* Phone ☐ Email ☐

* Required

MA082 ☐
 English Manual (STD)

MA089 ☐
 Spanish Manual

Taxes not included

MODEL		USD	
WC25	Veloce wheelchair - Helio Base (265 lb weight limit)	K0005	\$2,395
WC23	Carbon Fiber Frame		\$2,295
WC24	Fixed Front End		\$1,295
WC09	HD Heavy Duty package (350 lb weight capacity) ⁸		\$550
Frame			
SEAT WIDTH ³			
SW	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ² ☐ 19" ²		NC
	☐ 20" ² ☐ 21" ⁷ ☐ 22" ⁷	E2201	\$525
SEAT DEPTH			
SD	☐ 14" ¹ ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"		NC
	☐ 20"	E2203	\$575
FRAME LENGTH (Seat sling depth plus distance from seat sling to front of frame) ⁴			
SP01	☐ Standard frame length configuration ⁵		STD
SP02	☐ Frame length 1" longer than standard configuration		NC
SP03	☐ Frame length 2" longer than standard configuration		NC
SP04	☐ Frame length 1" shorter than standard configuration		NC
SP05	☐ Frame length 2" shorter than standard configuration		NC
TRANSIT TIE DOWN			
AC05-S	☐ WC-19 Transit Tie Down w/Q'straint belt ⁶ Short belt 20"-41"		\$365
AC05-L	☐ WC-19 Transit Tie Down w/Q'straint belt ⁶ Long belt 20"-46"		\$365
AC35	☐ WC-19 Transit Tie Down without belt ⁶		\$350

1- Only available with frame length 1" and 2" longer; 2- Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;

4- See chart on page 7; 5- 4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option;

7- Only available with HD option and the one piece adjustable angle flip-up; 8-HD fee is required for clients weighing 300-350 lbs


[3. Width calculator](#)

FRONT Seat-to-floor height

Caster sizes	Height available +/- 1/4"	
☐ 3" Caster	14" to 21"	Specify desired height: _____ (inches)
☐ 4" Caster	14.5" to 21"	
☐ 5" Caster	15.5" to 21"	
☐ 6" Caster	16.5" to 21" (N/A with 90° front frame)	

Frog Legs FRONT Seat-to-floor height

Caster Sizes	Height available +/- 1/4"	
☐ 3" Caster	15.25" to 21"	Specify desired height: _____ (inches)
☐ 4" Caster	15.75" to 21"	
☐ 5" Caster	16.75" to 21"	
☐ 6" Caster	17.75" to 21" (N/A with 90° front frame)	



FRONT CASTERS						
CASTERS				Availability		
CA01-03 to CA01-06	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		NC
CA10	☐	Pneumatic casters (1 ¼")		6"		\$165
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$95
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		\$215
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks			E1015	\$490

REAR Seat-to-floor Height

WHEEL SIZES		Height available	Specify desired height: (inches)
☐	22" (501) wheel	13 1/2" to 18 1/4"	
☐	24" (540) wheel	14 1/2" to 19 1/4"	
☐	25" (559) wheel	15" to 19 3/4"	
☐	26" (590) wheel	15 1/2" to 20 1/4"	

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

REAR WHEELS									
ONE ARM DRIVE ¹				Only available with Newton One rear wheels, permanent axles, 0° camber					
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐			E0958	\$1,495	
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐			E0958	\$1,600	
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)						NC	
REAR WHEELS					Available				
RW01	☐	Mag ³			22" 24"			\$57	
RW02	☐	Newton One - spoke wheel ⁴			22" 24"			NC	
RW03	☐	Newton Gravity - ultralight wheel			22", 24", 25", 26"			\$360	
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐	22", 24", 25", 26"			\$915	
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	22", 24", 25", 26"			\$1,145	
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	24", 25"			\$1,625	
TIRES (The brand and color may vary depending on parts availability)					Available				
AIRLESS									
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹			22" 24"		E2218	NC	
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{7,9}			22", 24", 25", 26"		E2218	\$97	
TI13	☐	Soft Urethane 1" - No tread, black - <i>Sentinel</i>			20", 22", 24", 25", 26"		E2218	\$118	
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸			22", 24", 26"		E2211/E2213	\$95	
PNEUMATIC									
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹			22", 24", 26"		E2211/E2212	NC	
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			22", 24", 25", 26"		E2211/E2212	\$105	
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"		E2211/E2212	\$150	
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰			22", 24", 25"		E2211/E2212	\$115	

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available;

3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additionnal 10 days lead-time;




6- Only available with 1" tires; 7- 22" is only available with Mags, Newton One rear wheels and Newton Gravity wheels

8 - Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels; 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;


HANDRIMS						
HANDRIM				Available		
HR01	☐	Aluminum anodized		22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²		22", 24", 25"		\$75
HR02	☐	Plastic coated ¹		22", 24", 25", 26"		\$110
HR03	☐	High friction coated ^{1,2,3}		22", 24", 25", 26"		\$375
HR07-HR18 and HR65- HR71	☐	Newton Air Grip ²	See page 6 for the color selection and price	22", 24", 25", 26"		
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"		\$450
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"		\$495
HR04	☐	Natural Fit ^{1,2,3,4}		22", 24", 25", 26"		\$395
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)		22", 24", 25", 26"		\$395
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"		\$445
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)		22", 24", 25", 26"		\$345
HR99	☐	Omit handrim				NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

4 - 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

HANDRIMS OPTION																																																			
ASSEMBLY POSITION																																																			
HR29	☐	Standard Handrim position	Full tab						STD																																										
HR030	☐	Narrow Handrim position no cut ⁶	Full tab						NC																																										
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions						\$30																																										
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐	Short tab 1 position ⁶ ☐					\$30																																										
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS																																																			
HR23	☐	Standard grip							\$95																																										
HR27	☐	Super grip							\$150																																										
HR099	☐	Omit thumb piece							NC																																										
5- Only available with Aluminum and Air Grip handrims; 6- Not available with the combination of 24" Newton Gravity, Spinerigy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;																																																			
AXLE ¹																																																			
AX01	☐	Quick release axle							STD																																										
AX03	☐	Quad release axle						K0108	\$130																																										
1- Not required with one arm drive																																																			
REAR WHEEL CAMBER ANGLE																																																			
CM00	☐	0° camber							STD																																										
CM03	☐	3° camber							NC																																										
CM06	☐	6° camber							\$165																																										
CENTER OF GRAVITY																																																			
		Choose position of rear wheel axle plate: Measure distance X from front of back post to center of axle receiver. <table border="1"> <thead> <tr> <th colspan="14">Most stable</th> </tr> </thead> <tbody> <tr> <td>X =</td> <td>1"</td> <td>1/4"</td> <td>1/2"</td> <td>3/4"</td> <td>2"</td> <td>1/4"</td> <td>1/2"</td> <td>3/4"</td> <td>3"</td> <td>1/4"</td> <td>1/2"</td> <td>3/4"</td> <td>4"</td> </tr> <tr> <td></td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> <td>☐</td> </tr> </tbody> </table> If not indicated, default setting is at 1.5" Possible conflict between a seat depth under 16" and a center of gravity over 2.5". Possible conflict between a 26" rear wheel and a center of gravity over 2.5". Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG								Most stable														X =	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Most stable																																																			
X =	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"																																						
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐																																						
WHEEL LOCK																																																			
WL08	☐	Newton Wheel lock - push to lock							STD																																										
WL02	☐	Pull to lock							NC																																										
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)						E0974	\$155																																										
WL13	☐	6" removable extension handle ¹						E0961	\$77																																										
WL14	☐																																																		
WL05	☐	Scissors Lock							\$155																																										
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐					\$330																																										
WL99	☐	Omit wheel locks							NC																																										
1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)																																																			
2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website: https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents																																																			
FRONT FRAME ANGLE																																																			
FA70	☐	70° Front frame angle								NC																																									
FA75	☐	75° Front frame angle								NC																																									
FA80	☐	80° Front frame angle								NC																																									
FA85	☐	85° Front frame angle								NC																																									
FA90	☐	90° Front frame angle								NC																																									
SEAT TO FOOTREST LENGTH																																																			
Set Length to ¹:		Length must be at least 1 1/4" shorter than front seat-to-floor height																																																	
1 - +/- 1/2"																																																			
																																																			

Footrest width / Taper selection			
Seat Width	Straight ¹		1" Taper
Outside Measurement	2 ½" Narrower than Seat Width		4 ½" Narrower than Seat Width
14"	11.5"	☐	9.5" ☐
15"	12.5"	☐	10.5" ☐
16"	13.5"	☐	11.5" ☐
17"	14.5"	☐	12.5" ☐
18"	15.5"	☐	13.5" ☐
19"	16.5"	☐	14.5" ☐
20"	17.5"	☐	15.5" ☐
21"	18.5"	☐	16.5" ☐
22"	19.5"	☐	17.5" ☐



1- Not compatible with auto-folding footplates

FOOTPLATES						
FP06	☐	One piece adjustable angle flip-up - Aluminum			K0040	NC
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				\$465
FP28	☐	One piece auto-folding adjustable angle - Aluminum			K0040	\$295
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				\$375
FOOTPLATE ACCESSORIES						
FP09	☐	Calf strap			K0038	\$39
FP30	☐	Bodypoint Aeromesh Padded calf strap			K0038	\$118

SEATING AND BACKREST						
BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				STD
BH08	☐	Adjustable tension Active-Mesh back upholstery			E2611/E2612	\$305
BH99	☐	Omit Back upholstery				NC
	☐	NXT back Support			See pricelist	
SEAT CUSHION						
SC01	☐	2" Cushion				\$170
SC02	☐	3" Cushion			E2601/E2602	\$190
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion series			See pricelist	
HEIGHT ADJUSTABLE WITH PUSH HANDLES OPTIONS (MUST select ONE push handle options)						
BC06-912	☐	Straight Carbon fiber back cane 9" - 12"	increment of 1/2"	Set at ² :		NC
BC06-1215	☐	Straight Carbon fiber back cane 12" - 15"	increment of 1/2"	Set at ² :		NC
BC06-1518	☐	Straight Carbon fiber back cane 15" - 18"	increment of 1/2"	Set at ² :		NC
BC06-1821	☐	Straight Carbon fiber back cane 18" - 21"	increment of 1/2"	Set at ² :		NC
BC08-1315	☐	8° Bend Carbon fiber back cane 13.5" - 15"	increment of 1/2"	Set at ² :		\$60
BC08-1518	☐	8° Bend Carbon fiber back cane 15" - 18"	increment of 1/2"	Set at ² :		\$60
BC08-1821	☐	8° Bend Carbon fiber back cane 18" - 21"	increment of 1/2"	Set at ² :		\$60
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC09-1215	☐	Straight Carbon fiber back cane 12" - 15"	increment of 1/2"	Set at ² :		\$215
BC09-1518	☐	Straight Carbon fiber back cane 15" - 18"	increment of 1/2"	Set at ² :		\$215
BC09-1821	☐	Straight Carbon fiber back cane 18" - 21"	increment of 1/2"	Set at ² :		\$215
BC10-1315	☐	8° bend Carbon fiber back cane 13.5" - 15"	increment of 1/2"	Set at ² :		\$270
BC10-1518	☐	8° bend Carbon fiber back cane 15" - 18"	increment of 1/2"	Set at ² :		\$270
BC10-1821	☐	8° bend Carbon fiber back cane 18" - 21"	increment of 1/2"	Set at ² :		\$270
		Select back cane angle: 85° ☐ 90° ☐ 95° STD ☐ 100° ☐ 105° ☐ 110° ☐				
PUSH HANDLES OPTIONS (ONE option at a time)						
PH01	☐	Push handles				NC
PH05	☐	Fold down push handles ^{3,5}				\$145
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" ^{3, 4, 5}				\$325
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" ^{3, 4, 5}				\$325
PH02	☐	Omit push handles				NC

1- Not available with 8 ° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes and 9-12" back canes; 5- Not compatible with headrest and attendant lock

ARMRESTS						
HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ¹						
AR01-812	☐	10" Desk length armpad 8" - 12"	Locked ☐ Unlocked STD ☐			\$355
AR01-1014	☐	10" Desk length armpad 10" - 14"	Locked ☐ Unlocked STD ☐			\$355
AR02-812	☐	14" Full length armpad 8" - 12"	Locked ☐ Unlocked STD ☐			\$355
AR02-1014	☐	14" Full length armpad 10" - 14"	Locked ☐ Unlocked STD ☐			\$355
HEIGHT ADJUSTABLE "T" ARMREST ¹						
AR13-610	☐	10" Desk length armpad ⁴ 6" - 10"			E0973	\$330
AR03-812	☐	10" Desk length armpad 8" - 12"			E0973	\$330
AR03-1014	☐	10" Desk length armpad 10" - 14"			E0973	\$330
AR14-610	☐	14" Full length armpad ⁴ 6" - 10"			E0973	\$330
AR04-812	☐	14" Full length armpad 8" - 12"			E0973	\$330
AR04-1014	☐	14" Full length armpad 10" - 14"			E0973	\$330
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²						
AR05	☐	Tubular armrest without plastic side guards - Short pad	8" - 11"		E0973	\$230
AR19	☐	Tubular armrest with plastic side guards - Short pad	8" - 11"		E0973	\$395
AR23	☐	Tubular armrest without plastic side guards - Long pad	8" - 11"		E0973	\$230
AR24	☐	Tubular armrest with plastic side guards - Long pad	8" - 11"		E0973	\$395
AR06	☐	10" Desk length flat armpad without plastic side guards	8" - 11"		E0973	\$230
AR20	☐	10" Desk length flat armpad with plastic side guards	8" - 11"		E0973	\$395
AR17	☐	14" Full length flat armpad without plastic side guards	8" - 11"		E0973	\$230
AR21	☐	14" Full length flat armpad with plastic side guards	8" - 11"		E0973	\$395
ARMREST OPTIONS						
AR12	☐	14" Full Length armpad - Gel Ovation ³				\$165
AR99	☐	Omit armrest				NC

1- May have a 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with fixed carbon fiber side guards

3 - Only available with full length armrest; Not available with swing away and desk length armrests; 4- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES									
SIDE GUARDS (not available with Flip Back and "T" Armrest)									
SG01	☐	Plastic side guards (removable)	Straight bracket - STD ☐ Offset bracket ¹ ☐					K0108	\$200
SG02	☐	Carbon fiber side guards (fixed)							\$270
SG03	☐	Carbon fiber side guards (removable)	Straight bracket - STD ☐ Offset bracket ¹ ☐					K0108	\$305
SG06	☐	Plastic fender side guards (removable)	Straight bracket - STD ☐ Offset bracket ¹ ☐					K0108	\$290
SG04	☐	Carbon fiber fender side guards (removable)	Straight bracket - STD ☐ Offset bracket ¹ ☐						\$550
SG99	☐	Omit side guard							NC

1- only compatible with 0" camber; Increase width by 1/2" per side for a total of 1"

POSITIONING BELT									
SB01-S	☐	Velcro adjustable - Small 35" ¹						E0978	\$46
SB01-L	☐	Velcro adjustable - Long 60"						E0978	\$46
SB02-S	☐	Auto buckle - Small 35" ¹						E0978	\$87
SB02-M	☐	Auto buckle - Medium 48"						E0978	\$87
SB02-L	☐	Auto buckle - Long 60"						E0978	\$87
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"						E0978	\$85
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"						E0978	\$85
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"						E0978	\$180
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"						E0978	\$180
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"						E0978	\$180
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"						E0978	\$195
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"						E0978	\$195
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"						E0978	\$195
SB03	☐	Hardware attachment for 2 extra points (1 pair)							\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"						K0108	\$305
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"						K0108	\$305
SB99	☐	Omit positioning belt							STD

1- not recommended for seat width of 18" and more



HEADREST ¹									
AC01-S	☐	Profiled Headrest - Small (3.5"X6")						E0955	\$185
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")						E0955	\$185
AC01-L	☐	Profiled Headrest - Large (5.5"X8")						E0955	\$185
AC28	☐	Headrest mounting hardware with stabilizer bar 1						E1028	\$300

1- Only available with push handles

ACCESSORIES									
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6" camber							\$36
AC03	☐	Anti-tippers						E0971	\$150
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐					\$190
AC27	☐	Tip assist ²							\$67
AC06	☐	Cane and crutch holder	Left ☐	Right ☐				E2207	\$200
AC08	☐	Spoke guards ³	Available: 22", 24", 25", 26"					K0065	\$118
AC10	☐	Seat pouch							\$82
AC26	☐	BodyPoint Mobility Bag							\$135
AC23	☐	Newton Backpack							\$130
AC30	☐	Newton multi-purpose pouch							\$61
AC12	☐	Velcro Closing Strap							\$57
AC13	☐	Tool kit							\$51
AC15	☐	Military Patches	Army ☐	Navy ☐	Air Force ☐	Marines ☐	Coast Guard ☐		\$96
		Combat Wounded Purple Heart ☐	Space Force ☐	American Flag ☐					

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate; 2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerjy; the 24" spoke guard will be installed on the 25" wheel

Colors																			
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)
																			
AirGrip	☐		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame (NC)	☐		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Highlight ²	☐		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color	Specify PPG Paintit color:															additional 10 days lead-time	
																			\$390

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- To see different highlight, visit our website: https://www.motioncomposites.com/en_us/products/folding-wheelchairs/veloce

3- The frame color will be confirmed with the quote conversion into an order (in absence of selection, the color will be black)

Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD	N/A						
1							N/A
2						N/A	N/A
-1	N/A	N/A					
-2	N/A	N/A	N/A				

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° camber		
One arm drive: not compatible with "T"		
One arm drive and flip-back armrest only available with 17 in depth frame		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4



CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :



FOLDING ON CHAIR – USA ORDER FORM



Use wheelchair configuration
Or use this size:
NXCAF
(W" x D")
HCPCS CODE E2607-E2608

DEPTH	WIDTH											
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
12"	\$468	\$468	\$468		N/A	N/A	N/A	N/A	N/A	N/A	N/A	
13"						N/A	N/A	N/A	N/A	N/A	N/A	
14"	\$468	\$468	\$468									
15"		\$468	\$468	\$468	\$468							
16"		\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	
17"	N/A			\$468	\$468	\$468	\$468	\$468	\$468			
18"	N/A			\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	
19"	N/A			\$468	\$468	\$468	\$468	\$468	\$468			
20"	N/A			\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	

Option: Change cover for incontinent/infection control outer cover ☐ ☐ N/C

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of \$160.00 and has a 2 to 4 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES

	MODEL #	Description	Fits Back Width	Retail	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$80	

HCPCS CODE E0992

THORACIC SHAPE	 Optima™ THORACIC BACK SUPPORT HCPCS CODE: E2620-E2621							
	WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
	14"	13"	NOACB1413	\$730	☐	NIB-2193V	\$130	☐
	14"	16"	NOACB1416	\$730	☐	NIB-2194V	\$130	☐
	15"	13"	NOACB1513	\$730	☐	NIB-2197V	\$130	☐
	15"	16"	NOACB1516	\$730	☐	NIB-2198V	\$130	☐
	16"	13"	NOACB1613	\$730	☐	NIB-2201V	\$130	☐
	16"	16"	NOACB1616	\$730	☐	NIB-2202V	\$130	☐
	16"	19"	NOACB1619	\$730	☐	NIB-2203V	\$130	☐
	17"	13"	NOACB1713	\$730	☐	NIB-2206V	\$130	☐
	17"	16"	NOACB1716	\$730	☐	NIB-2207V	\$130	☐
	18"	13"	NOACB1813	\$730	☐	NIB-2208V	\$130	☐
	18"	16"	NOACB1816	\$730	☐	NIB-2211V	\$130	☐
	18"	19"	NOACB1819	\$730	☐	NIB-2212V	\$130	☐
	19"	16"	NOACB1916	\$730	☐	NIB-2213V	\$130	☐
	19"	19"	NOACB1919	\$730	☐	NIB-2215V	\$130	☐
	20"	13"	NOACB2013	\$730	☐	NIB-2218V	\$130	☐
	20"	16"	NOACB2016	\$730	☐	NIB-2219V	\$130	☐
	20"	19"	NOACB2019	\$730	☐	NIB-2220V	\$130	☐
	21"	16"	NOACB2116	\$730	☐	NIB-2221V	\$130	☐
	21"	19"	NOACB2119	\$730	☐	NIB-2222V	\$130	☐
	22"	13"	NOACB2213	\$783	☐	NIB-2223V	\$130	☐
	22"	16"	NOACB2216	\$783	☐	NIB-2224V	\$130	☐
	22"	19"	NOACB2219	\$783	☐	NIB-2225V	\$130	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.

 Optima Carbon™ THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) HCPCS CODE: E2620-E2621							
WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOCFACB1413	\$972	☐	NIB-2193V	\$130	☐
14"	16"	NOCFACB1416	\$1,008	☐	NIB-2194V	\$130	☐
15"	13"	NOCFACB1513	\$972	☐	NIB-2197V	\$130	☐
15"	16"	NOCFACB1516	\$1,008	☐	NIB-2198V	\$130	☐
16"	13"	NOCFACB1613	\$972	☐	NIB-2201V	\$130	☐
16"	16"	NOCFACB1616	\$1,008	☐	NIB-2202V	\$130	☐
17"	13"	NOCFACB1713	\$972	☐	NIB-2206V	\$130	☐
17"	16"	NOCFACB1716	\$1,008	☐	NIB-2207V	\$130	☐
18"	13"	NOCFACB1813	\$972	☐	NIB-2211V	\$130	☐
18"	16"	NOCFACB1816	\$1,008	☐	NIB-2212V	\$130	☐
20"	13"	NOCFACB2013	\$972	☐	NIB-2218V	\$130	☐
20"	16"	NOCFACB2016	\$1,008	☐	NIB-2219V	\$130	☐

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

THORACIC SHAPE			Optima™ DEEP THORACIC BACK SUPPORT HCPCS CODE: E2620-E2621					
	WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
	14"	13"	NOPCB1413	\$730	☐	NIB-2246V	\$130	☐
	14"	16"	NOPCB1416	\$730	☐	NIB-2247V	\$130	☐
	15"	13"	NOPCB1513	\$730	☐	NIB-2250V	\$130	☐
	15"	16"	NOPCB1516	\$730	☐	NIB-2251V	\$130	☐
	16"	13"	NOPCB1613	\$730	☐	NIB-2254V	\$130	☐
	16"	16"	NOPCB1616	\$730	☐	NIB-2255V	\$130	☐
	16"	19"	NOPCB1619	\$730	☐	NIB-2256V	\$130	☐
	17"	13"	NOPCB1713	\$730	☐	NIB-2259V	\$130	☐
	17"	16"	NOPCB1716	\$730	☐	NIB-2260V	\$130	☐
	18"	13"	NOPCB1813	\$730	☐	NIB-2264V	\$130	☐
	18"	16"	NOPCB1816	\$730	☐	NIB-2265V	\$130	☐
	18"	19"	NOPCB1819	\$730	☐	NIB-2266V	\$130	☐
	19"	16"	NOPCB1916	\$730	☐	NIB-2267V	\$130	☐
	19"	19"	NOPCB1919	\$730	☐	NIB-2268V	\$130	☐
	20"	13"	NOPCB2013	\$730	☐	NIB-2271V	\$130	☐
	20"	16"	NOPCB2016	\$730	☐	NIB-2272V	\$130	☐
	20"	19"	NOPCB2019	\$730	☐	NIB-2273V	\$130	☐
	21"	16"	NOPCB2116	\$730	☐	NIB-2274V	\$130	☐
	21"	19"	NOPCB2119	\$730	☐	NIB-2275V	\$130	☐
	22"	13"	NOPCB2213	\$783	☐	NIB-2276V	\$130	☐
	22"	16"	NOPCB2216	\$783	☐	NIB-2277V	\$130	☐
	22"	19"	NOPCB2219	\$783	☐	NIB-2278V	\$130	☐

		Optima Carbon™ DEEP THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) HCPCS CODE: E2620-E2621					
WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOCFPCB1413	\$1,066	☐	NIB-2246V	\$130	☐
14"	16"	NOCFPCB1416	\$1,108	☐	NIB-2247V	\$130	☐
15"	13"	NOCFPCB1513	\$1,066	☐	NIB-2250V	\$130	☐
15"	16"	NOCFPCB1516	\$1,108	☐	NIB-2251V	\$130	☐
16"	13"	NOCFPCB1613	\$1,066	☐	NIB-2254V	\$130	☐
16"	16"	NOCFPCB1616	\$1,108	☐	NIB-2255V	\$130	☐
17"	13"	NOCFPCB1713	\$1,066	☐	NIB-2259V	\$130	☐
17"	16"	NOCFPCB1716	\$1,108	☐	NIB-2260V	\$130	☐
18"	13"	NOCFPCB1813	\$1,066	☐	NIB-2264V	\$130	☐
18"	16"	NOCFPCB1816	\$1,108	☐	NIB-2265V	\$130	☐
20"	13"	NOCFPCB2013	\$1,066	☐	NIB-2271V	\$130	☐
20"	16"	NOCFPCB2016	\$1,108	☐	NIB-2272V	\$130	☐

THORACIC SHAPE			Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT HCPCS CODE E2620-E2621					
	xtend™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
	NHAACB(L)1213/17	12"	13-17"	\$730	☐	NIB-2148V	\$130	☐
	NHAACB(L)1413/17	14"	13-17"	\$730	☐	NIB-2149V	\$130	☐
	NHAACB(L)1613/17	16"	13-17"	\$730	☐	NIB-2150V	\$130	☐
	NHAACB(L)1813/17	18"	13-17"	\$730	☐	NIB-2151V	\$130	☐
	NHAACB (L)2013/17	20"	13-17"	\$730	☐	NIB-2152V	\$130	☐
	NHAACB (L)2213/17	22"	13-17"	\$783	☐	NIB-2153V	\$130	☐

		Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT HCPCS CODE E2620-E2621					
xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NHAPCB1616/20	16"	16"-20"	\$730	☐	NIB-2143V	\$130	☐
NHAPCB1816/20	18"	16"-20"	\$730	☐	NIB-2144V	\$130	☐
NHAPCB2016/20	20"	16"-20"	\$730	☐	NIB-2145V	\$130	☐
NHAPCB2216/20	22"	16"-20"	\$783	☐	NIB-2146V	\$130	☐

PELVIC SHAPE			Classic™ PELVIC BACK SUPPORT HCPCS CODE E2620-E2621					
	Classic™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
	NACB1616	16"	16"	\$730	☐	NIB-2163V	\$130	☐
	NACB1618	16"	18"	\$730	☐	NIB-2164V	\$130	☐
	NACB1816	18"	16"	\$730	☐	NIB-2166V	\$130	☐
	NACB1818	18"	18"	\$730	☐	NIB-2167V	\$130	☐
	NACB1820	18"	20"	\$730	☐	NIB-2168V	\$130	☐
	NACB2020	20"	20"	\$730	☐	NIB-2171V	\$130	☐
	NACB2216	22"	16"	\$783	☐	NIB-2172V	\$130	☐
	NACB2218	22"	18"	\$783	☐	NIB-2173V	\$130	☐
	NACB2220	22"	20"	\$783	☐	NIB-2174V	\$130	☐

		Classic™ DEEP PELVIC BACK SUPPORT HCPCS CODE E2620-E2621					
Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NPCB1616	16"	16"	\$730	☐	NIB-2175V	\$130	☐
NPCB1618	16"	18"	\$730	☐	NIB-2176V	\$130	☐
NPCB1818	18"	18"	\$730	☐	NIB-2179V	\$130	☐
NPCB1820	18"	20"	\$730	☐	NIB-2180V	\$130	☐
NPCB2020	18"	20"	\$730	☐	NIB-2183V	\$130	☐
NPCB2216	22"	16"	\$783	☐	NIB-2184V	\$130	☐
NPCB2218	22"	18"	\$783	☐	NIB-2185V	\$130	☐
NPCB2220	22"	20"	\$783	☐	NIB-2186V	\$130	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$70	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

HCPCS CODE: E0955

	MODEL	Description	RETAIL	✓
	NIB-0784	10" x 5" Large Rectangular	\$205	
	NIB-0818	14" x 5" Extra-Large Rectangular	\$205	

MultiFit™ HEAD SUPPORT – FIXTURES

HCPCS CODE HARDWARE: E1033

	MODEL	Description	RETAIL	✓
	NIB-0747	Multi-Axis Hardware	\$150	
	NIB-0823	Head Support Receiver Bracket	\$60	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	NXCCBF_____V	\$130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	NXCCBF_____V	\$150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF_____V	\$130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF_____V	\$150	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com