



# Zwim

Pool Chair



Zwim offers a safe way to enter and exit a pool that has a ramp or for transportation from the locker room to a pool lift. A contoured seat (with drainage slot) provides comfort. Ideal for chlorinated environments.



330  
lbs

Weight capacity

# ZWIM STANDARD FEATURES



SPECIFICATIONS	ZWIM POOL CHAIR
WEIGHT CAPACITY	330 lbs
SEAT WIDTH	19.5"
SEAT DEPTH (FRONT OF BACK POST TO FRONT OF SEAT)	17"
BACK WIDTH	18.5"
BACK HEIGHT (FROM TOP OF SEAT TO TOP OF UPHOLSTERY)	17"
ARM SUPPORT HEIGHT ABOVE SEAT/ARM SEAT	9"
SEAT-TO-FLOOR HEIGHT* (2.5 CM [1"] INCREMENTS)	19" - 23"
OVERALL WIDTH/OVERALL WIDTH	26.5"
PRODUCT WEIGHT/PRODUCT WEIGHT	39.5"
FRAME MATERIAL	ALUMINUM