



Helio A7 - FiveDay

ORDER FORM - CANADA

P.O. number:* _____ Date: _____
 Company Name:* _____ Phone:* _____
 Ordered by:* _____
 Clinician/ATP/RTS: _____
 Email: _____
 Marked for: _____

Ship to:* _____
 Name: _____
 Address: _____
 City: _____ Prov: _____
 Postal Code: _____
 Phone: _____

* Required A7MA080 English Manual (STD) A7MA081 Manuel Français Taxes not included

MODEL								
1	A7WC02	Ultralight Helio A7 wheelchair 7000 series aluminum (265 lb weight capacity)				WM0001020	\$2 645	<input type="checkbox"/>
2	A7WC09	HD Heavy duty kit (350 lb weight capacity)				WAMF	\$775	<input type="checkbox"/>
3	A7WC10	E-Motion adaptor Kit (2 adaptors)					\$75	<input type="checkbox"/>
SEAT WIDTH ³								
4	A7SW	Seat width:	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" ¹ <input type="checkbox"/>	
5	A7SW	Seat width:	19" ¹ <input type="checkbox"/>	20" ¹ <input type="checkbox"/>	21" ² <input type="checkbox"/>	22" ² <input type="checkbox"/>		WAMA \$295 <input type="checkbox"/>
SEAT DEPTH								
6	A7SD	Seat depth:	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <input type="checkbox"/>	NC <input type="checkbox"/>
7	A7SD	Seat depth:	19" <input type="checkbox"/>	20" <input type="checkbox"/>				WAMB \$225 <input type="checkbox"/>

1 - Available with HD option 2 - Only available with HD option [3. Width calculator](#)

A7 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13	14	15	13 1/2	14 1/2	15 1/2												
	P2	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2									
5"	P1																		
	P2				14	15	16	14 1/2	15 1/2	16 1/2									
	P3	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2						
	P4	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
7"	P1										16 1/2	17 1/2	18 1/2	17	18	19			
	P2									16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	
	P3							17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2
	P4	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21
	P5	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21	19 1/2	20 1/2	21 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
8	A7CA01	Composite wheel w/ PU tire (1" width)	3" <input type="checkbox"/>	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>		NC	<input type="checkbox"/>	
9	A7CA10	Pneumatic casters 6" x 1 1/4"								\$160	<input type="checkbox"/>	
10	A7CA11	Pneumatic casters 8" x 1 1/4"								\$160	<input type="checkbox"/>	
11	A7CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 1/2"								\$98	<input type="checkbox"/>
12	A7CA15		5" x 1 1/2"								\$98	<input type="checkbox"/>
13	A7CA16		6" x 1 1/2"								\$98	<input type="checkbox"/>
14	A7CA17		4" x 1 1/2"								\$175	<input type="checkbox"/>
15	A7CA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 1/2"								\$185	<input type="checkbox"/>
16	A7CA19		6" x 1 1/2"								\$195	<input type="checkbox"/>
CASTER OPTIONS												
17	A7CA12	Caster Pin Locks (<i>Available with +1" or +2" Stem bolt ONLY</i>)							WAMS	\$255	<input type="checkbox"/>	
18	A7CA13	Carbon fiber Frog Legs suspension forks (<i>Not available with caster pin locks</i>)								\$489	<input type="checkbox"/>	

A7 - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height: _____
20" wheel <input type="checkbox"/>	12 1/2 to 17 1/4	
22" wheel <input type="checkbox"/>	13 1/2 to 18 1/4	
24" wheel <input type="checkbox"/>	14 1/2 to 19 1/4	
25" wheel <input type="checkbox"/>	15 to 19 3/4	
26" wheel <input type="checkbox"/>	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

SWING-IN/SWING-OUT FOOTREST											
19	A7FR03	60° Swing-in/swing-out footrest ¹	13.5 - 17" <input type="checkbox"/>	14.5 - 19.5" <input type="checkbox"/>	15.5 - 20.5" <input type="checkbox"/>				Set at:*	NC	<input type="checkbox"/>
20	A7FR04	70° Swing-in/swing-out footrest ¹	13.5 - 17" <input type="checkbox"/>	14.5 - 19.5" <input type="checkbox"/>	15.5 - 20.5" <input type="checkbox"/>				Set at:*	NC	<input type="checkbox"/>
21	A7FR15	80° Swing-in/swing-out footrest ¹	13.5 - 16" <input type="checkbox"/>	14.5 - 19.5" <input type="checkbox"/>	15.5 - 20.5" <input type="checkbox"/>				Set at:*	NC	<input type="checkbox"/>
FOOTREST OPTIONS											
22	A7FR09	High mount (<i>footplate attached to the hanger with clamps</i>)	5 - 13.5" <input type="checkbox"/>						Set at:*	\$90	<input type="checkbox"/>

1 - For a 5-13.5" range, see line 22

* +/- 1/2"

90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT (<i>Footplate mounted to hanger tube</i>)											
23	A7FR06	90° swing-in/swing-out with front mount - Standard (no extension) ²	7.5 - 11"						Set at:*	NC	<input type="checkbox"/>
24	A7FR07	90° swing-in/swing-out with front mount - Short extension ^{2,3}	9.5 - 13"						Set at:*	NC	<input type="checkbox"/>
25	A7FR08	90° swing-in/swing-out with front mount - Medium extension ^{4,2,3}	9.5 - 14.5"						Set at:*	NC	<input type="checkbox"/>
26	A7FR99	Omit Footrest								NC	<input type="checkbox"/>

1 - Available only with adjustable footplates 2 - May interfere with front casters

* +/- 1/2"

ELEVATING LEGREST (<i>unit</i>) The elevating legrest replace the footrest side selected											
27	A7FR10	Elevating legrest with calf pad, 14" - 19" Length adjustable	Left <input type="checkbox"/>		Right <input type="checkbox"/>				Set at:	WAMK	\$290 <input type="checkbox"/>

FOOTPLATES											
28	A7FP01	Standard - Composite								STD	<input type="checkbox"/>
29	A7FP02	Newton adjustable angle - Composite							WAMJ	\$145	<input type="checkbox"/>
30	A7FP03	Newton adjustable angle - Aluminum							WAMJ	\$165	<input type="checkbox"/>
31	A7FP04	Newton adjustable angle - Carbon fiber (<i>with rubber inserts</i>)							WAMJ	\$380	<input type="checkbox"/>
32	A7FP06	One piece adjustable angle flip-up - Aluminum						WAMJ	SEFND0015	\$206	<input type="checkbox"/>
33	A7FP07	One piece adjustable angle flip-up - Carbon fiber (<i>with rubber inserts</i>)						WAMJ	SEFND0015	\$450	<input type="checkbox"/>

FOOTPLATE ACCESSORIES											
34	A7FP08	Length adjustable heel loop (<i>not available with one piece footplate</i>)								NC	<input type="checkbox"/>
35	A7FP09	Calf strap							SEFND0010	\$52	<input type="checkbox"/>
36	A7FP30	Bodypoint Aeromesh Padded calf strap							SEFND0010	\$129	<input type="checkbox"/>

RESIDUAL LIMB SUPPORT (<i>Height and angle adjustable</i>) - (<i>unit</i>)											
37	A7FR23	Swing away residual limb support - 7" x 9"	Left <input type="checkbox"/>		Right <input type="checkbox"/>				SEFND0005	\$390	<input type="checkbox"/>
38	A7FR12	Swing away residual limb support - 8" x 11"	Left <input type="checkbox"/>		Right <input type="checkbox"/>				SEFND0005	\$390	<input type="checkbox"/>

Back canes with integrated fixed height push handle

FIXED ANGLE BACK										
39	A7BC01	8° bend aluminum back cane	16" <input type="checkbox"/>	18" <input type="checkbox"/>	20" <input type="checkbox"/>	STD			NC	<input type="checkbox"/>

Back canes, height adjustable, with option

FIXED ANGLE BACK										
40	A7BC11	Straight aluminum back cane	9-12" <input type="checkbox"/>	12 - 15" <input type="checkbox"/>	15 - 18" <input type="checkbox"/>	18 - 21" <input type="checkbox"/>	STD		\$75	<input type="checkbox"/>
41	A7BC12	8° bend aluminum back cane	13.5 - 15" <input type="checkbox"/>	15 - 18" <input type="checkbox"/>	18 - 21" <input type="checkbox"/>	STD			\$75	<input type="checkbox"/>

ADJUSTABLE ANGLE BACK											
42	A7BC13	Adjustable angle, straight aluminum back cane				15 - 18" <input type="checkbox"/>	18 - 21" <input type="checkbox"/>	STD		\$259	<input type="checkbox"/>
43	A7BC14	Adjustable angle, 8° bend aluminum back cane				15 - 18" <input type="checkbox"/>	18 - 21" <input type="checkbox"/>	STD		\$309	<input type="checkbox"/>
44		Select back cane angle:	85° ¹ <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>			

Set height to: ² _____ (measured from seat sling to top of back upholstery)

PUSH HANDLE OPTION										
45	A7PH01	Push handles							STD	<input type="checkbox"/>
46	A7PH05	Fold down push handles ^{3,5}							\$150	<input type="checkbox"/>
47	A7PH06	Clamp on height adjustable push handle ^{3,4,5}	9 - 20" STD <input type="checkbox"/>				9-23" <input type="checkbox"/>		\$325	<input type="checkbox"/>
48	A7PH02	Without push handles							NC	<input type="checkbox"/>

1 - Not available with a 8 degree bend cane 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit

4 - Not compatible with swing away armrests, 8° bend back canes and amputee axle plates.

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹ (only available with push handle)										
49	A7PH04	Stroller handle (not available with the HD option)						WAMW	\$210	<input type="checkbox"/>
50	A7PH03	Newton folding stabilizer bar							\$258	<input type="checkbox"/>

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY											
51	A7BH07	Slip-on Active-Mesh back upholstery							NC	<input type="checkbox"/>	
52	A7BH08	Adjustable tension Active-Mesh back upholstery						Product improved	WAMD/WAME	\$245	<input type="checkbox"/>
53	A7BH99	Omit Back upholstery							STD	<input type="checkbox"/>	

RIGID BACK										
54		nXt back series (See nXt back pricelist)				See pricelist				<input type="checkbox"/>

SEAT CUSHION										
55	A7SC01	2" Cushion							\$175	<input type="checkbox"/>
56	A7SC02	3" Cushion							\$195	<input type="checkbox"/>
57	A7SC99	Omit seat cushion							NC	<input type="checkbox"/>
58		nXt FIT cushion series (See nXt FIT pricelist)				See pricelist				<input type="checkbox"/>

ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims										
59	A7HR22	One arm drive ^{2,3}	24"	left <input type="checkbox"/>	right <input type="checkbox"/>				\$1 395	<input type="checkbox"/>
60	A7HR98	Omit handrim on opposite side of the One Arm Drive (OAD)							NC	<input type="checkbox"/>

1 - Offered with handrims on each side of the wheelchair 2 - Only standard handrim position available 3 - Non disponible avec l'option HD

REAR WHEELS											
61	A7RW01	Mag ¹	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				\$129	<input type="checkbox"/>	
62	A7RW02	Newton One - spoke wheel ²	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				NC	<input type="checkbox"/>	
63	A7RW03	Newton Gravity - ultralight wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>		\$360	<input type="checkbox"/>	
64	A7RW04	Spinergy Spox	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>			\$800	<input type="checkbox"/>	
65		Spinergy Spox spokes colors ³	Bik STD <input type="checkbox"/>	Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>		NC	<input type="checkbox"/>	
66	A7RW05	Spinergy LX	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>			\$995	<input type="checkbox"/>	
67		Spinergy LX spokes colors ³	Bik STD <input type="checkbox"/>	Pur <input type="checkbox"/>	Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>	Pink <input type="checkbox"/>	Grn <input type="checkbox"/>	Org <input type="checkbox"/>
68	A7RW07	Spinergy CLX Carbon Blade	24" <input type="checkbox"/>	25" <input type="checkbox"/>					\$1 540	<input type="checkbox"/>	
69		Spinergy CLX spokes colors ³	Bik STD <input type="checkbox"/>	Pur <input type="checkbox"/>	Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>	Pink <input type="checkbox"/>	Grn <input type="checkbox"/>	Org <input type="checkbox"/>

1 - Non compatible avec pneus HP ; Seulement disponible avec les cerceaux en aluminium ou plastifié en position standard 2 - Not compatible with Shox tires

3 - Colour spokes require an additionnal 10 days lead-time

TIRES			Available				
70	A7TI01	Soft Urethane 1 1/4 "	20", 22" 24"				STD <input type="checkbox"/>
71	A7TI02	Soft Urethane knobby 1 3/8 "	22", 24", 25", 26"				\$110 <input type="checkbox"/>
72	A7TI04	Soft Urethane - Shox 1"	22", 24", 25", 26"				\$130 <input type="checkbox"/>
73	A7TI05	Pneumatic 1 3/8 "	20", 22", 24", 25", 26"				NC <input type="checkbox"/>
74	A7TI08	Pneumatic with Airless insert 1 3/8"	22", 24"				\$215 <input type="checkbox"/>
75	A7TI06	Pneumatic HP - 1"	20", 22", 24", 25", 26"				\$110 <input type="checkbox"/>
76	A7TI07	Pneumatic HP - Schwalbe Marathon Plus 1"	24", 25", 26"				\$190 <input type="checkbox"/>

HANDRIM							
77	A7HR01	Aluminum anodized					STD <input type="checkbox"/>
78	A7HR07	Newton Air Grip ¹ Black <input type="checkbox"/>				WAMO	\$150 <input type="checkbox"/>
79	A7HR	Newton Air Grip ¹ (08) Sunkissed Orange <input type="checkbox"/> (09) Burgundy <input type="checkbox"/> (10) Sapphire Blue <input type="checkbox"/> (11) Monster Green <input type="checkbox"/> (12) Charcoal <input type="checkbox"/> (13) Fuchsia <input type="checkbox"/> (14) Steel Blue <input type="checkbox"/> (15) Ferrari Red <input type="checkbox"/> (16) Acid Green <input type="checkbox"/> (17) White <input type="checkbox"/> (18) Sky Blue <input type="checkbox"/>				WAMO	\$330 <input type="checkbox"/>
Assembly position for aluminum and Newton AirGrip handrims							
80	A7HR29	Standard Handrim position Full tab					STD <input type="checkbox"/>
81	A7HR30	Narrow Handrim position Short tab 2 positions					\$25 <input type="checkbox"/>
82	A7HR31	Super Narrow Handrim Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/>					\$25 <input type="checkbox"/>

HANDRIM							
83	A7HR02	Plastic coated ³				WAMO	\$160 <input type="checkbox"/>
84	A7HR03	Q-Grip High friction Neoprene coated ^{1,2,3}				WAMO	\$430 <input type="checkbox"/>
85	A7HR28	Surge ^{1,2,3} (1-1/2" Oval handrim w/ Gripton high friction strip)					\$470 <input type="checkbox"/>
86	A7HR05	Surge LT ^{1,2,3} (1-1/8" Oval handrim w/ Gripton high friction strip)					\$500 <input type="checkbox"/>
87	A7HR04	Natural Fit Standard Grip ^{1,2,3} (1-1/2" Oval handrim w/standard thumb piece)					\$500 <input type="checkbox"/>
88	A7HR99	Omit handrim					NC <input type="checkbox"/>
Thumb piece⁴ selection Natural Fit handrims							
89	A7HR23	Standard grip					\$190 <input type="checkbox"/>
90	A7HR27	Super grip					\$220 <input type="checkbox"/>
91	A7HR099	Omit thumb piece					NC <input type="checkbox"/>

1 - Not available with Mag wheels or OAD

2 - Not available in 20"

3 - Only standard handrim position available


4 - Compatible with Surge, Surge LT handrims

AXLE							
92	A7AX01	Quick release axle				WAMU	STD <input type="checkbox"/>
93	A7AX02	Permanent axle					NC <input type="checkbox"/>
94	A7AX03	Quad release axle				WAMU	\$170 <input type="checkbox"/>
AXLE PLATE							
95	A7AX04	Amputee adjustable axle plate ¹				WAMT	\$250 <input type="checkbox"/>

1 -May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER							
96	A7CM00	0° camber					STD <input type="checkbox"/>
97	A7CM03	3° camber					NC <input type="checkbox"/>
98	A7CM06	6° camber (not available with amputee axle plate)					\$80 <input type="checkbox"/>

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable								Most responsive			
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK									
99	A7WL08	Newton Wheel lock - push to lock						STD	<input type="checkbox"/>
100	A7WL02	Pull to lock					NC	<input type="checkbox"/>	
101	A7WL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>				WAMR	\$175	<input type="checkbox"/>	
102	A7WL13 A7WL14	6" removable extension handle					\$67	<input type="checkbox"/>	
103	A7WL05	Scissor Lock					\$175	<input type="checkbox"/>	
104	A7WL06	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>					\$230	<input type="checkbox"/>	
105	A7WL99	Omit wheel locks					NC	<input type="checkbox"/>	

FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ¹									
				Height Adjustment		Removable function			
106	A7AR01	10" Desk length arm pad		8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/>	STD	\$380 <input type="checkbox"/>
107	A7AR02	14" Full length arm pad		8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/>	STD	\$380 <input type="checkbox"/>

"T" ARMREST ¹									
				Height Adjustment					
108	A7AR03	10" Desk length arm pad		8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>				\$360 <input type="checkbox"/>
109	A7AR04	14" Full length arm pad		8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>				\$360 <input type="checkbox"/>

1 - May have a 1/2" adjustment restriction on both limits of the range

GEL ARM PAD									
110	A7AR12	14" Full Length Arm Pad - Gel Ovation					SEAND0030	\$150	<input type="checkbox"/>

SWING AWAY ARMREST ²									
				Height Adjustment					
111	A7AR05	Tubular arm pad		8 - 11"					\$295 <input type="checkbox"/>
112	A7AR06	10" Desk length arm pad		8 - 11"					\$295 <input type="checkbox"/>

113	A7AR99	Omit armrest						NC	<input type="checkbox"/>
-----	--------	--------------	--	--	--	--	--	----	--------------------------

2 - Could interfere with some rigid back rests

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>									
114	A7SG01	Plastic side guards <i>(removable)</i>					WAMV	\$225	<input type="checkbox"/>
115	A7SG02	Carbon fiber side guards <i>(fixed)</i> ¹					WAMV	\$255	<input type="checkbox"/>
116	A7SG03	Carbon fiber side guards <i>(removable)</i>					WAMV	\$295	<input type="checkbox"/>
117	A7SG06	Plastic fender side guards <i>(removable)</i> ²					WAMV	\$320	<input type="checkbox"/>
118	A7SG04	Carbon fiber fender side guards <i>(removable)</i>					WAMV	\$495	<input type="checkbox"/>
119	A7SG99	Omit side guard						NC	<input type="checkbox"/>

1 - Not available with back upholstery

2 - Not compatible with 20" and 22" rear wheels

POSITIONING BELT									
120	A7SB01	Velcro adjustable		Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>			\$45	<input type="checkbox"/>
121	A7SB02	Auto buckle		Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>			\$55	<input type="checkbox"/>
122	A7SB07	Hardware attachment for 2 point seatbelt (1 pair)						NC	<input type="checkbox"/>
123	A7SB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (1 pair)</i>						\$13	<input type="checkbox"/>
124	A7SB04	BodyPoint Hip Belt - 2 point not padded			Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>		\$82	<input type="checkbox"/>
125	A7SB05	BodyPoint Hip Belt - 2 point padded		Small 51" <input type="checkbox"/>	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>	SERND0030	\$170	<input type="checkbox"/>
126	A7SB06	BodyPoint Hip Belt - 4 point padded		Small 51" <input type="checkbox"/>	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>	SERND0020	\$195	<input type="checkbox"/>
127	A7SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps				Medium 18" <input type="checkbox"/>	Long 24" <input type="checkbox"/>	\$324	<input type="checkbox"/>
128	A7SB99	Omit positioning belt						STD	<input type="checkbox"/>

HEADREST ¹									
129	A7AC01	Profiled Headrest		Small (3.5"X6") <input type="checkbox"/>	Medium (4.5"X7") <input type="checkbox"/>	Large (5.5"X8") <input type="checkbox"/>		\$205	<input type="checkbox"/>
130	A7AC28	Headrest mounting hardware with removable non folding stabilizer bar						\$338	<input type="checkbox"/>


1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES														
131	A7AC03	Anti tippers								\$115	<input type="checkbox"/>			
132	A7AC27	Tip assist ¹								\$60	<input type="checkbox"/>			
133	A7AC04	Transit Tie Down (<i>Unoccupied</i>)								NC	<input type="checkbox"/>			
134	A7AC05	WC-19 Transit Tie Down w/Q'straint belt ²		Short belt 20-41" <input type="checkbox"/>		Long belt 20-46" <input type="checkbox"/>				\$345	<input type="checkbox"/>			
135	A7AC06	Cane carrier		Left <input type="checkbox"/>		Right <input type="checkbox"/>				\$105	<input type="checkbox"/>			
136	A7AC07	Oxygen tank holder - E Type		Left						WAMX	\$375	<input type="checkbox"/>		
<i>Not compatible with swing away armrest or amputee axle plate</i>														
137	A7AC08	Spoke guards ³		Available: 20", 22", 24", 26"						WAMP	\$145	<input type="checkbox"/>		
138	A7AC09	Footrest neoprene impact guard (pair) - Black								SEFND0035	\$80	<input type="checkbox"/>		
139	A7AC10	Seat pouch								\$95	<input type="checkbox"/>			
140	A7AC26	BodyPoint Mobility Bag								\$135	<input type="checkbox"/>			
141	A7AC23	Newton Backpack								\$120	<input type="checkbox"/>			
142	A7AC12	Velcro Closing Strap								\$50	<input type="checkbox"/>			
143	A7AC13	Tool kit								\$50	<input type="checkbox"/>			
144	A7AC16	Motion Composites Polo Shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>	M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>	M/XXXXL <input type="checkbox"/>	\$32	<input type="checkbox"/>

1 - Cannot be combined with Anti Tippers 2 - Not available with the swing away armrest ; Not available with HD option 3 - Not available with Mag wheels et Spinergy

FRAME COLOR											
145	A7FC01	Black <input type="checkbox"/>								NC	<input type="checkbox"/>
146	A7FC03	Sunkissed Orange <input type="checkbox"/>	Burgundy <input type="checkbox"/>	Sapphire Blue <input type="checkbox"/>	Monster Green <input type="checkbox"/>				additional 5 day lead-time	\$270	<input type="checkbox"/>
		Charcoal <input type="checkbox"/>	Fuchsia <input type="checkbox"/>	Steel Blue <input type="checkbox"/>	Ferrari Red <input type="checkbox"/>	Acid Green <input type="checkbox"/>					
147	A7FC99	Custom Color	Specify PPG Paintit color:						additional 15 days lead-time	\$410	<input type="checkbox"/>
TOTAL											

Fork and Stem bolt combinations



Notes :