



P.O. number: * _____ Date: _____

Company Name: * _____ Phone: * _____

Ordered by: * _____

Clinician/ATP/RTS: _____

Email: _____

Marked for: _____

Ship to: * _____

Name: _____

Address: _____

City: _____ Prov: _____

Postal Code: _____

Phone: _____

* Required MVMA080 English Manual (STD) MVMA081 Manuel Français Taxes not included

MODEL									
1	MVWC01	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)				WM0001196	\$1 395	☐	
2	MVWC09	Heavy duty kit (350 lb weight capacity)				WAMF	\$775	☐	
SEAT WIDTH ³									
3	MVSW	Seat width:	14" ☐	16" ☐	18" ¹ ☐			NC	☐
4	MVSW	Seat width:		20" ¹ ☐	22" ² ☐		WAMA	\$295	☐
SEAT DEPTH									
5	MVSD	Seat depth:		16" ☐	18" ☐			NC	☐
6	MVSD	Seat depth:		20" ☐			WAMB	\$225	☐

1 - Available with HD option

2 - Only available with HD option

[3. Width calculator](#)

Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS											
7	MVCA01	Composite wheel w/ PU tire (1" width)	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
8	MVCA10	Pneumatic 6" x 1 ¼"								\$160	☐
9	MVCA11	Pneumatic 8" x 1 ¼"								\$160	☐
10	MVCA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "							\$98	☐
11	MVCA15		5" x 1 ½ "							\$98	☐
12	MVCA16		6" x 1 ½ "							\$98	☐
13	MVCA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½"							\$175	☐
14	MVCA18		5" x 1 ½"							\$185	☐
15	MVCA19		6" x 1 ½"							\$195	☐

Move REAR Seat-to-floor Height - Please circle your choice

20"	22"	24"	25"	26"
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

SWING-AWAY FOOTREST									
16	MVFR01	60° swing away footrest ¹	13.5 - 17.5" ☐	15.5 - 19.5" ☐	17 - 21" ☐	Set at:*		\$250	☐
17	MVFR02	70° swing away footrest ¹	13.5 - 17.5" ☐	15.5 - 19.5" ☐	17 - 21" ☐	Set at:*		\$250	☐
		FOOTREST OPTIONS							
18	MVFR09	High mount (footplate attached to the hanger with clamps)	5 - 13.5" ☐			Set at:*		\$90	☐

1 - For a 5-13,5" range see line 18

* +/- 1/2"

ELEVATING LEGREST (unit)			The elevating legrest replace the footrest side selected					
19	MVFR10	Elevating legrest with calf pad, 14" - 19" Length adjustable	Left ☐	Right ☐	Set at:	WAMK	\$390	☐

FOOTPLATES									
20	MVFP01	Standard - Composite						STD	☐
21	MVFP02	Newton adjustable angle - Composite					WAMJ	\$145	☐
22	MVFP03	Newton adjustable angle - Aluminum					WAMJ	\$165	☐

FOOTPLATE ACCESSORIES									
23	MVFP08	Length adjustable heel loop						\$25	☐
24	MVFP09	Calf strap					SEFND0010	\$52	☐
25	MVFP30	Bodypoint Aeromesh Padded calf strap					SEFND0010	\$129	☐

RESIDUAL LIMB SUPPORT (Height and angle adjustable) (unit)									
26	MVFR23	Swing away residual limb support - 7" x 9"	Left ☐	Right ☐			SEFND0005	\$390	☐
27	MVFR12	Swing away residual limb support - 8" x 11"	Left ☐	Right ☐			SEFND0005	\$390	☐

Back canes with integrated fixed height push handle									
FIXED ANGLE BACK									
28	MVBC01	8° bend aluminum back cane	16" ☐	18" ☐	20" STD ☐			NC	☐
ADJUSTABLE ANGLE BACK									
29	MVBC02	8° bend aluminum back cane	16" ☐	20" STD ☐				\$309	☐
30		Select back cane angle:	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐		

STROLLER HANDLE AND STABILIZER BAR ¹									
31	MVPH04	Stroller handle (not available with the HD option)					WAMW	\$210	☐
32	MVPH03	Folding stabilizer bar						\$258	☐

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY						
33	MVBH05	Slip-on Nylon back upholstery			NC	☐
34	MVBH06	Adjustable tension Nylon back upholstery	Product improved	WAMD/WAME	\$245	☐
35	MVBH99	Omit back upholstery			STD	☐
RIGID BACK						
36		nXt back series (See nXt back pricelist)	See pricelist			☐

SEAT CUSHIONS						
37	MVSC01	2" Cushion			\$175	☐
38	MVSC02	3" Cushion			\$195	☐
39	MVSC99	Omit seat cushion			NC	☐
40		nXt FIT cushion series (See nXt FIT pricelist)	See pricelist			☐

ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims						
41	MVHR22	One arm drive ^{2,3} 24" left ☐ right ☐			\$1 395	☐
42	MVHR98	Omit handrim on opposite side of the One Arm Drive (OAD)			NC	☐

1- Offered with handrims on each side of the wheelchair 2- Only standard handrim position available 3- Not available with HD option

REAR WHEELS						
43	MVRW01	Mag ¹ 20" ☐ 22" ☐ 24" ☐			NC	☐
44	MVRW02	Newton One - spoke wheel ² 20" ☐ 22" ☐ 24" ☐			NC	☐
45	MVRW03	Newton Gravity - ultralight wheel 20" ☐ 22" ☐ 24" ☐ 25" ☐ 26" ☐			\$360	☐

1 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 2- Not compatible with Shox tires

TIRES			Available			
46	MVTI01	Soft Urethane 1 1/4 "	20", 22" 24"		STD	☐
47	MVTI02	Soft Urethane knobby 1 3/8 "	22", 24", 25", 26"		\$130	☐
48	MVTI04	Soft Urethane - Shox 1"	22", 24", 25", 26"		\$155	☐
49	MVTI05	Pneumatic 1 3/8 "	20", 22", 24", 25", 26"		\$130	☐
50	MVTI08	Pneumatic with Airless insert 1 3/8"	22", 24"		\$215	☐
51	MVTI06	Pneumatic HP - 1"	20", 22", 24", 25", 26"		\$150	☐
52	MVTI07	Pneumatic HP - Schwalbe Marathon Plus 1"	24", 25", 26"		\$230	☐

HANDRIM						
53	MVHR01	Aluminum anodized			STD	☐
54	MVHR07	Newton Air Grip ¹ Noir ☐		WAMO	\$150	☐
55	MVHR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐ (11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐ (15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐		WAMO	\$330	☐
Assembly position for aluminum and Newton AirGrip handrims						
56	MVHR29	Standard Handrim position Full tab			STD	☐
57	MVHR30	Narrow Handrim position Short tab 2 positions			\$25	☐
58	MVHR31	Super Narrow Handrim Short tab 2 positions STD ☐ Short tab 1 position ☐			\$25	☐


HANDRIM						
59	MVHR02	Plastic coated ³		WAMO	\$160	☐
60	MVHR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}		WAMO	\$430	☐
61	MVHR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)			\$470	☐
62	MVHR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)			\$500	☐
63	MVHR04	Natural Fit ^{1, 2, 3}			\$500	☐
64	MVHR99	Omit handrim			NC	☐
Thumb piece ⁴ selection Natural Fit handrims						
65	MVHR23	Standard grip			\$190	☐
66	MVHR27	Super grip			\$220	☐
67	MVHR099	Omit thumb piece			NC	☐

1 - Not available with Mag wheels or OAD 2 - Not available in 20" 3- Only standard handrim position available 4- Compatible with Surge, Surge LT handrims

AXLES (0° camber)					
68	MVAX01	Quick release axle		WAMU	STD ☐
69	MVAX02	Permanent axle			NC ☐
70	MVAX03	Quad release axle		WAMU	\$170 ☐
AXLE PLATE					
71	MVAX04	Amputee adjustable axle plate ¹		WAMT	\$250 ☐

¹ - May limit seat to floor height available

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																			</
-------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	----

WHEEL LOCKS					
72	MVWL08	Newton Wheel lock - push to lock			STD ☐
73	MVWL02	Pull to lock			NC ☐
74	MVWL04	Newton grade aid push to lock (<i>suggested with pneumatic or high pressure tire</i>)		WAMR	\$175 ☐
75	MVWL13 MVWL14	6" removable extension handle			\$67 ☐
76	MVWL06	Attendant Lock (<i>additional set of wheel lock mounted to the back</i>)			\$230 ☐
77	MVWL99	Omit wheel locks			NC ☐

FLIP BACK ARMREST (<i>Convertible to "T" armrest</i>) ¹			Height Adjustment		Removable function		
78	MVAR01	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ STD	\$380 ☐
79	MVAR02	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ STD	\$380 ☐
"T" ARMREST ¹			Height Adjustment				
80	MVAR03	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐			\$360 ☐
81	MVAR04	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐			\$360 ☐

¹ - May have a 1/2 " adjustment restriction on both limits of the range

GEL ARM PAD					
82	MVAR12	14" Full Length Arm Pad - Gel Ovation		SEAND0030	\$150 ☐
83	MVAR99	Omit armrest			NC ☐

POSITIONING BELTS					
84	MVSB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐	\$45 ☐
85	MVSB02	Auto buckle	Standard 35" ☐	Long 60" ☐	\$55 ☐
86	MVSB07	Hardware attachment for 2 point seatbelt (1 pair)			NC ☐
87	MVSB03	Additional hardware for 2 extra attachment points (<i>needed for the 4 point seat belt</i>) (1 pair)			\$13 ☐
88	MVSB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐	Long 61" ☐	\$82 ☐
89	MVSB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐ Long 61" ☐	SERND0030 \$170 ☐
90	MVSB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐ Long 61" ☐	SERND0020 \$195 ☐
91	MVSB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐	Long 24" ☐	\$324 ☐
92	MVSB99	Omit positioning belt			NC ☐

HEADREST ¹					
93	MVAC01	Profiled Headrest	Small (3.5"X6") ☐	Medium (4.5"X7") ☐	Large (5.5"X8") ☐
94	MVAC28	Headrest mounting hardware with stabiliser bar			\$338 ☐

¹ - Not available with a stroller handle or a folding stabiliser bar

