MOTION COMPOSITES
160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T.+1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



6.6

Helio A6

ORDER FORM - AADL #W332

1- Purchasing information

P.O. number:* Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

2- Dealer/Representative information

Dealer/Representative:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:*

Name:

Address:

City: Province:

Postal Code:

Phone:

** Required*

MA080 ☐ *English Manual (STD)* *MA081* ☐ *French Manual*

Taxes not included

MODEL					
WC19	☐	Ultralight folding Helio A6 Wheelchair - *6000 series aluminum (265lb weight capacity)*			**$2,445**
WC09	☐	**HD** Heavy duty kit *(350 lb weight capacity)*			$850
Frame					
SEAT WIDTH [3]					
SW	☐	14"☐ 15"☐ 16"☐ 17"☐ 18" [1]☐			NC
	☐	19" [1]☐ 20" [1]☐ 21" [2]☐ 22" [2]☐			$320
SEAT DEPTH					
SD	☐	14"☐ 15"☐ 16"☐ 17"☐ 18"☐			NC
	☐	19"☐ 20"☐			$280
TRANSIT TIE DOWN					
AC04	☐	Transit Tie Down *(Unoccupied)*			NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt [4]	Short belt 20"-41"	***Client cost***	$370
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt [4]	Long belt 20"-46"	***Client cost***	$370
AC35	☐	WC-19 Transit Tie Down without belt [4]		***Client cost***	$345

1- *Available with HD Kit (350 lb);* ***2-*** *Only with HD Kit (350 lb)* *3. Width calculator* ***4-*** *May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option*

A6 FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		**STD**	**+1"**	**+2"**	**STD**	**+1"**	**+2"**	**STD**	**+1"**	**+2"**	**STD**	**+1"**	**+2"**	**STD**	**+1"**	**+2"**	**STD**	**+1"**	**+2"**
4"	P1	12"☐	13"☐	14"☐	12 1/2"☐	13 1/2"☐	14 1/2"☐												
	P2	12 1/2"☐	13 1/2"☐	14 1/2"☐	13"☐	14"☐	15"☐		14 1/2"☐	15 1/2"☐									
5"	P1																		
	P2				13"☐	14"☐	15"☐	13 1/2"☐	14 1/2"☐	15 1/2"☐									
	P3	13"☐	14"☐	15"☐	13 1/2"☐	14 1/2"☐	15 1/2"☐	14"☐	15"☐	16"☐		15 1/2"☐	16 1/2"☐						
	P4	13 1/2"☐	14 1/2"☐	15 1/2"☐	14"☐	15"☐	16"☐	14 1/2"☐	15 1/2"☐	16 1/2"☐	15"☐	16"☐	17"☐		16 1/2"☐	17 1/2"☐			
7"	P1										15 1/2"☐	16 1/2"☐	17 1/2"☐	16"☐	17"☐	18"☐			
	P2							15 1/2"☐	16 1/2"☐	17 1/2"☐	16"☐	17"☐	18"☐	16 1/2"☐	17 1/2"☐	18 1/2"☐			
	P3							16"☐	17"☐	18"☐	16 1/2"☐	17 1/2"☐	18 1/2"☐	17"☐	18"☐	19"☐	17 1/2"☐	18 1/2"☐	19 1/2"☐
	P4	15 1/2"☐	16 1/2"☐	17 1/2"☐	16"☐	17"☐	18"☐	16 1/2"☐	17 1/2"☐	18 1/2"☐	17"☐	18"☐	19"☐	17 1/2"☐	18 1/2"☐	19 1/2"☐	18"☐	19"☐	20"☐
	P5	16"☐	17"☐	18"☐	16 1/2"☐	17 1/2"☐	18 1/2"☐	17"☐	18"☐	19"☐	17 1/2"☐	18 1/2"☐	19 1/2"☐	18"☐	19"☐	20"☐	18 1/2"☐	19 1/2"☐	20 1/2"☐

Height available +- 1/4" (see stem bolt and fork combination at the end of the order form); Orange boxes are not available with a seat slope of more than 1.5in.

FRONT CASTERS

CASTERS				Availability		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		$170
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		$120
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"	*Client cost*	$198
CASTER OPTIONS						
CA12	☐	Caster Pin Locks *(Available with +1" or +2" Stem bolt ONLY)* [1]			*Client cost*	$275
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)			*Client cost*	$525

1- Not available with WC-19

Helio A6 - REAR Seat-to-floor Height - Please select your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12" ☐	13" ☐	13 3/4" [1] ☐	14 1/4" ☐	14 3/4" ☐
12 1/2" ☐	13 1/2" ☐	14 1/4" ☐	14 3/4" ☐	15 1/4" ☐
13" ☐	14" ☐	14 3/4" ☐	15 1/4" ☐	15 3/4" ☐
13 1/2" ☐	14 1/2" ☐	15 1/4" ☐	15 3/4" ☐	16 1/4" ☐
14" ☐	15" ☐	15 3/4" ☐	16 1/4" ☐	16 3/4" ☐
14 1/2" ☐	15 1/2" ☐	16 1/4" ☐	16 3/4" ☐	17 1/4" ☐
15" ☐	16" ☐	16 3/4" ☐	17 1/4" ☐	17 3/4" ☐
15 1/2" ☐	16 1/2" ☐	17 1/4" ☐	17 3/4" ☐	18 1/4" ☐
16" ☐	17" ☐	17 3/4" ☐	18 1/4" ☐	18 3/4" ☐
16 1/2" ☐	17 1/2" ☐	18 1/4" ☐	18 3/4" ☐	19 1/4" ☐
17" ☐	18" ☐	18 3/4" ☐	19 1/4" ☐	19 3/4" ☐
17 1/2" ☐	18 1/2" ☐	19 1/4" [1] ☐	19 3/4" ☐	20 1/4" ☐

Configuration used: 1" seat slope with a center of gravity at 2" **1-** *Not available with One Arm Drive*

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

REAR WHEELS

ONE ARM DRIVE [1] *Only available with Newton One rear wheel, permanent axles and 0° camber*					
HR22	☐	One arm drive - aluminum inner handrim [2] Left☐ Right☐			$1,480
HR60	☐	One arm drive - plastic coated inner handrim [2] Left☐ Right☐			$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)			NC
REAR WHEELS			Available		
RW01	☐	Mag [3]	20", 22" 24"		$150
RW02	☐	Newton One - spoke wheel [4]	20", 22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel	20", 22", 24", 25", 26"	*Client cost*	$390
RW04	☐	Spinergy Spox **Blk STD**☐ Wht [5]☐ Blue [5]☐ Red [5]☐ Yel [5]☐	22", 24", 25", 26"	*Client cost*	$915
RW05	☐	Spinergy LX **Blk STD**☐ Pur [5]☐ Wht [5]☐ Blue [5]☐ Red [5]☐ Yel [5]☐ Pink [5]☐ Grn [5]☐ Org [5]☐	22", 24", 25", 26"	*Client cost*	$1,135
RW07	☐	Spinergy CLX [6] **Blk STD**☐ Pur [5]☐ Wht [5]☐ Blue [5]☐ Red [5]☐ Yel [5]☐ Pink [5]☐ Grn [5]☐ Org [5]☐	24", 25"	*Client cost*	$1,690
TIRES ***(The brand and color may vary depending on parts availability)***			Available		
		AIRLESS			
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black	20", 22" 24"		STD
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black [7]	22", 24", 25", 26"	*Client cost*	$120
TI13	☐	Soft Urethane 1" - No tread, black - *Sentinel*	20", 22", 24", 25", 26"	*Client cost*	$150
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread [8]	22", 24", 26"	*Client cost*	$235
		PNEUMATIC			
TI05	☐	Pneumatic 1 3/8" - Medium Tread	20", 22", 24", 26"		NC
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - *SpeedLite*	20", 22", 24", 25", 26"	*Client cost*	$130
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - *TrailBlazer*	24", 25", 26"	*Client cost*	$210

1- *Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide ;*

2- *Only standard handrim position available;* **3-** *Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position;* **4-** *Not compatible with Shox tires;*

5- *color spokes require an additionnal 10 days lead-time;* **6-** *Only available with 1" tires;* **7-** *22" is only available with Mags, Newton One and Newton Gravity rear wheels;*

8- *Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels*

HANDRIMS

HANDRIM			Available		
HR01	☐	Aluminum anodized	20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized [2]	20", 22", 24", 25"	*Client cost*	$95
HR02	☐	Plastic coated [1]	20", 22", 24", 25", 26"		$180
HR03	☐	High friction coated [1,2, 3]	22", 24", 25", 26"	*Client cost*	$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip [2] *See page 6 for the color selection and price*	20", 22", 24", 25", 26"	*Client cost*	
HR28	☐	Surge [1,2,3] *(1 1/2" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"	*Client cost*	$485
HR05	☐	Surge LT [1,2,3] *(1 1/8" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"	*Client cost*	$515
HR04	☐	Natural Fit [1,2,3]	22", 24", 25", 26"	*Client cost*	$320
HR62	☐	Simi H [1,2,3] *(1 3/8" Oval handrim with high friction strip)*	22", 24", 25", 26"	*Client cost*	$485
HR63	☐	Nova H [1,2,3] *(1 1/4" Oval handrim with high friction strip)*	22", 24", 25", 26"	*Client cost*	$515
HR64	☐	Optimum H [1,2] *(1 1/8" heart-shaped handrim)*	22", 24", 25", 26"	*Client cost*	$320
HR99	☐	Omit handrim			NC

1- *Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position;* **2-** *Not available with Mag wheels;* **3** *- Not available with one arm drive;*

HANDRIMS OPTION

ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	*Full tab*			STD
HR030	☐	Narrow Handrim position no cut [4]	*Full tab*			NC
HR30	☐	Narrow Handrim position with cut [3]	*Short tab 2 positions*			$29
HR31	☐	Super Narrow Handrim [3]	*Short tab 2 positions STD* ☐ *Short tab 1 position* [4] ☐			$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip			***Client cost***	$200
HR27	☐	Super grip			***Client cost***	$230
HR099	☐	Omit thumb piece				NC


3- *Only available with Aluminum and Air Grip handrims;* ***4-*** *Not available with the combination of 24" Newton Gravity, Spinergy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;*

AXLE					
AX01	☐	Quick release axle [1]			STD
AX02	☐	Permanent axle			NC
AX03	☐	Quad release axle [1]		***Clinical justification***	$180
AXLE PLATE					
AX04	☐	Amputee adjustable axle plate [2]			$280

1- *Not required with one arm drive;* ***2*** *- May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down and the swing away anti-tipper*

REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM03	☐	3° camber			NC

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 1 1/4"

Most stable					**Most responsive**
AMPUTEE[1]	**3/4"**	**1" 1/4**	**1" 3/4**	**2" 1/4**	**2" 3/4**
☐	☐	☐	☐	☐	☐

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1- Only available with amputee axle plate

WHEEL LOCK					
WL08	☐	Newton Wheel lock - push to lock			STD
WL02	☐	Pull to lock			NC
WL04	☐	Newton grade aid push to lock *(works best with pneumatic or high pressure tire)*			$200
WL13 WL14	☐	6" removable extension handle [1]			$75
WL05	☐	Scissors Lock			$200
WL06	☐	Attendant Lock *(additional set of wheel lock mounted to the back)*			$240
WL15	☐	Unilateral wheel lock - pull to lock[2] Left ☐ Right ☐			$335
WL99	☐	Omit wheel locks			NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilties and restriction, refer to the table at page 7 and the "How-to documents" page of our website:
https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

FOOTREST

SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at*		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at*		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" [5]		Set at*		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" [5]		Set at*		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" [5]		Set at*		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) [1]	7.5" - 11"	Set at*		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" [1, 2]	9.5"- 13"	Set at*		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" [1, 2]	9.5" - 14.5"	Set at*		NC
FOOTREST OPTIONS						
FR09	☐	High mount *(footplate attached to the hanger with clamps)* [3]	5" - 13.5"	Set at*		$105
FR10	☐	Elevating legrest with calf pad, Length adjustable [4] Left ☐ Right ☐	14" - 22"	Set at:		$229
FR99	☐	Omit Footrest				NC

1- *Available only with two piece adjustable footplates;* ***2-*** *May interfere with front casters;* ***3-*** *Not available with front mount and the elevating legrest; * +/- 1/2"*

4- *The elevating legrest replace the footrest side selected; Not available with one piece footplate;* ***5-*** *Not available with 8" caster*

FOOTPLATES							
FP01	☐	Standard - Composite[1]					STD
FP02	☐	Newton adjustable angle - Composite					$155
FP03	☐	Newton adjustable angle - Aluminum					$175
FP06	☐	One piece adjustable angle flip-up - Aluminum					$215
FOOTPLATE ACCESSORIES							
FP08	☐	Length adjustable heel loop [2]					NC
FP09	☐	Calf strap					$62
FP30	☐	Bodypoint Aeromesh Padded calf strap				***Seating benefit***	$145

1- *Not available with front mount;* ***2-*** *not available with one piece footplate*

SWING AWAY RESIDUAL LIMB SUPPORT							
FR27	☐	Residual limb support - flat 6"x8"	Left ☐	Right ☐			$505
FR28	☐	Residual limb support - flat 8"x10"	Left ☐	Right ☐			$505
FR29	☐	Residual limb support - flat 10"x10"	Left ☐	Right ☐			$505
FR30	☐	Residual limb support - flat 10"x14"	Left ☐	Right ☐			$505
FR32	☐	Residual limb support - contour 6"x8"	Left ☐	Right ☐			$505
FR33	☐	Residual limb support - contour 8"x10"	Left ☐	Right ☐			$505

SEATING AND BACKREST

BACK UPHOLSTERY							
BH05	☐	Slip-on Nylon back upholstery					NC
BH06	☐	Adjustable tension Nylon back upholstery				***Seating benefit***	$265
BH99	☐	Omit Back upholstery					STD
SEAT CUSHION							
SC01	☐	2" Cushion				***Seating benefit***	$185
SC02	☐	3" Cushion				***Seating benefit***	$205
SC99	☐	Omit seat cushion					STD
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE *(push handle included - no further selection needed)*							
BC58-915	☐	Straight aluminum back cane	9" - 15"	*increment of 1/2"*	Set at[2]:		NC
BC58-1521	☐	Straight aluminum back cane	15" - 21"	*increment of 1/2"*	Set at[2]:		NC
BC60	☐	Straight aluminum back cane	21" - 24"	*increment of 1"*	Set at[2]:		$100
BC01-16	☐	8° bend aluminum back cane	16"				NC
BC01-18	☐	8° bend aluminum back cane	18"				NC
BC01-20	☐	8° bend aluminum back cane	20"				NC
BC59	☐	8° bend aluminum back cane	21" - 24"	*increment of 1"*	Set at[2]:		$100
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*							
BC11-912	☐	Straight aluminum back cane	9" - 12"	*increment of 1/2"*	Set at[2]:		$80
BC11-1215	☐	Straight aluminum back cane	12" - 15"	*increment of 1/2"*	Set at[2]:		$80
BC11-1518	☐	Straight aluminum back cane	15" - 18"	*increment of 1/2"*	Set at[2]:		$80
BC11-1821	☐	Straight aluminum back cane	18" - 21"	*increment of 1/2"*	Set at[2]:		$80
BC12-1315	☐	8° bend aluminum back cane	13.5" - 15"	*increment of 1/2"*	Set at[2]:		$85
BC12-1518	☐	8° bend aluminum back cane	15" - 18"	*increment of 1/2"*	Set at[2]:		$85
BC12-1821	☐	8° bend aluminum back cane	18" - 21"	*increment of 1/2"*	Set at[2]:		$85
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*							
BC13-1518	☐	Straight aluminum back cane	15" - 18"	*increment of 1/2"*	Set at[2]:		$275
BC13-1821	☐	Straight aluminum back cane	18" - 21"	*increment of 1/2"*	Set at[2]:		$275
BC14-1518	☐	8° bend aluminum back cane	15" - 18"	*increment of 1/2"*	Set at[2]:		$330
BC14-1821	☐	8° bend aluminum back cane	18" - 21"	*increment of 1/2"*	Set at[2]:		$330
		Select back cane angle:	85° [1] ☐	90° **STD** ☐	95° ☐	100° ☐	105° ☐ 110° ☐
PUSH HANDLES OPTIONS *(ONE option at a time)*							
PH01	☐	Push handles					STD
PH05	☐	Fold down push handles [3, 5]				***Client cost***	$165
PH06-920	☐	Clamp on height adjustable push handles 9"-20" [3,4, 5]				***Client cost***	$345
PH06-923	☐	Clamp on height adjustable push handles 9"-23" [3,4, 5]				***Client cost***	$345
PH02	☐	Omit push handles					NC
STROLLER HANDLE AND STABILIZER BAR *(only available with push handle)*							
PH04	☐	Stroller handle *(not available with the HD option)*				***Client cost***	$220
PH03	☐	Newton folding stabilizer bar *(not available with stroller handle)*					$275

1- *Not compatible with 8 degree bent back cane;* ***2-*** *If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery;* ***3-*** *265lb weight limit*

4- *Not compatible with swing away armrests, 8° bend back canes, 21-24" back canes and amputee axle plates; 5- Not compatible with stroller handler, stabilizer bar, headrest and attendant lock;*

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST *(Convertible to "T" armrest)* [1,4]								
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐			$405
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐			$405
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐			$405
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐			$405
HEIGHT ADJUSTABLE "T" ARMREST [1]								
AR13-610	☐	10" Desk length armpad [5]	6" - 10"					$385
AR13-812	☐	10" Desk length armpad	8" - 12"					$385
AR13-1014	☐	10" Desk length armpad	10" - 14"					$385
AR14-610	☐	14" Full length armpad [5]	6" - 10"					$385
AR14-812	☐	14" Full length armpad	8" - 12"					$385
AR14-1014	☐	14" Full length armpad	10" - 14"					$385
SWING AWAY ARMREST [2]								
AR05	☐	Tubular armrest - Short pad	8" - 11"					$325
AR23	☐	Tubular armrest - long pad	8" - 11"					$325
AR06	☐	10" Desk length flat armpad	8" - 11"					$325
AR17	☐	14" Full length flat armpad	8" - 11"					$325
ARMREST OPTIONS								
AR12	☐	14'' Full Length armpad - Gel Ovation (pair) [3]						$160
AR99	☐	Omit armrest						NC

1- *May have 1/2" adjustment restriction on both limits of the range;* ***2-*** *Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;*

3- *Only available with full length armrest; Not available with swing away and desk lenght armrests;* ***4-*** *Not available with attendant lock;* ***5-*** *May have some restriction with seat-to-floor height and rear wheel*

ACCESSORIES

SIDE GUARDS ***(not available with Flip Back and "T" Armrest)***								
SG01	☐	Plastic side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐				$235
SG02	☐	Carbon fiber side guards *(fixed)*					***Client cost***	$285
SG03	☐	Carbon fiber side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			***Client cost***	$320
SG06	☐	Plastic fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			***Client cost***	$340
SG04	☐	Carbon fiber fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			***Client cost***	$535
SG99	☐	Omit side guard						NC

1- *only compatible with 0° camber; Increase distance between fender side guards by 1 1/8" ;*

POSITIONING BELT								
SB01-S	☐	Velcro adjustable - **Small 35"** [1]					***Seating benefit***	$50
SB01-L	☐	Velcro adjustable - **Long 60"**					***Seating benefit***	$50
SB02-S	☐	Auto buckle - **Small 35"** [1]						$22
SB02-M	☐	Auto buckle - **Medium 48"**						$22
SB02-L	☐	Auto buckle - **Long 60"**						$22
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - **Medium 56"**					***Seating benefit***	$105
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - **Long 61"**					***Seating benefit***	$105
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - **Small 51"**					***Seating benefit***	$195
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - **Medium 56"**					***Seating benefit***	$195
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - **Long 61"**					***Seating benefit***	$195
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - **Small 51"**					***Seating benefit***	$220
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - **Medium 56"**					***Seating benefit***	$220
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - **Long 61"**					***Seating benefit***	$220
SB03	☐	Hardware attachment for 2 extra points (1 pair)					***Seating benefit***	$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - **Medium 18"**					***Seating benefit***	$365
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - **Long 24"**					***Seating benefit***	$365
SB99	☐	Omit positioning belt						STD

1- Not recommended for seat width of 18" and more

HEADREST								
AC01-S	☐	Small contoured head support pad (9"x5")					***Seating benefit***	$215
AC01-M	☐	Large rectangular head support pad (10"x5")					***Seating benefit***	$215
AC01-L	☐	XL rectangular head support pad (14"x5")					***Seating benefit***	$215
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar [1]					***Seating benefit***	$365

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES					
AC03	☐	Anti-tippers			$125
AC29	☐	Swing-away anti-tipper[1] Left ☐ Right ☐		Client cost	$190
AC27	☐	Tip assist [2]		Client cost	$65
AC06	☐	Cane and crutch holder Left ☐ Right ☐		Client cost	$115
AC07	☐	Oxygen tank holder [4] - D and E Type - *Installed to the left*		Client cost	$395
AC08	☐	Spoke guards [3] Available: 20", 22", 24", 26"		Client cost	$165
AC09	☐	Footrest neoprene impact guard *(pair)* - Black		Client cost	$90
AC10	☐	Seat pouch		Client cost	$105
AC30	☐	Newton multi-purpose pouch		Client cost	$80
AC26	☐	BodyPoint Mobility Bag		Client cost	$150
AC23	☐	Newton Backpack		Client cost	$135
AC12	☐	Velcro Closing Strap		Client cost	$60
AC13	☐	Tool kit		Client cost	$60

1- *Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate;* ***2-*** *Cannot be combined with Anti Tippers*


3- *Not available with Mag wheels and Spinergy;* ***4-*** *Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down*

Colors																				
Options [1]	Black	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	☐ $160	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	$350
Frame [2]	☐ STD			☐																NC
FC99	☐	Custom Color *Specify PPG Paintit color:*																		*additional 10 days lead-time* $422

1- *Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;*

2- *with black hydroformed aluminum crossbrace*

Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.25 to 2.75
22 in	16	0.75 to 2.25
	17+	0.75 to 2.75
24 in	16	0.75 to 1.75
	17	0.75 to 2.25
	18+	0.75 to 2.75

Note :