160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T.+1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



6.6

Helio

ORDER FORM - AADL #W156

1- Purchasing information

P.O. number:* Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

2- Dealer/Representative information

Dealer/Representative:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:*

Name:

Address:

City: Province:

Postal Code:

Phone:

** Required*

MA080 ☐ *English Manual (STD)* *MA081* ☐ *Manuel Français*

Taxes not included

MODEL									
WC03	☐	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (265 lb weight capacity)							**$4,495**
WC09	☐	**HD** Heavy Duty package *(350 lb weight capacity)*							$205
FRAME									
SEAT WIDTH [3]									
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" [1] ☐			NC
	☐	19" [1] ☐	20" [1] ☐	21" [2] ☐	22" [2] ☐				$365
SEAT DEPTH									
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐			NC
	☐	19" ☐	20" ☐						$280
FRAME LENGTH [5]									
SP01	☐	Standard frame length configuration							STD
SP02	☐	Frame length 1" longer than standard configuration							NC
SP03	☐	Frame length 2" longer than standard configuration							NC
SP04	☐	Frame length 1" shorter than standard configuration							NC
SP05	☐	Frame length 2" shorter than standard configuration							NC
TRANSIT TIE DOWN									
AC04	☐	Transit Tie Down *(Unoccupied)*							NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt [4]	Short belt 20"-41"						$370
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt [4]	Long belt 20"-46"						$370
AC35	☐	WC-19 Transit Tie Down without belt [4]							$345

1- *Available with HD option;* ***2-*** *Only available with HD option*

4- *May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option*

5 *- See chart on page 7*

3. Width calculator

PACKAGES					
PK01	☐	Newton Ultralight Package *(include: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels)*		**Save $245**	$540

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13" ☐	14" ☐	15" ☐	13 1/2" ☐	14 1/2" ☐	15 1/2" ☐												
	P2	13 1/2" ☐	14 1/2" ☐	15 1/2" ☐	14" ☐	15" ☐	16" ☐	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐									
5"	P1																		
	P2				14" ☐	15" ☐	16" ☐	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐									
	P3	14" ☐	15" ☐	16" ☐	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐	15" ☐	16" ☐	17" ☐	15 1/2" ☐	16 1/2" ☐	17 1/2" ☐						
	P4	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐	15" ☐	16" ☐	17" ☐	15 1/2" ☐	16 1/2" ☐	17 1/2" ☐	16" ☐	17" ☐	18" ☐	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐			
7"	P1										16 1/2" ☐	17 1/2" ☐	18 1/2" ☐	17" ☐	18" ☐	19" ☐			
	P2							16 1/2" ☐	17 1/2" ☐	18 1/2" ☐	17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐			
	P3							17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐	18" ☐	19" ☐	20" ☐	18 1/2" ☐	19 1/2" ☐	20 1/2" ☐
	P4	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐	17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐	18" ☐	19" ☐	20" ☐	18 1/2" ☐	19 1/2" ☐	20 1/2" ☐	19" ☐	20" ☐	21" ☐
	P5	17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐	18" ☐	19" ☐	20" ☐	18 1/2" ☐	19 1/2" ☐	20 1/2" ☐	19" ☐	20" ☐	21" ☐	19 1/2" ☐	20 1/2" ☐	21 1/2" ☐
Froglegs	P1				14 1/2" ☐	15 1/2" ☐	16 1/2" ☐												
	P2	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐	15" ☐	16" ☐	17" ☐	15 1/2" ☐	16 1/2" ☐	17 1/2" ☐									
	P3	15" ☐	16" ☐	17" ☐	15 1/2" ☐	16 1/2" ☐	17 1/2" ☐	16" ☐	17" ☐	18" ☐	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐						

Height available +- 1/4" (see stem bolt and fork combination at the end of the order form)

FRONT CASTERS

CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		$170
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		$120
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		$198
CASTER OPTIONS						
CA12	☐	Caster Pin Locks *(Available with +1" or +2" Stem bolt ONLY)* [1]				$275
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				$525

1- *Not available with WC-19*

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available			
☐	20" (451) wheel	12 1/2	to	17 1/4	**Specify desired height:** _____
☐	22" (501) wheel	13 1/2	to	18 1/4	
☐	24" (540) wheel	14 1/2	to	19 1/4	
☐	25" (559) wheel	15	to	19 3/4	
☐	26" (590) wheel	15 1/2	to	20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

REAR WHEELS

REAR WHEELS			Available		
RW01	☐	Mag [1]	20", 22" 24"		$150
RW02	☐	Newton One - spoke wheel [2]	20", 22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel	20", 22", 24", 25", 26"		$390
RW04	☐	Spinergy Spox **Blk STD** ☐ Wht [3] ☐ Blue [3] ☐ Red [3] ☐ Yel [3] ☐	22", 24", 25", 26"		$915
RW05	☐	Spinergy LX **Blk STD** ☐ Pur [3] ☐ Wht [3] ☐ Blue [3] ☐ Red [3] ☐ Yel [3] ☐ Pink [3] ☐ Grn [3] ☐ Org [3] ☐	22", 24", 25", 26"		$1,135
TIRES ***(The brand and color may vary depending on parts availability)***			**Available**		
		AIRLESS			
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black [6]	20", 22" 24"		STD
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black [4,6]	22", 24", 25", 26"		$120
TI13	☐	Soft Urethane 1" - No tread, black - *Sentinel*	20",22", 24", 25", 26"		$150
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread [5]	22", 24", 26"		$235
		PNEUMATIC			
TI05	☐	Pneumatic 1 3/8" - Medium Tread [6]	20", 22", 24", 26"		NC
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - *SpeedLite*	20", 22", 24", 25", 26"		$130
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - *TrailBlazer*	24", 25", 26"		$210

1- *Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position;* ***2-*** *Not compatible with Shox tires;* ***3-*** *color spokes require an additionnal 10 days lead-time;*

4- *22" is only available with Mags, Newton One and Newton Gravity wheels*

5 *- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels* ***6-*** *Not available with flip-back armrest when the 6° camber is selected;*

HANDRIMS

HANDRIM			Available		
HR01	☐	Aluminum anodized	20", 22", 24", 25", 26"		STD
HR02	☐	Plastic coated [1]	20", 22", 24", 25", 26"		$180
HR03	☐	High friction coated [1,2,3]	22", 24", 25", 26"		$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip [2] *See on page 6 for color choice and price*	20", 22", 24", 25", 26"		
HR28	☐	Surge [1,2,3,4] *(1 1/2" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"		$485
HR05	☐	Surge LT [1,2,3,4] *(1 1/8" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"		$515
HR04	☐	Natural Fit [1,2,3,4]	22", 24", 25", 26"		$320
HR99	☐	Omit handrim			NC

1- *Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position;* ***2-*** *Not available with Mag wheels;* ***3*** *- Not available with one arm drive;*
4 *- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece*

HANDRIMS OPTION

ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	*Full tab*			STD
HR030	☐	Narrow Handrim position no cut [6]	*Full tab*			NC
HR30	☐	Narrow Handrim position with cut [5]	*Short tab 2 positions*			$29
HR31	☐	Super Narrow Handrim [5]	*Short tab 2 positions STD* ☐ *Short tab 1 position* [6] ☐			$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				$200
HR27	☐	Super grip				$230
HR099	☐	Omit thumb piece				NC


5- *Only available with Aluminum and Air Grip handrims; 6- Not available with the combination of 24" Newton Gravity, Spinergy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;*

AXLE [1]						
AX01	☐	Quick release axle				STD
AX02	☐	Permanent axle				NC
AX03	☐	Quad release axle				$180
AXLE PLATE						
AX04	☐	Amputee adjustable axle plate [2]				$280

1- *Not required with one arm drive;* ***2*** *- May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down and the swing-away anti-tipper*

REAR WHEEL CAMBER						
CM00	☐	0° camber				STD
CM03	☐	3° camber				NC
CM06	☐	6° camber *(not available with amputee axle plate)*				$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =

Most stable														
AMPUTEE[1] ☐	**1"** ☐	1/4" ☐	1/2" ☐	3/4" ☐	**2"** ☐	1/4" ☐	1/2" ☐	3/4" ☐	**3"** ☐	1/4" ☐	1/2" ☐	3/4" ☐	**4"** ☐

If not indicated, default setting is at 1.5''

Possible conflict between a seat depth under 16'' and a center of gravity over 2.5''.

Possible conflict between a 26'' rear wheel and a center of gravity over 2.5''.

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK						
WL08	☐	Newton Wheel lock - push to lock				STD
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock *(works best with pneumatic or high pressure tire)*				$200
WL13 WL14	☐	6" removable extension handle				$75
WL05	☐	Scissors Lock				$200
WL06	☐	Attendant Lock *(additional set of wheel lock mounted to the back)*				$240
WL99	☐	Omit wheel locks				NC

FOOTREST

SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" [5]		Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" [5]		Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" [5]		Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) [1]	7.5" - 11"	Set at*:		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" [1, 2]	9.5"- 13"	Set at*:		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" [1, 2]	9.5" - 14.5"	Set at*:		NC
FOOTREST OPTIONS [3]						
FR09	☐	High mount *(footplate attached to the hanger with clamps)* [3]	5" - 13.5"	Set at*:		$105
FR10	☐	Elevating legrest with calf pad, Length adjustable [4] Left ☐ Right ☐	14" - 22"	Set at:		$229
FR99	☐	Omit Footrest				NC

***1-** Available only with two piece adjustable footplates; **2-** May interfere with front casters; **3-** Not available with front mount and the elevating legrest; * +/- 1/2"*
***4-** The elevating legrest replace the footrest side selected; Not available with one piece footplate; **5-** Not available with 8" caster*

FOOTPLATES						
FP01	☐	Standard - Composite [1]				STD
FP02	☐	Newton adjustable angle - Composite				$155
FP03	☐	Newton adjustable angle - Aluminum				$175
FP04	☐	Newton adjustable angle - Carbon fiber *(with rubber inserts)*				$395
FP06	☐	One piece adjustable angle flip-up - Aluminum				$215
FP07	☐	One piece adjustable angle flip-up - Carbon fiber *(with rubber inserts)*				$480
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop [2]				NC
FP09	☐	Calf strap				$62
FP10	☐	Standard Padded Calf strap				$80
FP30	☐	Bodypoint Aeromesh Padded calf strap				$145

***1-** Not available with front mount; **2-** not available with one piece footplate*

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 6"x8"	Left ☐	Right ☐		$505
FR28	☐	Residual limb support - flat 8"x10"	Left ☐	Right ☐		$505
FR29	☐	Residual limb support - flat 10"x10"	Left ☐	Right ☐		$505
FR30	☐	Residual limb support - flat 10"x14"	Left ☐	Right ☐		$505
FR32	☐	Residual limb support - contour 6"x8"	Left ☐	Right ☐		$505
FR33	☐	Residual limb support - contour 8"x10"	Left ☐	Right ☐		$505

SEATING AND BACKREST

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery				$265
BH99	☐	Omit Back upholstery				STD
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion				$185
SC02	☐	3" Cushion				$205
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*						
BC06-912	☐	Straight Carbon fiber back cane	9" - 12"	*increment of 1/2"*	Set at [2]	NC
BC06-1215	☐	Straight Carbon fiber back cane	12" - 15"	*increment of 1/2"*	Set at [2]	NC
BC06-1518	☐	Straight Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at [2]	NC
BC06-1821	☐	Straight Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at [2]	NC
BC08-1315	☐	8° Bend Carbon fiber back cane	13.5" - 15"	*increment of 1/2"*	Set at [2]	$170
BC08-1518	☐	8° Bend Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at [2]	$170
BC08-1821	☐	8° Bend Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at [2]	$170
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*						
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	*increment of 1/2"*	Set at [2]	$330
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at [2]	$330
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at [2]	$330
BC10-1315	☐	8° bend Carbon fiber back cane	13.5" - 15"	*increment of 1/2"*	Set at [2]	$380
BC10-1518	☐	8° bend Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at [2]	$380
BC10-1821	☐	8° bend Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at [2]	$380
		Select back cane angle:	85° [1] ☐	90° ☐	95° **STD** ☐	100° ☐	105° ☐	110° ☐

PUSH HANDLES OPTIONS *(ONE option at a time)*							
PH01	☐	Push handles					STD
PH05	☐	Fold down push handles [3, 5]					$165
PH06-920	☐	Clamp on height adjustable push handles 9"-20" [3,4, 5]					$345
PH06-923	☐	Clamp on height adjustable push handles 9"-23" [3,4, 5]					$345
PH02	☐	Omit push handles					NC
STROLLER HANDLE AND STABILIZER BAR *(only available with push handle)*							
PH04	☐	Stroller handle *(not available with the HD option)*					$220
PH03	☐	Newton folding stabilizer bar *(not available with stroller handle)*					$275

1- *Not available with 8 degree bend cane;* ***2-*** *If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery;* ***3-*** *265 lb weight limit;*
4- *Not compatible with swing away armrests, 8° bend back canes, 9-12" back canes and amputee axle plates;* ***5-*** *Not compatible with stroller handle, stabilizer bar, headrest and attendant lock*

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST *(Convertible to "T" armrest)* [1, 4]							
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		$405
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		$405
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		$405
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		$405
HEIGHT ADJUSTABLE "T" ARMREST [1]							
AR13-610	☐	10" Desk length armpad [5]	6" - 10"				$385
AR03-812	☐	10" Desk length armpad	8" - 12"				$385
AR03-1014	☐	10" Desk length armpad	10" - 14"				$385
AR14-610	☐	14" Full length armpad [5]	6" - 10"				$385
AR04-812	☐	14" Full length armpad	8" - 12"				$385
AR04-1014	☐	14" Full length armpad	10" - 14"				$385
HEIGHT ADJUSTABLE SWING AWAY ARMREST [2]							
AR05	☐	Tubular armrest - short pad	8" - 11"				$325
AR23	☐	Tubular armrest - long pad	8" - 11"				$325
AR06	☐	10" Desk length flat armpad	8" - 11"				$325
AR17	☐	14" Full length flat armpad	8" - 11"				$325
ARMRESTS OPTIONS							
AR12	☐	14'' Full Length armpad - Gel Ovation (pair) [3]					$160
AR99	☐	Omit armrest					NC

1- *May have 1/2" adjustment restriction on both limits of the range;* ***2-*** *Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;*
3- *Only available with full length armrest; Not available with swing away and desk lenght armrests;* ***4-*** *Not available with attendant lock;* ***5-*** *May have some restriction with seat-to-floor height and rear wheel*

ACCESSORIES

SIDE GUARDS *(not available with Flip Back and "T" Armrest)*							
SG01	☐	Plastic side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$235
SG02	☐	Carbon fiber side guards *(fixed)*					$285
SG03	☐	Carbon fiber side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$320
SG06	☐	Plastic fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$340
SG04	☐	Carbon fiber fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$535
SG99	☐	Omit side guard					NC

1- *only compatible with 0° camber;Increase width by 1/2" per side for a total of 1" ;*

POSITIONING BELT							
SB01-S	☐	Velcro adjustable - **Small 35"** [1]					$50
SB01-L	☐	Velcro adjustable - **Long 60"**					$50
SB02-S	☐	Auto buckle - **Small 35"** [1]					$22
SB02-M	☐	Auto buckle - **Medium 48"**					$22
SB02-L	☐	Auto buckle - **Long 60"**					$22
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - **Medium 56"**					$105
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - **Long 61"**					$105
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - **Small 51"**					$195
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - **Medium 56"**					$195
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - **Long 61"**					$195
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - **Small 51"**					$220
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - **Medium 56"**					$220
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - **Long 61"**					$220
SB03	☐	Hardware attachment for 2 extra points (1 pair)					$15
SB99	☐	Omit positioning belt					STD

1- Not recommended for seat width of 18" and more

HEADREST [1]					
AC01-S	☐	Small contoured head support pad (9"x5")			$215
AC01-M	☐	Large rectangular head support pad (10"x5")			$215
AC01-L	☐	XL rectangular head support pad (14"x5")			$215
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar [1]			$365

***1-** Only available with push handles - Not available with a stroller handle or a folding stabilizer bar*

ACCESSORIES					
WC10	☐	E-Motion adaptor Kit *(2 adaptors) - Only necessary with 6° camber*			$35
AC03	☐	Anti-tippers			$125
AC06	☐	Cane and crutch holder Left ☐ Right ☐			$115
AC07	☐	Oxygen tank holder [1] - D and E Type - installed on the left			$395
AC08	☐	Spoke guards [2] Available: 20", 22", 24", 25", 26"			$165
AC09	☐	Footrest neoprene impact guard *(pair)* - Black			$90
AC10	☐	Seat pouch			$105
AC23	☐	Newton Backpack			$135
AC12	☐	Velcro Closing Strap			$60
AC13	☐	Tool kit			$60

1 - Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down ***2-** Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel ;*


Colors

Options [1]	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	☐ $160		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	$350
Frame[2]		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC
FC99	☐	Custom Color	Specify PPG Paintit color:																*additional 10 days lead-time*		$422

***1-** Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;*

***2-** With carbon fiber crossbrace;*

Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
22 in	16	1.0 to 2.75
	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Frame length chart

	Seat Depth						
Frame length	**14"**	**15"**	**16"**	**17"**	**18"**	**19"**	**20"**
STD							
+1							*N/A*
+2						*N/A*	*N/A*
-1	*N/A*						
-2	*N/A*	*N/A*					

Notes :