MOTION COMPOSITES
160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T.+1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



Move MVQS

Quick Ship Order Form

Ontario 6.5

P.O. number:* | Ship to:*

Company Name:* | Name:

Purchasing contact:* Phone:* | Address:

Clinician/Sales rep: | City: Province

Email: | Postal Code:

Marked for: | Phone:

** Required*

Taxes not included

MOVE

Ultralight MOVE wheelchair aluminum	Basalt Black	WM0001196		**$1,325**	☐
Included: seating upholstery / **Excluded** back upholstery and cushion	265 lbs. weight capacity				

Front Seat-to-Floor Height Configuration

13"					14"					16 ½"				
M6	16x16		NC	☐	**M1**	18x16		NC	☐	**M11**	16x18		NC	☐
M7	16x18		NC	☐	**M2**	18x18		NC	☐	**M12**	18x18		NC	☐
M9	18x18		NC	☐	**M3**	18x20	WAMB	$280	☐	**M13**	18x20	WAMB	$280	☐
M10	20x18	WAMA	$320	☐	**M4**	20x20	WAMA WAMB	$600	☐	**M14**	20x18	WAMA	$320	☐
WHEELS Front, Rear & Handrim					**WHEELS** Front, Rear & Handrim					**WHEELS** Front, Rear & Handrim				
4"x1" Composites caster w/PU tire			NC	■	5"x1" Composites caster w/PU tire			NC	■	6"x1" Composites caster w/PU tire			NC	■
5" fork + Std stem bolt			NC	■	5" fork + Std stem bolt			NC	■	7" fork + Std stem bolt			NC	■
22" Newton One rear wheel			NC	■	22" Newton One rear wheel			NC	■	24" Newton One rear wheel			NC	■
Soft urethane tire $1^{3/8}$"			NC	■	Soft urethane tire $1^{3/8}$"			NC	■	Soft urethane tire $1^{3/8}$"			NC	■
Aluminum anodized handrim			NC	■	Aluminum anodized handrim			NC	■	Aluminum anodized handrim			NC	■
SWING-AWAY FOOTREST					**SWING-AWAY FOOTREST**					**SWING-AWAY FOOTREST**				
60° Swing-away footrests			$270	■	70° Swing-away footrests			$270	■	70° Swing-away footrests			$270	■
Newton Composite Adjustable angle footplate w/ heel loop		WAMJ	$181	■	Newton Composite Adjustable angle footplate w/ heel loop		WAMJ	$181	■	Newton Composite Adjustable angle footplate w/ heel loop		WAMJ	$181	■
This configuration:					**This configuration:**					**This configuration:**				
Front STF Height	Rear STF Height				Front STF Height	Rear STF Height				Front STF Height	Rear STF Height			
13" - 14"	**13" - 18"**				**13 ¾" - 14 ¾"**	**13" - 18"**				**15 ½" - 17 ½"**	**14" - 19"**			

Possible configurations:

Front STF Height **12 ½" - 19 ½"** Rear STF Height **12 ¼" - 20"**

PARTS & ACCESSORIES on all configurations

8° BEND BACK (with fixed height push handles)				
Aluminum back cane	20" height		NC	■
FLIP BACK ARMREST (Convertible to "T" armrest)				
14" Full Length, 8 - 12" height adjustable			$405	■
AXLE (with 0° camber)				
Quick release axle		WAMU	$100	■
WHEEL LOCK				
Push to lock with 6" removable extension handle			$75	■
POSITIONING BELT				
Additional attachment hardware for 4 point seatbelt			$15	■
ACCESSORIES				
Anti tippers			$125	■

Quick Ship Parts and Accessories Program

Motion Composites' Quick Ship Parts and Accessories Program adds flexibility to the convenience of our Quick Ship wheelchairs

These specially discounted parts will be added to your Quick Ship wheelchair order for you to install and customize upon receipt

Invoice code						
QP0004	Aluminum straight back cane 18-21"			PR	$120	☐
QP0005	Aluminum back cane 8° bend 18-21"			PR	$140	☐
QP0006	Angle adjustable - Aluminium straight back cane 18-21"			PR	$240	☐
QP0098	Push handles			PR	$110	☐
QP0134	Slip-On Nylon Back UPH, 16x18-21in			EA	$67	☐
QP0135	Slip-On Nylon Back UPH, 18x18-21in			EA	$67	☐
QP0136	Slip-On Nylon Back UPH, 20x18-21in			EA	$67	☐
QP0137	Adj Tension Nylon Back UPH, 16x18-21in			EA	$265	☐
QP0138	Adj Tension Nylon Back UPH, 18x18-21in			EA	$265	☐
QP0139	Adj Tension Nylon Back UPH, 20x18-21in			EA	$265	☐
QP0167	Swing Out residual limb support flat cushion - L	6''x8''	*E1020*	EA	$505	☐
QP0168	Swing Out residual limb support flat cushion - R	6''x8''	*E1020*	EA	$505	☐
QP0169	Swing Out residual limb support flat cushion - L	8''x10''	*E1020*	EA	$505	☐
QP0170	Swing Out residual limb support flat cushion - R	8''x10''	*E1020*	EA	$505	☐
QP0171	Swing Out residual limb support flat cushion - L	10"x10	*E1020*	EA	$505	☐
QP0172	Swing Out residual limb support flat cushion - R	10"x10	*E1020*	EA	$505	☐
QP0173	Swing Out residual limb support flat cushion - L	10"x14	*E1020*	EA	$505	☐
QP0174	Swing Out residual limb support flat cushion - R	10"x14	*E1020*	EA	$505	☐
QP0175	Swing Out residual limb support contour cushion - L	6"x8"	*E1020*	EA	$505	☐
QP0176	Swing Out residual limb support contour cushion - R	6"x8"	*E1020*	EA	$505	☐
QP0177	Swing Out residual limb support contour cushion - L	8"x10"	*E1020*	EA	$505	☐
QP0178	Swing Out residual limb support contour cushion - R	8"x10"	*E1020*	EA	$505	☐
QP0149	Swing Away ELR/calf pad/std FP, Left & heel loop (14-16in)			EA	$410	☐
QP0150	Swing AwayELR/calf pad/std FP, Left & heel loop (17-18in)			EA	$410	☐
QP0151	SwingAway ELR/calf pad/std FP, Left & heel loop (19in +)			EA	$410	☐
QP0152	Swing Away ELR/calf pad/std FP, Right & heel loop (14-16in)			EA	$410	☐
QP0153	Swing Away ELR/calf pad/std FP, Right & heel loop (17-18in)			EA	$410	☐
QP0154	Swing Away ELR/calf pad/std FP, Right & heel loop (19in +)			EA	$410	☐
QP0013	Calf strap (wheelchair width: 15 to 17 in)			EA	$62	☐
QP0014	Calf strap (wheelchar width: 18 to 20 in)			EA	$62	☐
QP0125	Bodypoint Aeromesh Padded calf strap (wheelchair Width: 12 to 16 in)			EA	$145	☐
QP0126	Bodypoint Aeromesh Padded calf strap (wheelchar width: 17 to 20 in)			EA	$145	☐
QP0127	Bodypoint Aeromesh Padded calf strap (wheelchar width: 21 to 24 in)			EA	$145	☐
QP0017	Composite caster 3" x 1"			PR	$85	☐
QP0018	Composite caster 4" x 1"			PR	$85	☐
QP0019	Composite caster 5" x 1"			PR	$100	☐
QP0020	Composite caster 6" x 1"			PR	$100	☐
QP0021	Composite caster 7" x 1"			PR	$110	☐
QP0022	Composite caster 8" x 1"			PR	$110	☐
QP0077	Newton Ultracaster Composite caster 4" x 1½" Soft-Roll			PR	$120	☐
QP0078	Newton Ultracaster Composite caster 5" x 1½" Soft-Roll			PR	$120	☐
QP0079	Newton Ultracaster Composite caster 6" x 1½" Soft-Roll			PR	$120	☐
QP0080	Newton Ultracaster Aluminum caster 4" x 1½" Soft-Roll			PR	$198	☐
QP0081	Newton Ultracaster Aluminum caster 5" x 1½" Soft-Roll			PR	$198	☐
QP0082	Newton Ultracaster Aluminum caster 6" x 1½" Soft-Roll			PR	$198	☐
QP0031	Pneumatic caster 6" x 1¼"			PR	$170	☐
QP0032	Pneumatic caster 8" x 1¼"			PR	$170	☐

Invoice code				
QP0033	Standard stem bolt (anti-flutter)	PR	$21	☐
QP0034	+1" stem bolt (anti-flutter)	PR	$21	☐
QP0035	+2" stem bolt (anti-flutter)	PR	$21	☐
QP0076	Length adjustable heel loop	PR	$41	☐
QP0096	Newton Airgrip 22'' black handrim	PR	$160	☐
QP0097	Newton Airgrip 24'' black handrim	PR	$160	☐
QP0036	4" fork (2 holes)	PR	$62	☐
QP0037	5" fork (4 holes)	PR	$62	☐
QP0038	7" fork (5 holes)	PR	$62	☐
QP0114	Newton One wheels w/Urethane tire & Alu Pushrim - 22"	PR	$450	☐
QP0115	Newton One wheels w/Urethane tire & Alu Pushrim - 24"	PR	$450	☐
QP0039	Mag wheels w/Urethane tire & Alu Pushrim - 20"	PR	$450	☐
QP0040	Mag wheels w/Urethane tire & Alu Pushrim - 22"	PR	$450	☐
QP0041	Mag wheels w/Urethane tire & Alu Pushrim - 24"	PR	$450	☐
QP0044	Wheel lock - pull to lock	PR	$105	☐
QP0045	Wheel lock - pull to lock with removable extension handle	PR	$175	☐
QP0046	Wheel lock - Newton grade aid push to lock	PR	$305	☐
QP0048	Positioning belt - velcro adjustable - 35"	EA	$50	☐
QP0049	Positioning belt - auto buckle - 35"	EA	$80	☐
QP0195	Positioning belt - auto buckle - 48"	EA	$80	☐
QP0050	Bodypoint 2 point Non Padded - medium (56")	EA	$105	☐
QP0051	Bodypoint 2 point Non Padded - large (61")	EA	$105	☐
QP0052	Bodypoint 2 point Padded - small (51")	EA	$195	☐
QP0053	Bodypoint 2 point Padded - medium (56")	EA	$195	☐
QP0054	Bodypoint 2 point Padded - large (61")	EA	$195	☐
QP0055	Bodypoint 4 point Padded - small (51")	EA	$220	☐
QP0056	Bodypoint 4 point Padded - medium (56")	EA	$220	☐
QP0057	Bodypoint 4 point Padded - large (61")	EA	$220	☐
QP0058	Cane carrier - Left	EA	$115	☐
QP0059	Cane carrier - Right	EA	$115	☐
QP0060	Spoke guards - 22"	PR	$165	☐
QP0061	Spoke guards - 24"	PR	$165	☐
QP0062	Footrest neoprene impact guard (pr) - Black	PR	$90	☐
QP0063	Seat pouch	EA	$105	☐
QP0103	Newton Backpack	EA	$135	☐
QP0071	Tool Kit	EA	$90	☐
QP0192	Unoccupied transit tie-down	pr	$120	☐

Table 1

SLIP ON BACK UPHOLSTERY		LENGTH inches (cm)
		18-21
		(45.7-53.3)
Width Inches (cm)	16 (40.6)	9001832
	18 (45.7)	9001842
	20 (50.8)	9001852

Table 2

ADJUSTABLE TENSION BACK UPHOLSTERY		LENGTH inches (cm)
		18-21
		(45.7-53.3)
Width Inches (cm)	16 (40.6)	9001777
	18 (45.7)	9001787
	20 (50.8)	9001797

Move REAR Seat-to-floor Height

20"	22"	24"	25"	26"
12 1/4"	13"	14"	14 1/4"	15"
13 1/4"	14"	15"	15 1/4"	16"
14 1/4"	15"	16"	16 1/4"	17"
15 1/4"	16"	17"	17 1/4"	18"
16 1/4"	17"	18"	18 1/4"	19"
17 1/4"	18"	19"	19 1/4"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Move FRONT Seat-to-floor Height

		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
Fork/Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15		14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16		15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17		16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +- 1/4"

Notes :