RAZ-AT ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date:_____

Dealer / Location: _____ T:_____ F:_____

PO# / Client:_____ Contact:_____ Email: _____

Raz-AT – Attendant Tilt Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height-adjustable stainless-steel frame, Fixed Back, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. User weight capacity 407 lbs.
For help with justification, please see our Evaluation and Medical Justification Form.

Need help with this order form? Call customer service at (866) 875-2877 or (519) 942 8441

1. Standard Raz-AT Frame (18" W) – If you require a special width, depth and / or height see [Sec. 2]

☐	Z30018	Raz-AT w/ 18" W base frame – Standard base frame width, will fit over most toilets	$3,825

2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.7]

☐	Z30016	Raz-AT w/ 16" W frame width[1] (base chair price – custom upcharge below will be added)	$3,825
	Z3SW16	Upcharge for custom 16" W seat and base frame – will not fit over most toilets	$690
☐	Z30020	Raz-AT w/ 20" W frame width[2] (base chair price – custom upcharge below will be added)	$3,825
	Z3SW20	Upcharge for custom 20" W seat and base frame	$690
☐	Z30022	Raz-AT w/ 22" W frame width (base chair price – custom upcharge below will be added)	$3,825
	Z3SW22	Upcharge for custom 22" W seat and base frame	$805
☐	Z30024	Raz-AT w/ 24" W frame width (base chair price – custom upcharge below will be added)	$3,825
	Z3SW24	Upcharge for custom 24" W seat and base frame	$910
Custom Frame Configurations (Upcharges below can be added to all available frame widths)			
☐	Z3LNGSF3	Upcharge for custom 3” Longer Seat Frame[3] (Required for Custom Seat depths ≥ 21")	$690
☐	Z3LNGBF3	Upcharge for custom 3" Longer Base Frame[4] (Recommended for Custom Seat depths ≥ 21")	$690
☐	Z3LOWF5	Upcharge for Super-Low Frame[5] (16"-20" STF, 5" lower than standard) Note: Tilt is limited to 15°	$690

[1] 18" W Visco-Foam Interface seats can be modified to fit.

[2] 18" W seats can be selected but must also select Seat Support Kit ZSSKT [Sec.3]

[3] Seat frame is lengthened forward of the tilt pivot in order to provide better support for the front of the seat

[4] Base frame is lengthened forward of the tilt pivot in order to reduce front caster loading

[5] Forward Offset Foot Support Receiver ZFOR [Sec. 11] is required when selecting seats that are 20" or deeper.

3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 16" D and 19" D, are now designated 17" D and 20" D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Molded Seats (thickness ~2.5”-3”))		
☐ ZMS17FRO	18" W x 17" D – Front Access	Standard
☐ ZMS17BRI	18" W x 17" D – Front Bridged	N/C
☐ ZMS20FRO	18" W x 20" D – Front Access	$155
☐ ZMS20BRI	18" W x 20" D – Front Bridged	$155
Specialty Seats		
☐ ZSOM	Shower-Only Seat – Molded 18" W x 16”-21” D	$460
☐ ZEJRM	E&J Replica Seat – Molded – 18" W x 17" D	$225
☐ ZEJRV	E&J Replica Seat – Visco Foam – 18" W x 17" D	$970

Visco-Foam Interface Seats (thickness ~3.5")		
☐ ZVFI17BRI	18" W x 17" D – Front Bridged	$495
☐ ZVFI17FRO	18" W x 17" D – Front Access	$495
☐ ZVFI17LTO	18" W x 17" D – Left Access	$495
☐ ZVFI17RTO	18" W x 17" D – Right Access	$495
☐ ZVFI17RRO	18" W x 17" D – Rear Access	$495
☐ ZVFI20BRI	18" W x 20" D – Front Bridged	$570
☐ ZVFI20FRO	18" W x 20" D – Front Access	$570
☐ ZVFI20LTO	18" W x 20" D – Left Access	$570
☐ ZVFI20RTO	18" W x 20" D – Right Access	$570
☐ ZVFI20RRO	18" W x 20" D – Rear Access	$570

Seat Support Kit - These are used to support the sides of seats that are narrower than seat frames on which they are mounted. A typical example would be when a 20" W chair is needed but an 18" W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.

☐ ZSSKT	Seat Support Kit for Raz-AT	$165
☐ ZR00092	4" Link	$100
☐ ZCM0456	Green Adjustment Knob	$10

Headrest Post is 6" long. For longer tubes, order additional links (ZR00092 + ZCM0456)

4. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. See Custom Seat Assessment Form – Resources tab at razdesigninc.com

☐ ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below) Select if both width and depth are ≤ 21"	$1190
☐ ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	$1360
☐ ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if width and/or depth is ≥ 22"	$1,545

Custom Specifications (for above custom seats)

Seat Frame Width:	_ _	16” / 18” / 20" / 22" / 24"
Custom Seat Width:	_ _	16” / 18” / 20” / 22" / 24" (cannot be more than 2” wider than seat frame width)
Custom Seat Depth:	_ _	15” / 16” / 17” / 18” / 19” / 20” / 21” / 22" / 23" / 24" [≥ 21” requires Longer Frame (Z3LNGSF3)]
Access Opening:	_ _	FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
Access Opening Width:	_ _	2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)
Aperture Width:	_ _	3.0"-8.0" (6.5” recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

5. Fixed Back Frames (18" W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded, vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (90°/97°) with cane-style posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

☐ ZFB1818	Fixed Back – 18”W x 18”H	Standard
☐ ZFBU1818	Tension-Adjustable Back Upholstery for Fixed Back	$180

6. Adjusta-Back Frames – Back frame widths will be the same as the seat frame widths unless otherwise selected in Sec 8.

Includes height-and angle-adjustable arms, tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (90°/97°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. Adjusta-Back – 24"H includes straight posts, which do not allow the back upholstery to be raised. Includes the Armrest Lock Kit (ZALK).

☐ ZAB _ _ 15*	Adjusta-Back[1] – 15" H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	$285
☐ ZAB _ _ 17*	Adjusta-Back[1] – 17" H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	$285
☐ ZAB _ _ 20*	Adjusta-Back[1] – 20" H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	$370
☐ ZAB _ _ 24*	Adjusta-Back[1] – 24" H * insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	$420
☐ ZDTL	Dual Triggers – both on left side	N/C
☐ ZDTR	Dual Triggers – both on right side	N/C

Reclining Backs – 30° (Dual triggers are standard when recline option is added) Reclining back frame width must be the same as the seat frame width. 2-Point Arms (Z2PLT / Z2PST) are included with Reclining Backs.

☐ ZRBT _ _ 20*	Reclining Adjusta-Back with 20" H back * insert width to complete the code (16" / 18" / 20" / 22" / 24")	$1,320
☐ ZRBT _ _ 24*	Reclining Adjusta-Back with 24" H back * insert width to complete the code (16" / 18" / 20" / 22" / 24")	$1,385

[1] Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation.

7. Custom Back Frame Widths Different From Seat Frame Widths – Must choose back frame [Sec.7]

(configurations below are not available for recliner backs)

☐ ZNB_ _ _ _*	Custom Back – 2" Narrower than Seat Frame	$395	Seat Frame Width	16"	18"	20"	22"	24"
			Back Frame Width	14"[1]	16"[2]	18"[3]	20"[4]	22"[4]
* Enter base and back frame widths to complete code. eg. ZNB2018 specifies a 20"W chair frame with an 18"W back frame								
☐ ZNB_ _ _ _*	Custom Back – 4" Narrower than Seat Frame – requires custom seat	$395	Seat Frame Width	N/A	18"	20"	22"	24"
			Back Frame Width	N/A	14"[1]	16"[3]	18"[4]	20"[4]
* Enter seat and back frame widths to complete code. eg. ZNB2016 specifies a 20"W seat frame with a 18"W back frame								
☐ ZWB_ _ _ _*	Custom Back – 5" Wider than Seat Frame	$460	Seat Frame Width	16"	18"	20"	22"	N/A
			Back Frame Width	18"	20"	22"[3]	24"[4]	N/A
* Enter seat and back frame widths to complete code. eg. ZWB1820 specifies an 18"W seat frame with a 20"W back frame			Consider also using Arm Spacer Kit (ZASK05/10) as another solution for providing more width between the arms.					

[1] Requires Custom Seat

[2] Compatible with all 18" W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

[3] Custom Seat recommended. Also compatible with all 18" W seats, requires Seat Support Kit ZSSKT [Sec.3]

[4] Custom Seat recommended.

8. Back Support Options – These back supports replace the molded vinyl upholstery on a Fixed Back or tension-adjustable upholstery on an Adjusta-Back. Symphony and Harmony Backs are recommended for infection-control environments.

☐ ZSB14	14" W Symphony Back[1] – recessed solid back support	$340
☐ ZSB16	16" W Symphony Back[1] – recessed solid back support	$340
☐ ZSB18	18" W Symphony Back[1] – recessed solid back support	$340
☐ ZSB20	20" W Symphony Back[1] – recessed solid back support	$460
☐ ZSB22	22" W Symphony Back[1] – recessed solid back support	$525
☐ ZSB24	24" W Symphony Back[1] – recessed solid back support	$555
☐ ZHB	Harmony Back[2] (contoured solid back – 18"W x 13"H) Fits 17" or 20" back heights, 18"W back frames only	$340

[1] Does not fit 15" H backs

[2] Not recommended for reclining backs.

9. Arm Supports and Accessories

☐ ZARMTL	Long Arm Support Tubes (pr)	Standard
☐ ZARMTS	Short Arm Support Tubes (pr)[1]	N/C
☐ ZARMLL	Long Arm Support Pad 13" L – Left	Standard
☐ ZARMLR	Long Arm Support Pad 13" L – Right	Standard
☐ ZARMSL	Short Arm Support Pad 10" L – Left	N/C
☐ ZARMSR	Short Arm Support Pad 10" L – Right	N/C
☐ ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	$90
☐ Z2PLT	2-Point Arm Supports[2] (pr) with 13" L pads (std with reclining backs)	$525
☐ Z2PST	2-Point Arm Supports[2] (pr) with 10" L pads (std with reclining backs)	$525
☐ ZFAPLL	Long Molded Flat Armpad[2] – Left (14"L x 3.25"W)	$95
☐ ZFAPLR	Long Molded Flat Armpad[2] – Right (14"L x 3.25"W)	$95
☐ ZFAPSL	Short Molded Flat Armpad[3] – Left (10" L x 2.75" W – fits ZARM10)	$95
☐ ZFAPSR	Short Molded Flat Armpad[3] – Right (10" L x 2.75" W – fits ZARM10)	$95
☐ ZTRSL	Arm Trough - Small Left (16"L x4.25"W) *NEW*	$130
☐ ZTRSR	Arm Trough - Small Right (16"L x 4.25"W) *NEW*	$130
☐ ZTRML	Arm Trough - Medium Left (18.5"L × 5"W)	$130
☐ ZTRMR	AArm Trough - Medium Right (18.5"L x 5"W)	$130
☐ ZTRLL	Arm Trough - Large Left (21"L x5.25"W) *NEW*	$130
☐ ZTRLR	Arm Trough - Large Right (21"L x 5.25"W) *NEW*	$130
☐ ZPAML	Pivoting Arm Mount[5] – Left (for ZFAPLL / ZFAPLL / ZTRSL / ZTRML / ZTRLL)	$190
☐ ZPAMR	Pivoting Arm Mount[5] – Right (for ZFAPLR / ZFAPLR / ZTRSR / ZTRMR / ZTRLR)	$190
☐ ZASK05	Arm Support Spacer Kit[6] – 0.5" per side for all Adjusta-Backs (pr)	$70
☐ ZASK10	Arm Support Spacer Kit[6] – 1" per side for all Adjusta-Backs (pr)	$135
☐ ZPHGL	Pivoting Hand Grips (pr) for 13" L arms[7]	$460
☐ ZPHGS	Pivoting Hand Grips (pr) for 10" L arms[7]	$460
☐ ZAPSB_ _*	Anterior Posture Support Bar[8,9] * Enter frame width to complete code (14" / 16" / 18" / 20")	$225
☐ ZAPSB_ _*	Anterior Posture Support Bar[8,9] * Enter frame width to complete code (22" / 24")	$260

Reclining Back Arm Support Accessories

☐ Z2PLTR	2-Point Flip-up Arm Supports (Long) (pr)	N/C
☐ Z2PSTR	2-Point Flip-up Arm Supports (Short) (pr)	N/C
☐ ZASKR05	Arm Support Spacers 0.5" per side (pr) for Reclining Back Support	$70
☐ ZASKR10	Arm Support Spacers 1" per side (pr) for Reclining Back Support	$70

[1] Not compatible with ZPAM / ZFAPL / ZFAPS / ZTRM / ZTRL
[2] For chairs with same back width and frame width only
[3] May require ZPAM for Arm Support to flip up with Flat Arm Pad
[4] ZPAM is required for Arm Support to flip up with Arm Trough
[5] Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered.
[6] ZASK05 or ZASK10 is required for Standard Arm Support to flip up with lateral pad in forward position
[7] ZPHGL ZPHGS combined with 2 Point Arms requires a custom up charge
[8] Not compatible with ZFAP / ZTR
[9] Special length required with ZASK05 and ZASK10 (contact customer service)

10. Foot and Leg Supports

☐ ZFSL	Foot Support Left – 16" - 20.5"	Standard
☐ ZFSR	Foot Support Right – 16" - 20.5"	Standard
☐ ZFET2L	Footrest Extension Tube[1] Left – 2" longer	$30
☐ ZFET2R	Footrest Extension Tube[1] Right – 2" longer	$30
☐ ZMFX8L	MFX8 Foot Support Left – 8" - 12"	$45
☐ ZMFX8R	MFX8 Foot Support Right – 8" - 12"	$45
☐ ZMFX12L	MFX12 Foot Support Left – 12" - 16.5"	$55
☐ ZMFX12R	MFX12 Foot Support Right – 12" - 16.5"	$55
☐ ZADABS18L	Angle/Depth Adj ABS Footplate Left For 18" W frames	$110
☐ ZADABS18R	Angle/Depth Adj ABS Footplats Right For 18" W frames	$110
☐ ZADML	Adj Angle/Depth Molded Footplate[2] – Left	$110
☐ ZADMR	Adj Angle/Depth Molded Footplate[2] – Right	$110
☐ ZFET2AL	Footrest Extension Tube Left for ZADML – 2" longer	$30
☐ ZFET2AR	Footrest Extension Tube Right for ZADMR – 2" longer	$30
☐ ZVS_ _*	V-Style Foot Support (1-piece) * Enter frame width to complete code (16” / 18”)	$340
☐ ZVS_ _*	V-Style Foot Support (1-piece) * Enter frame width to complete code (20” / 22” / 24")	$375
☐ ZFBT	Flip-Back Foot Support (1-piece) with calf pads	$525
☐ ZELRL	Elevating Leg Support – Left (includes ZADML)	$570
☐ ZELRR	Elevating Leg Support – Right (includes ZADMR)	$570
☐ ZRLSL	Universal Residual Limb Support – Left	$225
☐ ZRLSR	Universal Residual Limb Support – Right	$225
☐ ZFOR	Forward Offset Foot Support Receiver (ea) Moves Foot support forward 2.5" to accommodate deeper seats (19" or deeper). Specify 1 or 2 __ ** Need to purchase ZFOR when seats are deeper then 19” d**	$140
Foot Support Options for 16"W Frames		
☐ ZFS16L	Foot Support for 16" W frames Left – (16" - 20.5")	N/C
☐ ZFS16R	Foot Support for 16" W frames Right – (16" - 20.5")	N/C
☐ ZADABS16L	Angle/Depth Adj ABS Footplate Left for 16" W frames	$110
☐ ZADABS16R	Angle/Depth Adj ABS Footplate Right for 16" W frames	$110
☐ ZMFX816L	MFX8 Foot Support for 16" W Frame Left – 8" - 12"	$45
☐ ZMFX816R	MFX8 Foot Support for 16" W Frame Right – 8" - 12"	$45
☐ ZMFX1216L	MFX12 Foot Support for 16" W Frame Left – 12" - 16.5"	$55
☐ ZMFX1216R	MFX12 Foot Support for 16" W Frame Right –12" - 16.5"	$55

11. Foot and Leg Support Accessories

☐ ZFSLK	Foot Support Lock (pr)	$110
☐ ZCS_ _*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (16" / 18" / 20" / 22"/ 24")	$155
☐ ZCSIC_ _*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (16" / 18" / 20")	$145
☐ ZCSIC_ _*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (22" / 24")	$165
☐ ZCPL	Calf Panel – Bodypoint Aeromesh[1] Fits frame widths (16” / 18” / 20”)	$210
☐ ZHS1820	H-Strap – for 18” / 20"W frames – hook and loop attachment	$105
☐ ZHS2224	H-Strap – for 22"/ 24"W frames – hook and loop attachment	$130
☐ ZAHL	Adjustable Heel Loops (pr) – not available for ZVS / ZFBT	$70
☐ ZNC	Neoprene Cover for Standard Footplate (Ea)	$55
☐ ZNCA	Neoprene Cover for Angle/Depth Adjustable (Each)	$55
☐ ZNC16	Neoprene Covers for Narrow Molded Footplates – 16" W frame (Each)	$55
☐ ZNCABS16	Neoprene Cover for Narrow Molded ABS Angle/Depth Adjustable (Each)	$55
☐ ZNCVS	Neoprene Cover for V-Style Foot Support	$55
☐ ZNCFB	Neoprene Cover for Flip-Back Foot Support	$55
☐ ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 __	$135
☐ ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 __	$135
☐ ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 6.2" to accomodate deeper seats and reduced knee flexion range. Specify 1 or 2 __	$140

[1] Not compatible with ZMFX12 / ZFBT / ZVS

12. Positioning Pads

Head Supports		
☐ ZHSP1	Head Support Pad – Small (10" x 4")	Standard
☐ ZHSP2	Head Support Pad – Large (12" x 12")	$60
☐ ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	$90
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 2.75" per side. For up to 4.3" medial offset, specify ZLAEP extension plate. (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
☐ ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left	$380
☐ ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right	$380
☐ ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left	$380
☐ ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right	$380
Leg Adductors (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
☐ ZLGT1L	Leg Adductor w/ Small Pad – Left	$380
☐ ZLGT1R	Leg Adductor w/ Small Pad – Right	$380
☐ ZLGT2L	Leg Adductor w/ Large Pad – Left	$380
☐ ZLGT2R	Leg Adductor w/ Large Pad – Right	$380
Hip Adductors (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
☐ ZHP1L	Hip Adductor w/ Small Pad – Left	$380
☐ ZHP1R	Hip Adductor w/ Small Pad – Right	$380
☐ ZHP2L	Hip Adductor w/ Large Pad – Left	$380
☐ ZHP2R	Hip Adductor w/ Large Pad – Right	$380
Hip Adductors For Reclining Backs (Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3")		
☐ ZHPR1L	Hip Adductor for Reclining Back w/ Small Pad – Left	$380
☐ ZHPR1R	Hip Adductor for Reclining Back w/ Small Pad – Right	$380
☐ ZHPR2L	Hip Adductor for Reclining Back w/ Large Pad – Left	$380
☐ ZHPR2R	Hip Adductor for Reclining Back w/ Large Pad – Right	$380
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors These plates allow for an additional 4cm medial offset and can be used on Lateral Thoracic Support, Hip and Leg Adductors		
☐ ZLAEPL	Lateral/Adductor Extension Plate – Left	$30
☐ ZLAEPR	Lateral/Adductor Extension Plate – Right	$30
Lateral Thoracic Support	☐ Left ☐ Right	
Leg Adductor	☐ Left ☐ Right	
Hip Adductor	☐ Left ☐ Right	
Total Quantity ___		
Medial Thigh Supports (Pommel) Pad sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3"		
☐ ZMTST1	Medial Thigh Support w/ Small Pads	$515
☐ ZMTST2	Medial Thigh Support w/ Large Pads	$515

13. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			
☐	ZPBS	Pelvic Belt – Small (frame widths 14-16" , 2-piece w/ side-release buckle)	$50
☐	ZPBM	Pelvic Belt – Medium (frame widths 16" - 18" , 2-piece w/ side-release buckle)	$60
☐	ZPBL	Pelvic Belt – Large (frame widths 18" - 22" , 2-piece w/ side-release buckle)	$55
☐	ZPBXL	Pelvic Belt – XLarge (seat frame widths 22" - 26", 2-piece w/ side-release buckle)	$75
☐	ZPBICS	Infection Control Pelvic Belt[1] – Small (frame widths 14" to 16" , 2-piece w/ side-release buckle)	$110
☐	ZPBICM	Infection Control Pelvic Belt[1] – Medium (frame widths 16" - 18" , 2-piece w/ side-release buckle)	$110
☐	ZPBICL	Infection Control Pelvic Belt[1] – Large (frame widths 18" -22" , 2-piece w/ side-release buckle)	$110
☐	ZPBICXL	Infection Control Pelvic Belt[1] – XLarge (seat frame widths 22" - 26", 2-piece w/ side-release buckle)	$125
☐	ZPBBPM	Pelvic Belt – Bodypoint Aeromesh – Med	$295
☐	ZPBBPL	Pelvic Belt – Bodypoint Aeromesh – Large	$295
☐	ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh – XLarge	$175
☐			
☐			

Chest Belts			
☐	ZCBS	Chest Belt[2] – Small (back widths 14-16" , 1-piece w/ hook and loop closure)	$40
☐	ZCBM	Chest Belt[2] – Medium (back widths 16" - 18" , 1-piece w/ hook and loop closure)	$45
☐	ZCBL	Chest Belt[2] – Large (back widths 18" -22" , 1-piece w/ hook and loop closure)	$55
☐	ZCBXL	Chest Belt – XLarge[2] (back frame widths 22" - 26", 1-piece w/ hook and loop closure)	$70
☐	ZCBICS	Infection Control Chest Belt[1,2] – Small (back widths 14" to 16" , 1-piece w/ side-release buckle)	$90
☐	ZCBICM	Infection Control Chest Belt[1,2] – Medium (back widths 16" -18" , 1-piece w/ side-release buckle)	$90
☐	ZCBICL	Infection Control Chest Belt[1,2] – Large (back widths 18" -22" , 1-piece w/ side-release buckle)	$100
☐	ZCBICXL	Infection Control Chest Belt[1,2] – XLarge (back frame widths 22" - 26", 1-piece w/ side-release buckle)	$125
☐	ZCBBPM	Chest Belt – Bodypoint Aeromesh – Med	$310
☐	ZCBBPL	Chest Belt – Bodypoint Aeromesh – Large	$310
☐	ZCBBPXL	Chest Belt – Bodypoint Aeromesh – XLarge	$310

Recommended sizes for Bodypoint belts

Pelvic belt	ZPBBPS	14-16" W seat frame	ZPBBPM	16" /18" W seat frame	ZPBBPL	18" -22" W seat frame	ZPBBPXL	22" - 26" W seat frame
Chest belt	ZCBBPS	14-16" W back frame	ZCBBPM	16" /18" W back frame	ZCBBPL	18" -20" W back frame	ZCBBPXL	22" - 26"W back frame

[1] Non-absorbent belts – recommended for infection-control environments

[2] Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM / ZCBL / ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

14. Urine Management

☐	ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	$75
☐	ZSS	SplaShield – Molded urine deflector (must be held in place during use)	$170
☐	ZSC17	Splash Curtain[1] for 17" D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	$345
☐	ZSC20	Splash Curtain[1] for 20" D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	$390
☐	ZSCFR	Splash Curtain for 17" / 20" D (at front of aperture only, not for Front Open Seats)	$405
☐	ZSCC	Custom Splash Curtain for Custom Seat – Raises seat 1/2"	$420
☐	ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	$390
☐	ZBUDBH	Bard urine drainage bag holder (available for all frame widths)	$190

[1] Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for order codes and pricing.

15. Other

☐ ZDCLT	Directional Caster Lock – fits on one rear caster (facilitates steering)	$195
☐ ZATTRR	Anti-Tippers[1] – Rear (pr)	$395
☐ ZATTFR	Anti-Tippers[1] – Front (pr)	$395
☐ ZCUT4	4" Caster Upgrade (4x)	$140
☐ ZCUT6	6" Caster Upgrade (4x)	$400
☐ ZTLS	Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx 5°. Can be stacked for additional tilt limiting. Quantity ___	N/C

[1] Anti tippers for 4" / 5" / 6" casters available.

☐ ZPSK1	Commode Pan Spacer Kit[2] – 1" lower	$225
☐ ZPSK2	Commode Pan Spacer Kit[2] – 2" lower	$245
☐ ZPSK3	Commode Pan Spacer Kit[2] – 3" lower	$285
☐ ZPSK4	Commode Pan Spacer Kit[2] – 4" lower	$285
☐ ZSPKT22	AT Self-Propel Kit with 22" wheels – incl anti-tippers (tilt limited to 25°)	$1,765
☐ ZSPKT24	AT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°)	$1,765

[2] Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.

16. Credits and Omissions

☐ Z905	Omit Commode Pan and Mounting Brackets	–$20	☐ Z911	Omit Molded Head Support Pad	–$20
☐ Z907	Omit Foot Supports – if no other Foot / Leg Support Option is selected	–$45	☐ Z912	Omit Complete Head Support Assembly	–$50

Some assembly of Raz products is required.

17. Notes:

Specifications	Raz-AT (Attendant Tilt)
Weight Capacity	407 lbs
Overal Width	Seat width plus 4.5" with Standard Arm Pads
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports - at 40°	43" 54"
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat) - with pan brackets - without pan brackets	 17"-21" (1" increments) 17.75"-21.75 (1" increments)
Toilet Clearance Width	15.75" for standard 18"W seat frame. For other widths, seat frame width minus 2.25"
Arm Support Height - above molded seat (1/2" increments)	
- **18"H Fixed Back - standard** - 15"H Adjusta-Back - 17"H Adjusta-Back - 20"H Adjusta-Back - 24"H Adjusta-Back - 20"H Reclining Adjusta-Back - 24"H Reclining Adjusta-Back	**8.5"** 5.5"-8.5" (1/2" increments) 6.5"-9.5" (1/2" increments) 7.5"-10.5" (1/2" increments) 7.5"-10.5" (1/2" increments) 7.5"-10.5" (1/2" increments) 7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	21"-25" in 1" increments
Product Weight (with 18"W seat)	61 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)