160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T.+1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



6.8

VELOCE
ORDER FORM - INTERNATIONAL

1- Purchasing information

P.O. number:* Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

2- ATP/Clinician information

ATP/Clinician:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐ ** Required*

3- Shipping information

Ship to:*

Name:

Address:

City: State:

Zip Code:

Phone:

MA083 ☐ *English Manual (STD)* *MA084* ☐ *German Manual* *MA085* ☐ *Italian Manual* *MA086* ☐ *Dutch Manual* *MA087* ☐ *Spanish Manual* *MA088* ☐ *French Manual* *MA089* ☐ *Portuguese Manual*

MODEL								USD (Taxes not included)
WC20	☐	Ultralight **VELOCE** Wheelchair *(120 kg - 265 lb weight limit)*						**$5,295**
WC09	☐	**HD** Heavy Duty package *(159 kg - 350 lb weight capacity)*						$540
Frame								
SEAT WIDTH [3] - *cm (inches)*								
SW	☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") [2] ☐	48.26 (19") [2] ☐	NC
	☐			50.80 (20") [2] ☐	53.34 (21") [7] ☐	55.88 (22") [7] ☐		$535
SEAT DEPTH - *cm (inches)*								
SD	☐	35.56 (14") [1] ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ☐	48.26 (19") ☐	NC
	☐			50.80 (20") ☐				$515
FRAME LENGTH *(Seat sling depth plus distance from seat sling to front of frame)* [4]								
SP01	☐	Standard frame length configuration [5]						STD
SP02	☐	Frame length 2.54 cm (1") longer than standard configuration						NC
SP03	☐	Frame length 5.08 cm (2") longer than standard configuration						NC
SP04	☐	Frame length 2.54 cm (1") shorter than standard configuration						NC
SP05	☐	Frame length 5.08 cm (2") shorter than standard configuration						NC
TRANSIT TIE DOWN								
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt [6]	Short belt 50-104 cm (20" - 41")					$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt [6]	Long belt 50-116 cm (20" - 46")					$365
AC35	☐	WC-19 Transit Tie Down without belt [6]						$335

1- *Only available with frame length 1" and 2" longer;* ***2-*** *Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;*

3. Width calculator

4- *See chart on page 7;* ***5-*** *4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option;*


7- *Only available with HD option and the one piece adjustable angle flip-up*

FRONT Seat-to-floor height

Caster sizes		Height available +- 1/4"	
☐	3" Caster	35.56cm to 53.34cm *(14" to 21")*	**Specify desired height:**
☐	4" Caster	36.83cm to 53.34cm *(14.5" to 21")*	
☐	5" Caster	39.37cm to 53.34cm *(15.5" to 21")*	(inches)
☐	6" Caster	41.91cm to 53.34cm *(16.5" to 21")* (N/A with 90° front frame)	

Frog Legs FRONT Seat-to-floor height

Caster Sizes		Height available +- 1/4"	
☐	3" Caster	38.7cm to 53.3cm *(15.25'' to 21'')*	**Specify desired height:**
☐	4" Caster	40.0cm to 53.3cm *(15.75'' to 21'')*	
☐	5" Caster	42.54cm to 53.3cm *(16.75'' to 21'')*	(inches)
☐	6" Caster	45.08cm to 53.3cm *(17.75'' to 21'')* (N/A with 90° front frame)	



FRONT CASTERS

CASTERS				Availability		
CA01-03 to CA01-06	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		NC
CA10	☐	Pneumatic casters (1 ¼")		6"		$175
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		$105
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		$210
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks				$515

REAR Seat-to-floor Height

	WHEEL SIZES	Height available	
☐	22'' (501) wheel	34.29cm to 46.35cm *(13.5" to 18.25")*	**Specify desired height:**
☐	24'' (540) wheel	36.83cm to 48.89cm *(14.5" to 19.25")*	
☐	25'' (559) wheel	38.10cm to 50.16cm *(15" to 19.75")*	_____ (inches)
☐	26'' (590) wheel	39.37cm to 51.43cm *(15.5" to 20.25")*	

Configuration used: 2.5 cm (1") seat slope with a center of gravity at 5cm (2")

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (½") of selected Seat to Floor height"

REAR WHEELS

ONE ARM DRIVE [1] *Only available with Newton One rear wheels, permanent axles, 0° camber*							
HR22	☐	One arm drive - aluminum inner handrim [2]	Left ☐	Right ☐			$1,510
HR60	☐	One arm drive - plastic coated inner handrim [2]	Left ☐	Right ☐			$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)					NC
REAR WHEELS					**Available**		
RW01	☐	Mag [3]			22" 24"		$57
RW02	☐	Newton One - spoke wheel [4]			22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel			22", 24", 25", 26"		$350
RW04	☐	Spinergy Spox	**Blk STD** ☐	Wht [5] ☐ Blue [5] ☐ Red [5] ☐ Yel [5] ☐	22", 24", 25", 26"		$895
RW05	☐	Spinergy LX	**Blk STD** ☐	Pur [5] ☐ Wht [5] ☐ Blue [5] ☐ Red [5] ☐ Yel [5] ☐ Pink [5] ☐ Grn [5] ☐ Org [5] ☐	22", 24", 25", 26"		$1,125
RW07	☐	Spinergy CLX [6]	**Blk STD** ☐	Pur [5] ☐ Wht [5] ☐ Blue [5] ☐ Red [5] ☐ Yel [5] ☐ Pink [5] ☐ Grn [5] ☐ Org [5] ☐	24", 25"		$1,600
TIRES (*The brand and color may vary depending on parts availability)*					**Available**		
		AIRLESS					
TI01	☐	Soft Urethane 3.49cm (1 3/8") - Low Tread, black			22" 24"		NC
TI02	☐	Soft Urethane 3.49cm (1 3/8") - Medium Tread, black [7]			22", 24", 25", 26"		$105
TI13	☐	Soft Urethane 2.54cm (1") - No tread, black - *Sentinel*			20",22", 24", 25", 26"		$125
TI08	☐	Pneumatic Airless insert 3.49cm (1 3/8") - Medium Tread [8]			22", 24", 26"		$90
		PNEUMATIC					
TI05	☐	Pneumatic 3.49cm (1 3/8") Medium Tread			22", 24", 26"		NC
TI06	☐	Pneumatic HP 2.54cm (1") - Low Tread, Puncture resistant, black - *SpeedLite*			22", 24", 25", 26"		$110
TI07	☐	Pneumatic HP 2.54cm (1") - Medium Tread, Puncture resistant, black - *TrailBlazer*			24", 25", 26"		$145
TI03	☐	Pneumatic 5.08cm (2") - Large Tread, black			22", 24", 25"		$115

***1-** Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; **2-** Only standard handrim position available;*
***3-** Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; **4-** Not compatible with Shox tires; **5-** color spokes require an additionnal 10 days lead-time;*
***6-** Only available with 1" tires; **7-** 22" is only available with Mags, Newton One rear wheels and Newton Gravity wheels*
***8** - Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels*

HANDRIMS

HANDRIM				Available		
HR01	☐	Aluminum anodized		22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized [2]		22", 24", 25"		$80
HR02	☐	Plastic coated [1]		22", 24", 25", 26"		$145
HR03	☐	High friction coated [1,2,3]		22", 24", 25", 26"		$380
HR07-HR18 and HR65-HR71	☐	Newton Air Grip [2]	*See page 6 for the color selection and price*	22", 24", 25", 26"		
HR28	☐	Surge [1,2,3,4]	*(1 1/2" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"		$425
HR05	☐	Surge LT [1,2,3,4]	*(1 1/8" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"		$465
HR04	☐	Natural Fit [1,2,3,4]		22", 24", 25", 26"		$395
HR62	☐	Simi H [1,2,3]	*(1 3/8" Oval handrim with high friction strip)*	22", 24", 25", 26"		$395
HR63	☐	Nova H [1,2,3]	*(1 1/4" Oval handrim with high friction strip)*	22", 24", 25", 26"		$425
HR64	☐	Optimum H [1,2,3]	*(1 1/8" heart-shaped handrim)*	22", 24", 25", 26"		$325
HR99	☐	Omit handrim				NC

***1-** Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; **2-** Not available with Mag wheels; **3** - Not available with one arm drive;*
***4** - 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece*

HANDRIMS OPTION

ASSEMBLY POSITION							
HR29	☐	Standard Handrim position	*Full tab*				STD
HR030	☐	Narrow Handrim position no cut [6]	*Full tab*				NC
HR30	☐	Narrow Handrim position with cut [5]	*Short tab 2 positions*				$32
HR31	☐	Super Narrow Handrim [5]	*Short tab 2 positions STD* ☐	*Short tab 1 position* [6] ☐			$32
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS							
HR23	☐	Standard grip					$95
HR27	☐	Super grip					$150
HR099	☐	Omit thumb piece					NC


***5-** Only available with Aluminum and Air Grip handrims; 6- Not available with the combination of 24" Newton Gravity, Spinergy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;*

AXLE [1]					
AX01	☐	Quick release axle			STD
AX03	☐	Quad release axle			$135

1- Not required with one arm drive

REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM03	☐	3° camber			NC
CM06	☐	6° camber			$165

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =

Most stable												
2,5	3,2	3,8	4,4	**5,1**	5,7	6,3	7	**7,6**	8,3	8,9	9,5	**10,2**
1"	1/4"	1/2"	3/4"	**2"**	1/4"	1/2"	3/4"	**3"**	1/4"	1/2"	3/4"	**4"**
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 3,8 cm/1.5''

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6,3 cm (2.5").

Possible conflict between a 26" rear wheel and a center of gravity over 6,3 cm (2.5").

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG


WHEEL LOCK					
WL08	☐	Newton Wheel lock - push to lock			STD
WL02	☐	Pull to lock			NC
WL04	☐	Newton grade aid push to lock *(works best with pneumatic or high pressure tire)*			$165
WL13 WL14	☐	6" removable extension handle [1]			$83
WL05	☐	Scissors Lock			$155
WL15	☐	Unilateral wheel lock - pull to lock [2] Left ☐ Right ☐			$335
WL99	☐	Omit wheel locks			NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilties and restriction, refer to the table at page 6 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents


FRONT FRAME ANGLE					
FA70	☐	70° Front frame angle			NC
FA75	☐	75° Front frame angle			NC
FA80	☐	80° Front frame angle			NC
FA85	☐	85° Front frame angle			NC
FA90	☐	90° Front frame angle			NC




SEAT TO FOOTREST LENGTH

Set Length to[1]:	*Length must be at least 4,5 cm (1 ¾") shorter than front seat-to-floor height*

1 - +/- 1.27 cm (1/2")



Footrest width / Taper selection		
Seat Width	**Straight [1]**	**1" Taper**
Outside Measurement	**6,4 cm Narrower than Seat Width**	**11,4 cm Narrower than Seat Width**
35.56 *14"*	29.21cm/11.5'' ☐	24.13cm/9.5'' ☐
38.10 *15"*	31.75cm/12.5'' ☐	26.67cm/10.5'' ☐
40.64 *16"*	34.29cm/13.5'' ☐	29.21cm/11.5'' ☐
43.18 *17"*	36.83cm/14.5'' ☐	31.75cm/12.5'' ☐
45.72 *18"*	39.37cm/15.5'' ☐	34.29cm/13.5'' ☐
48.26 *19"*	41.91cm/16.5'' ☐	36.83cm/14.5'' ☐
50.80 *20"*	44.45cm/17.5'' ☐	39.37cm/15.5'' ☐
53.34 *21"*	46.99cm/18.5'' ☐	41.91cm/16.5'' ☐
55.88 *22"*	49.53cm/19.5'' ☐	44.45cm/17.5'' ☐



1- Not compatible with auto-folding footplate

FOOTPLATES							
FP06	☐	One piece adjustable angle flip-up - Aluminum					NC
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts					$470
FP28	☐	One piece auto-folding adjustable angle - Aluminum					$310
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts					$380
FOOTPLATE ACCESSORIES							
FP09	☐	Calf strap					$43
FP30	☐	Bodypoint Aeromesh Padded calf strap					$125

SEATING AND BACKREST

BACK UPHOLSTERY							
BH07	☐	Slip-on Active-Mesh back upholstery					NC
BH08	☐	Adjustable tension Active-Mesh back upholstery					$190
BH99	☐	Omit Back upholstery					NC
	☐	NXT back Support				See price list	
SEAT CUSHION							
SC01	☐	5.08 cm (2") Cushion					$194
SC02	☐	7.62 cm (3") Cushion					$200
SC99	☐	Omit seat cushion					STD
	☐	NXT cushion				See price list	
HEIGHT ADJUSTABLE WITH PUSH HANDLES *(MUST select ONE push handle options)*							
BC06-912	☐	Straight Carbon fiber back cane	22.86 - 30.48cm (9-12")	*increment of 1/2"*	Set at [2]:		NC
BC06-1215	☐	Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	*increment of 1/2"*	Set at [2]:		NC
BC06-1518	☐	Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	*increment of 1/2"*	Set at [2]:		NC
BC06-1821	☐	Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	*increment of 1/2"*	Set at [2]:		NC
BC08-1315	☐	8° Bend Carbon fiber back cane	34.29 - 38.10cm (13 ½-15")	*increment of 1/2"*	Set at [2]:		$60
BC08-1518	☐	8° Bend Carbon fiber back cane	38.10 - 45.72 (15-18")	*increment of 1/2"*	Set at [2]:		$60
BC08-1821	☐	8° Bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	*increment of 1/2"*	Set at [2]:		$60
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*							
BC09-1215		Straight Carbon fiber back cane	30.48 - 38.10cm (12-15")	*increment of 1/2"*	**Not available until fall 2026**		$220
BC09-1518		Straight Carbon fiber back cane	38.10 - 45.72cm (15-18")	*increment of 1/2"*	**Not available until fall 2026**		$220
BC09-1821		Straight Carbon fiber back cane	45.72 - 53.37cm (18-21")	*increment of 1/2"*	**Not available until fall 2026**		$220
BC10-1315		8° bend Carbon fiber back cane	34.29 - 38.10cm (13 ½-15")	*increment of 1/2"*	**Not available until fall 2026**		$275
BC10-1518		8° bend Carbon fiber back cane	38.10 - 45.72 (15-18")	*increment of 1/2"*	**Not available until fall 2026**		$275
BC10-1821		8° bend Carbon fiber back cane	45.72 - 53.34cm (18-21")	*increment of 1/2"*	**Not available until fall 2026**		$275
		Select back cane angle:	85° [1] □	90° □	95° **STD** □	100° □	105° □ 110° □
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*							
BC13-1518	☐	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	*increment of 1/2"*	Set at[2]:		$275
BC13-1821	☐	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	*increment of 1/2"*	Set at[2]:		$275
BC14-1518	☐	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	*increment of 1/2"*	Set at[2]:		$275
BC14-1821	☐	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	*increment of 1/2"*	Set at[2]:		$275
		Select back cane angle:	85° [1] ☐	90° ☐	95° **STD** ☐	100° ☐	105° ☐ 110° ☐
PUSH HANDLES OPTIONS *(ONE option at a time)*							
PH01	☐	Push handles					NC
PH05	☐	Fold down push handles [3, 5]					$150
PH06-920	☐	Clamp on height adjustable push handles 22.86 -50.80cm (9" - 20") [3,4, 5]					$330
PH06-923	☐	Clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") [3,4, 5]					$330
PH02	☐	Omit push handles					NC

1- *Not available with a 8° bend back cane:* ***2-*** *If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery:* ***3-*** *120 kg/265lb weight limit*
4- *Not compatible with swing away armrests, 8° bend back canes and 9-12" back canes;* ***5-*** *Not compatible with headrest and attendant lock*

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST ***(convertible to "T" armrest)*** [1]							
AR01-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐		$390
AR01-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐		$390
AR02-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐		$390
AR02-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐		$390
HEIGHT ADJUSTABLE "T" ARMREST [1]							
AR13-610	☐	25.4cm (10") Desk length armpad [4]	15.25 - 25cm (6" - 10")				$360
AR03-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")				$360
AR03-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")				$360
AR14-610	☐	35.56cm (14") Full length armpad [4]	15.25 - 25cm (6" - 10")				$360
AR04-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")				$360
AR04-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")				$360
HEIGHT ADJUSTABLE SWING AWAY ARMREST [2]							
AR05	☐	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")				$250
AR23	☐	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")				$250
AR06	☐	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")				$250
AR17	☐	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")				$250
ARMREST OPTIONS							
AR12	☐	35.56cm (14") Full Length Armpad - *Gel Ovation (pair)* [3]					$165
AR99	☐	Omit armrest					NC

1- *May have a 1/2 " adjustment restriction on both limits of the range;* ***2-*** *Could interfere with some rigid back rests; Not available with fixed carbon fiber side guards*
3- *Only available with full length armrest; Not available with swing away and desk lenght armrests;* ***4-*** *May have some restriction with seat-to-floor height and rear wheel*

ACCESSORIES

SIDE GUARDS *(not available with Flip Back and "T" Armrest)*							
SG01	☐	Plastic side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$210
SG02	☐	Carbon fiber side guards *(fixed)*					$275
SG03	☐	Carbon fiber side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$315
SG06	☐	Plastic fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$295
SG04	☐	Carbon fiber fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐			$550
SG99	☐	Omit side guard					NC

1- *only compatible with 0° camber;Increase width by 1/2" per side for a total of 1"*

POSITIONING BELT							
SB01-S	☐	Velcro adjustable - **Small 88.9cm (35")** [1]					$50
SB01-L	☐	Velcro adjustable - **Long 152.40cm (60")**					$50
SB02-S	☐	Auto buckle - **Small 88.9cm (35")** [1]					$95
SB02-M	☐	Auto buckle - **Medium 121.92cm (48")**					$95
SB02-L	☐	Auto buckle - **Long 152.40cm (60")**					$95
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - **Medium 142.24cm (56")**					$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - **Long 154.94cm (61")**					$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - **Small129.54cm (51")**					$185
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - **Medium 142.24cm (56")**					$185
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - **Long 154.94cm (61")**					$185
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - **Small 129.54cm (51")**					$205
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - **Medium 142.24cm (56")**					$205
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - **Long 154.94cm (61")**					$205
SB03	☐	Hardware attachment for 2 extra points (1 pair)					$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - **Medium 73.66 (18")**					$310
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - **Long 93.98cm (24")**					$310
SB99	☐	Omit positioning belt					STD

1- Not recommended for seat width of 45.72cm (18") and more

HEADREST [1]							
AC01-S	☐	Small contoured head support pad (22.9x12.7cm/9"x5")					$195
AC01-M	☐	Large rectangular head support pad (25.4x12.7 cm/10"x5")					$195
AC01-L	☐	XL rectangular head support pad (35.6x12.7cm/14"x5")					$195
AC28	☐	Headrest mounting hardware with tension bar					$320

1- *Only available with push handles*

ACCESSORIES							
WC10	☐	E-Motion adaptor Kit (2 adaptors) - *Only necessary with 6° camber*					$35
AC03	☐	Anti-tippers					$155
AC29	☐	Swing-away anti-tipper[1]	Left ☐	Right ☐			$165
AC27	☐	Tip assist [2]					$67
AC06	☐	Cane and crutch holder	Left ☐	Right ☐			$130
AC08	☐	Spoke guards [3]	Available: 22", 24", 25", 26"				$125
AC10	☐	Seat pouch					$83
AC30	☐	Newton multi-purpose pouch					$62
AC26	☐	BodyPoint Mobility Bag					$135
AC23	☐	Newton Backpack					$130
AC12	☐	Velcro Closing Strap					$57
AC13	☐	Tool kit					$50

1- *Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); Not available with amputee axle plate* ***2-*** *Cannot be combined with Anti Tippers*
3- *Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel*

Colors

Options [1]	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
																					
AirGrip	☐ $125		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	$265
Frame (NC)		☐ (STD)	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	TBD[3] ☐
																					
Highlight [2]	☐		☐					☐		☐	☐			☐			☐				NC
FC99	☐	Custom Color	Specify PPG Paintit color:																	*additional 10 days lead-time*	$390

1- *Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;*
2- To see different highlight, visit our website: *https://www.motioncomposites.com/en_intl/products/folding-wheelchairs/veloce*

3- *The frame color will be confirmed with the quote conversion into an order (in absence of selection, the color will be black)*

Frame length chart

Frame length	Seat Depth						
	35.56 (14")	38.10 (15")	40.64 (16")	43.18 (17")	45.72 (18")	48.26 (19")	50.80 (20")
STD	N/A						
1							N/A
2						N/A	N/A
-1	N/A	N/A					
-2	N/A	N/A	N/A				

Unilateral wheel lock restrictions

Seat depth: shall be 40.64 cm and greater

Rear wheel: shall be 24 in and less

"T" or flip-back armrest only compatible with 0° camber

One arm drive: not compatible with "T"

One arm drive and flip-back armrest only available with 17 in depth frame

"T" or flip-back armrest may require to move the rear wheels outwards

Not compatible with wheel lock lever extension

Rear wheel diameter	Seat depth (cm)	COG Compatibility Range (cm)
20 in	40.64+	4.44 to 10.16
22 in	40.64+	6.35 to 10.16
24 in	40.64+	2.54 to 10.16

EVO CUSTOM

CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's subtle tweaks for comfort and mobility improvements or custom colors that make your accessories uniquely yours, we're here to support you.

Explore the **EVO Custom Catalog**.

You can include your request or reference an existing customization from our catalog in the notes section below.

Notes :

FOLDING ON CHAIR – INTL ORDER FORM

BioFit™ CUSHION



☐ Use wheelchair configuration

Or use this size:

NXCABF _____ (W" x D")

DEPTH \ WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
12"	$425	$425	$425		N/A	N/A	N/A	N/A	N/A	N/A	N/A
13"						N/A	N/A	N/A	N/A	N/A	N/A
14"	$425	$425	$425								
15"		$425	$425	$425	$425						
16"		$425	$425	$425	$425	$425	$425	$460	$460	$460	$460
17"	N/A			$425	$425	$425	$425	$460	$460		
18"	N/A			$425	$425	$425	$425	$460	$460	$460	$460
19"	N/A			$425	$425	$425	$425	$460	$460		
20"	N/A			$425	$425	$425	$460	$460	$460	$460	$460

Option: Change cover for incontinent/infection control outer cover ☐ N/C

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of $160.00 and has a 2 to 4 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES



MODEL #	Description	Fits Back Width	Retail	QTY
RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	$77	
RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	$77	
RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	$77	

THORACIC SHAPE



Optima™

THORACIC BACK SUPPORT

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOACB1413	$695	☐	NIB-2193V	$125	☐
14"	16"	NOACB1416	$695	☐	NIB-2194V	$125	☐
15"	13"	NOACB1513	$695	☐	NIB-2197V	$125	☐
15"	16"	NOACB1516	$695	☐	NIB-2198V	$125	☐
16"	13"	NOACB1613	$695	☐	NIB-2201V	$125	☐
16"	16"	NOACB1616	$695	☐	NIB-2202V	$125	☐
16"	19"	NOACB1619	$695	☐	NIB-2203V	$125	☐
17"	13"	NOACB1713	$695	☐	NIB-2206V	$125	☐
17"	16"	NOACB1716	$695	☐	NIB-2207V	$125	☐
18"	13"	NOACB1813	$695	☐	NIB-2208V	$125	☐
18"	16"	NOACB1816	$695	☐	NIB-2211V	$125	☐
18"	19"	NOACB1819	$695	☐	NIB-2212V	$125	☐
19"	16"	NOACB1916	$695	☐	NIB-2213V	$125	☐
19"	19"	NOACB1919	$695	☐	NIB-2215V	$125	☐
20"	13"	NOACB2013	$695	☐	NIB-2218V	$125	☐
20"	16"	NOACB2016	$695	☐	NIB-2219V	$125	☐
20"	19"	NOACB2019	$695	☐	NIB-2220V	$125	☐
21"	16"	NOACB2116	$695	☐	NIB-2221V	$125	☐
21"	19"	NOACB2119	$695	☐	NIB-2222V	$125	☐
22"	13"	NOACB2213	$745	☐	NIB-2223V	$125	☐
22"	16"	NOACB2216	$745	☐	NIB-2224V	$125	☐
22"	19"	NOACB2219	$745	☐	NIB-2225V	$125	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™

THORACIC BACK SUPPORT
CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover &	Retail	✓
14"	13"	NOCFACB1413	$900	☐	NIB-2193V	$125	☐
14"	16"	NOCFACB1416	$930	☐	NIB-2194V	$125	☐
15"	13"	NOCFACB1513	$900	☐	NIB-2197V	$125	☐
15"	16"	NOCFACB1516	$930	☐	NIB-2198V	$125	☐
16"	13"	NOCFACB1613	$900	☐	NIB-2201V	$125	☐
16"	16"	NOCFACB1616	$930	☐	NIB-2202V	$125	☐
17"	13"	NOCFACB1713	$900	☐	NIB-2206V	$125	☐
17"	16"	NOCFACB1716	$930	☐	NIB-2207V	$125	☐
18"	13"	NOCFACB1813	$900	☐	NIB-2211V	$125	☐
18"	16"	NOCFACB1816	$930	☐	NIB-2212V	$125	☐
20"	13"	NOCFACB2013	$900	☐	NIB-2218V	$125	☐
20"	16"	NOCFACB2016	$930	☐	NIB-2219V	$125	☐

THORACIC SHAPE



Optima™ DEEP

THORACIC BACK SUPPORT

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOPCB1413	$695		NIB-2246V	$125	
14"	16"	NOPCB1416	$695		NIB-2247V	$125	
15"	13"	NOPCB1513	$695		NIB-2250V	$125	
15"	16"	NOPCB1516	$695		NIB-2251V	$125	
16"	13"	NOPCB1613	$695		NIB-2254V	$125	
16"	16"	NOPCB1616	$695		NIB-2255V	$125	
16"	19"	NOPCB1619	$695		NIB-2256V	$125	
17"	13"	NOPCB1713	$695		NIB-2259V	$125	
17"	16"	NOPCB1716	$695		NIB-2260V	$125	
18"	13"	NOPCB1813	$695		NIB-2264V	$125	
18"	16"	NOPCB1816	$695		NIB-2265V	$125	
18"	19"	NOPCB1819	$695		NIB-2266V	$125	
19"	16"	NOPCB1916	$695		NIB-2267V	$125	
19"	19"	NOPCB1919	$695		NIB-2268V	$125	
20"	13"	NOPCB2013	$695		NIB-2271V	$125	
20"	16"	NOPCB2016	$695		NIB-2272V	$125	
20"	19"	NOPCB2019	$695		NIB-2273V	$125	
21"	16"	NOPCB2116	$695		NIB-2274V	$125	
21"	19"	NOPCB2119	$695		NIB-2275V	$125	
22"	13"	NOPCB2213	$745		NIB-2276V	$125	
22"	16"	NOPCB2216	$745		NIB-2277V	$125	
22"	19"	NOPCB2219	$745		NIB-2278V	$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™ DEEP

THORACIC BACK SUPPORT

CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover &	Retail	✓
14"	13"	NOCFPCB1413	$985		NIB-2246V	$125	
14"	16"	NOCFPCB1416	$1,025		NIB-2247V	$125	
15"	13"	NOCFPCB1513	$985		NIB-2250V	$125	
15"	16"	NOCFPCB1516	$1,025		NIB-2251V	$125	
16"	13"	NOCFPCB1613	$985		NIB-2254V	$125	
16"	16"	NOCFPCB1616	$1,025		NIB-2255V	$125	
17"	13"	NOCFPCB1713	$985		NIB-2259V	$125	
17"	16"	NOCFPCB1716	$1,025		NIB-2260V	$125	
18"	13"	NOCFPCB1813	$985		NIB-2264V	$125	
18"	16"	NOCFPCB1816	$1,025		NIB-2265V	$125	
20"	13"	NOCFPCB2013	$985		NIB-2271V	$125	
20"	16"	NOCFPCB2016	$1,025		NIB-2272V	$125	

THORACIC SHAPE



Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

xtend™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NHAACB(L)1213/17	12"	13-17"	$655		NIB-2148V	$125	
NHAACB(L)1413/17	14"	13-17"	$655		NIB-2149V	$125	
NHAACB(L)1613/17	16"	13-17"	$655		NIB-2150V	$125	
NHAACB(L)1813/17	18"	13-17"	$655		NIB-2151V	$125	
NHAACB (L)2013/17	20"	13-17"	$655		NIB-2152V	$125	
NHAACB (L)2213/17	22"	13-17"	$825		NIB-2153V	$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.



Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPC

xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover &	Retail	✓
NHAPCB1616/20	16"	16"-20"	$655		NIB-2143V	$125	
NHAPCB1816/20	18"	16"-20"	$655		NIB-2144V	$125	
NHAPCB2016/20	20"	16"-20"	$655		NIB-2145V	$125	
NHAPCB2216/20	22"	16"-20"	$825		NIB-2146V	$125	

PELVIC SHAPE



Classic™ PELVIC BACK SUPPORT

Classic™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NACB1616	16"	16"	$655		NIB-2163V	$125	
NACB1618	16"	18"	$655		NIB-2164V	$125	
NACB1816	18"	16"	$655		NIB-2166V	$125	
NACB1818	18"	18"	$655		NIB-2167V	$125	
NACB1820	18"	20"	$655		NIB-2168V	$125	
NACB2020	20"	20"	$655		NIB-2171V	$125	
NACB2216	22"	16"	$825		NIB-2172V	$125	
NACB2218	22"	18"	$825		NIB-2173V	$125	
NACB2220	22"	20"	$825		NIB-2174V	$125	



Classic™ DEEP PELVIC BACK SUPPORT

Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover &	Retail	✓
NPCB1616	16"	16"	$655		NIB-2175V	$125	
NPCB1618	16"	18"	$655		NIB-2176V	$125	
NPCB1818	18"	18"	$655		NIB-2179V	$125	
NPCB1820	18"	20"	$655		NIB-2180V	$125	
NPCB2020	18"	20"	$655		NIB-2183V	$125	
NPCB2216	20"	16"	$825		NIB-2184V	$125	
NPCB2218	22"	18"	$825		NIB-2185V	$125	
NPCB2220	22"	20"	$825		NIB-2186V	$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN USD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT ACCESSORIES



MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	$160	
NIB-1938	Extra Small Privacy Flap for 11” - 13” wide back supports	11” - 13” (27.9 - 33.0)	$65	
NIB-1078	Small Privacy Flap for 14”-16” wide back supports	14” - 16” (35.6-40.6 cm)	$65	
NIB-1079	Medium Flap for 17”-19” wide back supports	17” - 19” (43.2-48.3 cm)	$65	
NIB-1080	Large Privacy Flap for 20”-22” wide back supports	20” - 22” (50.8-55.9 cm)	$65	
NIB-1081	Extra Large Privacy Flap for 23” - 26” wide back supports	23” - 26” (58.4 - 66.0 cm)	$65	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

	MODEL	Description	RETAIL	✓
	NIB-0784	10” x 5” Large Rectangular	$195	
	NIB-0818	14” x 5” Extra-Large Rectangular	$195	

MultiFit™ HEAD SUPPORT – FIXTURES

	MODEL	Description	RETAIL	✓
	NIB-0747	Multi-Axis Hardware	$145	
	NIB-0823	Head Support Receiver Bracket	$58	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx^{3D} Cover – For a cover up to 18” x 18”	NXCCBF_____V	$125	
Extra Outer smartx^{3D} Cover – For a cover 18” x 20” and above	NXCCBF_____V	$145	
Extra Outer Incontinent Cover – For a cover up to 18” x 18”	NXCCAF_____V	$125	
Extra Outer Incontinent Cover – For a cover 18” x 20” and above	NXCCAF_____V	$145	

ALL PRICES QUOTED IN USD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com