160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T.+1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



6.8

VELOCE

ORDER FORM - ONTARIO

1- Purchasing information

P.O. number:* Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

2- Dealer/Representative information

Dealer/Representative:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:*

Name:

Address:

City: Province:

Postal Code:

Phone:

** Required*

MA080 ☐ *English Manual (STD)* *MA081* ☐ *Manuel Français*

MODEL									*Taxes not included*
WC20	☐	Ultralight VELOCE Wheelchair *(265 lb weight limit)*					WAMZ	WM0001598	**$5,195**
WC09	☐	**HD** Heavy Duty package *(350 lb weight capacity)*						WAMF	$205
Frame									
SEAT WIDTH [3]									
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐			NC
	☐	19" [2] ☐	20" [2] ☐	21" [7] ☐	22" ☐			WAMA	$365
SEAT DEPTH									
SD	☐	14" [1] ☐	15" ☐	16" ☐	17" ☐	18" ☐			NC
	☐	19" ☐	20" ☐					WAMB	$280
FRAME LENGTH *(Seat sling depth plus distance from seat sling to front of frame)* [4]									
SP01	☐	Standard frame length configuration [5]							STD
SP02	☐	Frame length 1" longer than standard configuration							NC
SP03	☐	Frame length 2" longer than standard configuration							NC
SP04	☐	Frame length 1" shorter than standard configuration							NC
SP05	☐	Frame length 2" shorter than standard configuration							NC
TRANSIT TIE DOWN									
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt [6] Short belt 20"-41"							$370
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt [6] Long belt 20"-46"							$370
AC35	☐	WC-19 Transit Tie Down without belt [6]							$345

1- *Only available with frame length 1" and 2" longer;* ***2-*** *Available with HD option; When HD option is selected, only the one piece angle adjustable flip-up footplate is available;*


4- *See chart on page 7; 5- 4.375" from seat sling to front frame; 6- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option;*

7- *Only available with HD option and the one piece adjustable angle flip-up*

3. Width calculator

FRONT Seat-to-floor height

	Caster sizes	Height available +- 1/4"	
☐	3" Caster	14" to 21"	**Specify desired height:**
☐	4" Caster	14.5" to 21"	
☐	5" Caster	15.5" to 21"	(inches)
☐	6" Caster	16.5" to 21" (N/A with 90° front frame)	



Frog Legs FRONT Seat-to-floor height

	Caster Sizes	Height available +- 1/4"	
☐	3" Caster	15.25'' to 21''	**Specify desired height:**
☐	4" Caster	15.75'' to 21''	
☐	5" Caster	16.75'' to 21''	(inches)
☐	6" Caster	17.75'' to 21'' (N/A with 90° front frame)	

FRONT CASTERS

CASTERS				Availability		
CA01-03 to CA01-06	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		NC
CA10	☐	Pneumatic casters (1 ¼")		6"		$170
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		$120
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		$198
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks				$525

Veloce - REAR Seat-to-floor Height

	WHEEL SIZES	Height available	
☐	22'' (501) wheel	13 1/2" to 18 1/4"	**Specify desired height:**
☐	24'' (540) wheel	14 1/2" to 19 1/4"	
☐	25'' (559) wheel	15" to 19 3/4"	(inches)
☐	26'' (590) wheel	15 1/2" to 20 1/4"	

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

REAR WHEELS

ONE ARM DRIVE [1] *Only available with Newton One rear wheels, permanent axles, 0° camber*						
HR22	☐	One arm drive - aluminum inner handrim [2] Left ☐ Right ☐		WAML		$1,480
HR60	☐	One arm drive - plastic coated inner handrim [2] Left ☐ Right ☐		WAML		$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS			**Available**			
RW01	☐	Mag [3]	22" 24"			$150
RW02	☐	Newton One - spoke wheel [4]	22" 24"			SF
RW03	☐	Newton Gravity - ultralight wheel	22", 24", 25", 26"			$390
RW04	☐	Spinergy Spox **Blk STD** ☐ Wht [5] ☐ Blue [5] ☐ Red [5] ☐ Yel [5] ☐	22", 24", 25", 26"			$915
RW05	☐	Spinergy LX **Blk STD** ☐ Pur [5] ☐ Wht [5] ☐ Blue [5] ☐ Red [5] ☐ Yel [5] ☐ Pink [5] ☐ Grn [5] ☐ Org [5] ☐	22", 24", 25", 26"			$1,135
RW07	☐	Spinergy CLX [6] **Blk STD** ☐ Pur [5] ☐ Wht [5] ☐ Blue [5] ☐ Red [5] ☐ Yel [5] ☐ Pink [5] ☐ Grn [5] ☐ Org [5] ☐	24", 25"			$1,690
TIRES *(The brand and color may vary depending on parts availability)*			**Available**			
		AIRLESS				
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black [9]	22" 24"			NC
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black [7,9]	22", 24", 25", 26"			$120
TI13	☐	Soft Urethane 1" - No tread, black - *Sentinel*	20",22", 24", 25", 26"			$150
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread [8]	22", 24", 26"			$235
		PNEUMATIC				
TI05	☐	Pneumatic 1 3/8" - Medium Tread [9]	22", 24", 26"			NC
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - *SpeedLite*	22", 24", 25", 26"			$130
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - *TrailBlazer*	24", 25", 26"			$210
TI03	☐	Pneumatic 2" - Large Tread, black [10]	22", 24", 25"			$140

1- *Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide;* ***2-*** *Only standard handrim position available;*
3- *Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position;* ***4-*** *Not compatible with Shox tires;* ***5-*** *color spokes require an additionnal 10 days lead-time;*
6- *Only available with 1" tires;* ***7-*** *22" is only available with Mags, Newton One and Newton Gravity wheels*
8 *- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels;* ***9-*** *Not available with flip-back armrest when the 6° camber is selected;* ***10-*** *Not available with armrests when the 3° or 6° camber are selected;*

HANDRIMS

HANDRIM					
HR01	☐	Aluminum anodized	22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized [2]	22", 24", 25"		$95
HR02	☐	Plastic coated [1]	22", 24", 25", 26"	WAMO	$180
HR03	☐	High friction coated [1,2,3]	22", 24", 25", 26"	WAMO	$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip [2] *See page 6 for the color selection and price*	22", 24", 25", 26"	WAMO	
HR28	☐	Surge [1,2,3,4] *(1 1/2" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"		$485
HR05	☐	Surge LT [1,2,3,4] *(1 1/8" Oval handrim w/ Gripton high friction strip)*	22", 24", 25", 26"		$515
HR04	☐	Natural Fit [1,2,3,4]	22", 24", 25", 26"		$320
HR62	☐	Simi H [1,2,3] *(1 3/8" Oval handrim with high friction strip)*	22", 24", 25", 26"		$485
HR63	☐	Nova H [1,2,3] *(1 1/4" Oval handrim with high friction strip)*	22", 24", 25", 26"		$515
HR64	☐	Optimum H [1,2,3] *(1 1/8" heart-shaped handrim)*	22", 24", 25", 26"		$320
HR99	☐	Omit handrim			NC

1- *Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position;* ***2-*** *Not available with Mag wheels;* ***3*** *- Not available with one arm drive;*
4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

HANDRIMS OPTION

ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	*Full tab*			STD
HR030	☐	Narrow Handrim position no cut [6]	*Full tab*			NC
HR30	☐	Narrow Handrim position with cut [5]	*Short tab 2 positions*			$29
HR31	☐	Super Narrow Handrim [5]	*Short tab 2 positions STD* ☐	*Short tab 1 position* [6] ☐		$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				$200
HR27	☐	Super grip				$230
HR099	☐	Omit thumb piece				NC


5- *Only available with Aluminum and Air Grip handrims;* *6- Not available with the combination of 24" Newton Gravity, Spinergy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;*

AXLE [1]					
AX01	☐	Quick release axle		WAMU	$100
AX03	☐	Quad release axle		WAMU	$280

1- *Not required with one arm drive*

REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM03	☐	3° camber			NC
CM06	☐	6° camber			$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =

Most stable												
1" ☐	1/4" ☐	1/2" ☐	3/4" ☐	**2"** ☐	1/4" ☐	1/2" ☐	3/4" ☐	**3"** ☐	1/4" ☐	1/2" ☐	3/4" ☐	**4"** ☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16'' and a center of gravity over 2.5''.


Possible conflict between a 26'' rear wheel and a center of gravity over 2.5''.

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK						
WL08	☐	Newton Wheel lock - push to lock				STD
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock *(works best with pneumatic or high pressure tire)*			WAMR	$200
WL13 WL14	☐	6" removable extension handle [1]				$75
WL05	☐	Scissors Lock				$185
WL15	☐	Unilateral wheel lock - pull to lock [2]	Left ☐	Right ☐	WAMN	$335
WL99	☐	Omit wheel locks				NC


1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilties and restriction, refer to the table at page 7 and the "How-to documents" page of our website:
https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents


FRONT FRAME ANGLE						
FA70	☐	70° Front frame angle				NC
FA75	☐	75° Front frame angle				NC
FA80	☐	80° Front frame angle				NC
FA85	☐	85° Front frame angle				NC
FA90	☐	90° Front frame angle				NC

SEAT TO FOOTREST LENGTH	
Set Length to: [1]	*Length must be at least 1 ¾" shorter than front seat-to-floor height*

1- +/- 1/2"



Footrest width / Taper selection		
Seat Width	**Straight** [1]	**1'' Taper**
Outside Measurement	**2 ½'' Narrower than Seat Width**	**4 ½'' Narrower than Seat Width**
14''	11.5'' ☐	9.5'' ☐
15''	12.5'' ☐	10.5'' ☐
16''	13.5'' ☐	11.5'' ☐
17''	14.5'' ☐	12.5'' ☐
18''	15.5'' ☐	13.5'' ☐
19''	16.5'' ☐	14.5'' ☐
20''	17.5'' ☐	15.5'' ☐
21''	18.5'' ☐	16.5'' ☐
22''	19.5'' ☐	17.5'' ☐



1- Not compatible with auto-folding footplate

FOOTPLATES						
FP06	☐	One piece adjustable angle flip-up - Aluminum		WAMJ	SEFND0015	STD
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts		WAMJ	SEFND0015	$480
FP28	☐	One piece auto-folding adjustable angle - Aluminum				$380
FP29	☐	One piece auto-folding adjustable angle - Carbon fiber with rubber inserts				$470
FOOTPLATE ACCESSORIES						
FP09	☐	Calf strap			SEFND0010	$62
FP30	☐	Bodypoint Aeromesh Padded calf strap				$145

SEATING AND BACKREST

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery			WAMD/WAME	$265
BH99	☐	Omit Back upholstery				STD
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion				$185
SC02	☐	3" Cushion				$205
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	
HEIGHT ADJUSTABLE WITH PUSH HANDLES OPTIONS *(MUST select ONE push handle options)*						
BC06-912	☐	Straight Carbon fiber back cane	9" - 12"	*increment of 1/2"*	Set at[2]:	NC
BC06-1215	☐	Straight Carbon fiber back cane	12" - 15"	*increment of 1/2"*	Set at[2]:	NC
BC06-1518	☐	Straight Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at[2]:	NC
BC06-1821	☐	Straight Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at[2]:	NC
BC08-1315	☐	8° Bend Carbon fiber back cane	13.5" - 15"	*increment of 1/2"*	Set at[2]:	$175
BC08-1518	☐	8° Bend Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at[2]:	$175
BC08-1821	☐	8° Bend Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at[2]:	$175
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS *(MUST select ONE push handle options)*						
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	*increment of 1/2"*	Set at[2]:	$330
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at[2]:	$330
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at[2]:	$330
BC10-1315	☐	8° bend Carbon fiber back cane	13.5" - 15"	*increment of 1/2"*	Set at[2]:	$380
BC10-1518	☐	8° bend Carbon fiber back cane	15" - 18"	*increment of 1/2"*	Set at[2]:	$380
BC10-1821	☐	8° bend Carbon fiber back cane	18" - 21"	*increment of 1/2"*	Set at[2]:	$380
		Select back cane angle:	85° [1] ☐ 90° ☐ 95° **STD** ☐ 100° ☐ 105° ☐ 110° ☐			
PUSH HANDLES OPTIONS *(ONE option at a time)*						
PH01	☐	Push handles				STD
PH05	☐	Fold down push handles [3,5]				$165
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" [3, 4, 5]				$345
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" [3, 4, 5]				$345
PH02	☐	Omit push handles				NC

1- *Not available with a 8° bend back cane;* ***2-*** *If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery;* ***3-*** *265lb weight limit*
4- *Not compatible with swing away armrests and 8° bend back canes and 9-12" back canes;* ***5-*** *Not compatible with headrest and attendant lock*

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST *(Convertible to "T" armrest)* [1]							
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		$405
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		$405
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		$405
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		$405
HEIGHT ADJUSTABLE "T" ARMREST [1]							
AR13-610	☐	10" Desk length armpad [4]	6" - 10"				$385
AR03-812	☐	10" Desk length armpad	8" - 12"				$385
AR03-1014	☐	10" Desk length armpad	10" - 14"				$385
AR14-610	☐	14" Full length armpad [4]	6" - 10"				$385
AR04-812	☐	14" Full length armpad	8" - 12"				$385
AR04-1014	☐	14" Full length armpad	10" - 14"				$385
HEIGHT ADJUSTABLE SWING AWAY ARMREST [2]							
AR05	☐	Tubular armrest - Short pad	8" - 11"				$325
AR23	☐	Tubular armrest - Long pad	8" - 11"				$325
AR06	☐	10" Desk length flat armpad	8" - 11"				$325
AR17	☐	14" Full length flat armpad	8" - 11"				$325
ARMRESTS OPTIONS							
AR12	☐	14" Full Length armpad [3] - *Gel Ovation*				SEAND0030	$160
AR99	☐	Omit armrest					NC

1- *May have a 1/2 " adjustment restriction on both limits of the range;* ***2-*** *Could interfere with some rigid backrests; Not available with fixed carbon fiber side guards*
3- *Only available with full length armrest; Not available with swing away and desk lenght armrests;* ***4-*** *May have some restriction with seat-to-floor height and rear wheel*

ACCESSORIES

SIDE GUARDS *(not available with Flip Back and "T" Armrest)*							
SG01	☐	Plastic side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐		WAMV	$235
SG02	☐	Carbon fiber side guards *(fixed)*				WAMV	$285
SG03	☐	Carbon fiber side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐		WAMV	$320
SG06	☐	Plastic fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐		WAMV	$340
SG04	☐	Carbon fiber fender side guards *(removable)*	Straight bracket - STD ☐	Offset bracket[1] ☐		WAMV	$535
SG99	☐	Omit side guard					NC

1- *only compatible with 0° camber;Increase width by 1/2" per side for a total of 1" ;*

POSITIONING BELT							
SB01-S	☐	Velcro adjustable - **Small 35"** [1]					$50
SB01-L	☐	Velcro adjustable - **Long 60"**					$50
SB02-S	☐	Auto buckle - **Small 35"** [1]					$80
SB02-M	☐	Auto buckle - **Medium 48"**					$80
SB02-L	☐	Auto buckle - **Long 60"**					$80
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - **Medium 56"**					$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - **Long 61"**					$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - **Small 51"**				SERND0030	$195
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - **Medium 56"**				SERND0030	$195
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - **Long 61"**				SERND0030	$195
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - **Small 51"**				SERND0020	$220
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - **Medium 56"**				SERND0020	$220
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - **Long 61"**				SERND0020	$220
SB03	☐	Hardware attachment for 2 extra points (1 pair)					$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - **Medium 18"**					$365
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - **Long 24"**					$365
SB99	☐	Omit positioning belt					STD

1- not recommended for seat width of 18" and more

HEADREST [1]							
AC01-S	☐	Small contoured head support pad (9"x5")					$215
AC01-M	☐	Large rectangular head support pad (10"x5")					$215
AC01-L	☐	XL rectangular head support pad (14"x5")					$215
AC28	☐	Headrest mounting hardware with stabilizer bar [1]					$365

1- *Only available with push handles*

ACCESSORIES							
WC10	☐	E-Motion adaptor Kit (2 adaptors) - *Only necessary with 6° camber*					$35
AC03	☐	Anti-tippers					$125
AC29	☐	Swing-away anti-tipper[1]	Left ☐	Right ☐			$190
AC27	☐	Tip assist [2]					$70
AC06	☐	Cane and crutch holder	Left ☐	Right ☐			$115
AC08	☐	Spoke guards [3]	Available: 22", 24", 25", 26"			WAMP	$165
AC10	☐	Seat pouch					$105
AC26	☐	BodyPoint Mobility Bag					$150
AC23	☐	Newton Backpack					$135
AC30	☐	Newton multi-purpose pouch					$80
AC12	☐	Velcro Closing Strap					$60
AC13	☐	Tool kit					$60

1- *Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate;* ***2-*** *Cannot be combined with Anti Tippers*

3- *Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel*

Colors

Options [1]	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	☐ $160		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	$350
Frame (NC)		☐ (STD)	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	TBD[3] ☐
Highlight [2]	☐		☐					☐		☐	☐			☐			☐				NC
FC99	☐	Custom Color	Specify PPG Paintit color:																	additional 10 days lead-time	$422

1- *Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;*

2- *To see different highlight, visit our website:* https://www.motioncomposites.com/en_ca/products/folding-wheelchairs/veloce

3- *The frame color will be confirmed with the quote conversion into an order (in absence of selection, the color will be black)*

Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD	N/A						
1							N/A
2						N/A	N/A
-1	N/A	N/A					
-2	N/A	N/A	N/A				

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater

Rear wheel: shall be 24 in and less

"T" or flip-back armrest only compatible with 0° camber

One arm drive: not compatible with "T"

One arm drive and flip-back armrest only available with 17 in depth frame

"T" or flip-back armrest may require to move the rear wheels outwards

Not compatible with wheel lock lever extension

Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.75 to 4
22 in	16+	2.5 to 4
24 in	16+	1 to 4

EVO CUSTOM

CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's subtle tweaks for comfort and mobility improvements or custom colors that make your accessories uniquely yours, we're here to support you.

Explore the **EVO Custom Catalog**.

You can include your request or reference an existing customization from our catalog in the notes section below.

Notes :

FOLDING ON CHAIR – CDN ORDER FORM

BioFit™ CUSHION



☐ Use wheelchair configuration

Or use this size:

NXCABF _____ (W" x D")

DEPTH \ WIDTH	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
12"	$750	$750	$750		N/A	N/A	N/A	N/A	N/A	N/A	N/A
13"						N/A	N/A	N/A	N/A	N/A	N/A
14"	$750	$750	$750								
15"		$750	$750	$750	$750						
16"		$750	$750	$750	$750	$750	$750	$786	$786	$786	$786
17"	N/A			$750	$750	$750	$750	$786	$786		
18"	N/A			$750	$750	$750	$750	$786	$786	$786	$786
19"	N/A			$750	$750	$750	$750	$786	$786		
20"	N/A			$750	$750	$750	$786	$786	$786	$786	$786

Option: Change cover for incontinent/infection control outer cover ☐ N/C

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of $160.00 and has a 2 to 4 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES



MODEL #	Description	Fits Back Width	Retail	QTY
RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	$96	
RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	$96	
RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	$96	

THORACIC SHAPE



Optima™

THORACIC BACK SUPPORT

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOACB1413	$870	☐	NIB-2193V	$156	☐
14"	16"	NOACB1416	$870	☐	NIB-2194V	$156	☐
15"	13"	NOACB1513	$870	☐	NIB-2197V	$156	☐
15"	16"	NOACB1516	$870	☐	NIB-2198V	$156	☐
16"	13"	NOACB1613	$870	☐	NIB-2201V	$156	☐
16"	16"	NOACB1616	$870	☐	NIB-2202V	$156	☐
16"	19"	NOACB1619	$870	☐	NIB-2203V	$156	☐
17"	13"	NOACB1713	$870	☐	NIB-2206V	$156	☐
17"	16"	NOACB1716	$870	☐	NIB-2207V	$156	☐
18"	13"	NOACB1813	$870	☐	NIB-2208V	$156	☐
18"	16"	NOACB1816	$870	☐	NIB-2211V	$156	☐
18"	19"	NOACB1819	$870	☐	NIB-2212V	$156	☐
19"	16"	NOACB1916	$870	☐	NIB-2213V	$156	☐
19"	19"	NOACB1919	$870	☐	NIB-2215V	$156	☐
20"	13"	NOACB2013	$870	☐	NIB-2218V	$156	☐
20"	16"	NOACB2016	$870	☐	NIB-2219V	$156	☐
20"	19"	NOACB2019	$870	☐	NIB-2220V	$156	☐
21"	16"	NOACB2116	$906	☐	NIB-2221V	$156	☐
21"	19"	NOACB2119	$906	☐	NIB-2222V	$156	☐
22"	13"	NOACB2213	$906	☐	NIB-2223V	$156	☐
22"	16"	NOACB2216	$906	☐	NIB-2224V	$156	☐
22"	19"	NOACB2219	$906	☐	NIB-2225V	$156	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™

THORACIC BACK SUPPORT
CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover &	Retail	✓
14"	13"	NOCFACB1413	$1,380	☐	NIB-2193V	$156	☐
14"	16"	NOCFACB1416	$1,380	☐	NIB-2194V	$156	☐
15"	13"	NOCFACB1513	$1,380	☐	NIB-2197V	$156	☐
15"	16"	NOCFACB1516	$1,380	☐	NIB-2198V	$156	☐
16"	13"	NOCFACB1613	$1,380	☐	NIB-2201V	$156	☐
16"	16"	NOCFACB1616	$1,380	☐	NIB-2202V	$156	☐
17"	13"	NOCFACB1713	$1,380	☐	NIB-2206V	$156	☐
17"	16"	NOCFACB1716	$1,380	☐	NIB-2207V	$156	☐
18"	13"	NOCFACB1813	$1,380	☐	NIB-2211V	$156	☐
18"	16"	NOCFACB1816	$1,380	☐	NIB-2212V	$156	☐
20"	13"	NOCFACB2013	$1,380	☐	NIB-2218V	$156	☐
20"	16"	NOCFACB2016	$1,380	☐	NIB-2219V	$156	☐

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.
For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

THORACIC SHAPE



Optima™ DEEP

THORACIC BACK SUPPORT

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOPCB1413	$990		NIB-2246V	$156	
14"	16"	NOPCB1416	$990		NIB-2247V	$156	
15"	13"	NOPCB1513	$990		NIB-2250V	$156	
15"	16"	NOPCB1516	$990		NIB-2251V	$156	
16"	13"	NOPCB1613	$990		NIB-2254V	$156	
16"	16"	NOPCB1616	$990		NIB-2255V	$156	
16"	19"	NOPCB1619	$990		NIB-2256V	$156	
17"	13"	NOPCB1713	$990		NIB-2259V	$156	
17"	16"	NOPCB1716	$990		NIB-2260V	$156	
18"	13"	NOPCB1813	$990		NIB-2264V	$156	
18"	16"	NOPCB1816	$990		NIB-2265V	$156	
18"	19"	NOPCB1819	$990		NIB-2266V	$156	
19"	16"	NOPCB1916	$990		NIB-2267V	$156	
19"	19"	NOPCB1919	$990		NIB-2268V	$156	
20"	13"	NOPCB2013	$990		NIB-2271V	$156	
20"	16"	NOPCB2016	$990		NIB-2272V	$156	
20"	19"	NOPCB2019	$990		NIB-2273V	$156	
21"	16"	NOPCB2116	$1,032		NIB-2274V	$156	
21"	19"	NOPCB2119	$1,032		NIB-2275V	$156	
22"	13"	NOPCB2213	$1,032		NIB-2276V	$156	
22"	16"	NOPCB2216	$1,032		NIB-2277V	$156	
22"	19"	NOPCB2219	$1,032		NIB-2278V	$156	



Optima Carbon™ DEEP

THORACIC BACK SUPPORT

CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover &	Retail	✓
14"	13"	NOCFPCB1413	$1,500		NIB-2246V	$156	
14"	16"	NOCFPCB1416	$1,500		NIB-2247V	$156	
15"	13"	NOCFPCB1513	$1,500		NIB-2250V	$156	
15"	16"	NOCFPCB1516	$1,500		NIB-2251V	$156	
16"	13"	NOCFPCB1613	$1,500		NIB-2254V	$156	
16"	16"	NOCFPCB1616	$1,500		NIB-2255V	$156	
17"	13"	NOCFPCB1713	$1,500		NIB-2259V	$156	
17"	16"	NOCFPCB1716	$1,500		NIB-2260V	$156	
18"	13"	NOCFPCB1813	$1,500		NIB-2264V	$156	
18"	16"	NOCFPCB1816	$1,500		NIB-2265V	$156	
20"	13"	NOCFPCB2013	$1,500		NIB-2271V	$156	
20"	16"	NOCFPCB2016	$1,500		NIB-2272V	$156	

THORACIC SHAPE



Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

xtend™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NHAACB(L)1213/17	12"	13-17"	$870		NIB-2148V	$156	
NHAACB(L)1413/17	14"	13-17"	$870		NIB-2149V	$156	
NHAACB(L)1613/17	16"	13-17"	$870		NIB-2150V	$156	
NHAACB(L)1813/17	18"	13-17"	$870		NIB-2151V	$156	
NHAACB (L)2013/17	20"	13-17"	$870		NIB-2152V	$156	
NHAACB (L)2213/17	22"	13-17"	$906		NIB-2153V	$156	



Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover &	Retail	✓
NHAPCB1616/20	16"	16"-20"	$870		NIB-2143V	$156	
NHAPCB1816/20	18"	16"-20"	$870		NIB-2144V	$156	
NHAPCB2016/20	20"	16"-20"	$870		NIB-2145V	$156	
NHAPCB2216/20	22"	16"-20"	$906		NIB-2146V	$156	

PELVIC SHAPE



Classic™ PELVIC BACK SUPPORT

Classic™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NACB1616	16"	16"	$870		NIB-2163V	$156	
NACB1618	16"	18"	$870		NIB-2164V	$156	
NACB1816	18"	16"	$870		NIB-2166V	$156	
NACB1818	18"	18"	$870		NIB-2167V	$156	
NACB1820	18"	20"	$870		NIB-2168V	$156	
NACB2020	20"	20"	$870		NIB-2171V	$156	
NACB2216	22"	16"	$906		NIB-2172V	$156	
NACB2218	22"	18"	$906		NIB-2173V	$156	
NACB2220	22"	20"	$906		NIB-2174V	$156	



Classic™ DEEP PELVIC BACK SUPPORT

Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NPCB1616	16"	16"	$996		NIB-2175V	$156	
NPCB1618	16"	18"	$996		NIB-2176V	$156	
NPCB1818	18"	18"	$996		NIB-2179V	$156	
NPCB1820	18"	20"	$996		NIB-2180V	$156	
NPCB2020	18"	20"	$996		NIB-2183V	$156	
NPCB2216	20"	16"	$1,032		NIB-2184V	$156	
NPCB2218	22"	18"	$1,032		NIB-2185V	$156	
NPCB2220	22"	20"	$1,032		NIB-2186V	$156	

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT ACCESSORIES



MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	$246	
NIB-1938	Extra Small Privacy Flap for 11” - 13” wide back supports	11” - 13” (27.9 - 33.0)	$84	
NIB-1078	Small Privacy Flap for 14”-16” wide back supports	14” - 16” (35.6-40.6 cm)	$84	
NIB-1079	Medium Flap for 17”-19” wide back supports	17” - 19” (43.2-48.3 cm)	$84	
NIB-1080	Large Privacy Flap for 20”-22” wide back supports	20” - 22” (50.8-55.9 cm)	$84	
NIB-1081	Extra Large Privacy Flap for 23” - 26” wide back supports	23” - 26” (58.4 - 66.0 cm)	$84	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

	MODEL	Description	RETAIL	✓
	NIB-0784	10” x 5” Large Rectangular	$246	
	NIB-0818	14” x 5” Extra-Large Rectangular	$246	

MultiFit™ HEAD SUPPORT – FIXTURES

	MODEL	Description	RETAIL	✓
	NIB-0747	Multi-Axis Hardware	$180	
	NIB-0823	Head Support Receiver Bracke	$72	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx^{3D} Cover – For a cover up to 18” x 18”	NXCCBF_____V	$156	
Extra Outer smartx^{3D} Cover – For a cover 18” x 20” and above	NXCCBF_____V	$180	
Extra Outer Incontinent Cover – For a cover up to 18” x 18”	NXCCAF_____V	$156	
Extra Outer Incontinent Cover – For a cover 18” x 20” and above	NXCCAF_____V	$180	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com